InstaFlo

US Bowel Catheter System – NON-STERILE

The InstaFlo Bowel Catheter is constructed primarily of silicone materials. All system components are latex-free.
Single patient use only.

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
REF 33004

USA: Rx Only


CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician or other healthcare practitioner licensed under state law to order this product.

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INTENDED USE

For diversion of liquid or semi-liquid stool to facilitate the collection of fecal matter in patients with little or no bowel control.

PRODUCT DESCRIPTION

Contains two main parts: the catheter and the collection bag. The insertion end of the catheter contains a retention cuff and the other end of the catheter has a twist lock fitting to attach the collection bag. The retention cuff leads to a drain tube that allows stool to drain directly from the rectum into the bag (**See Figure 1**). There are two catheter connectors attached to the catheter. The BLUE connector (**See Figure 2, Inset, Item #10**) is used to inflate and deflate the retention cuff. The CLEAR connector (**See Figure 2, Inset, Item #11**) is used only to irrigate the device when needed.

See Figure 2 and the following table for configuration and a description of the InstaFlo Bowel Catheter System.

No.	Description
1	Low-pressure Retention Cuff
1a	Collapse-Resistant Annulus
2	Collapsible Transsphincteric Zone
3	Coated Drain Tube
4	Catheter Connectors (See below)
5	Flush/Sampling Port
6	Drain Cap
7	Twist Lock Fitting
8	Closed Collection Bag
9	60 mL Luer Syringe

Catheter Connectors (see inset Figure 2)

No.	Color	Purpose	Label
10	BLUE	Inflation and deflation of the low-pressure retention cuff	CUFF (☒)= 35-40 mL H ₂ O
11	CLEAR with white-capped removable rigid Luer locked fitting and silicone funnel fitting. Both fittings have closure caps.	Allows for irrigation of the device	IRRIG

CONTRAINDICATIONS

- Do not use in patients with known sensitivities or allergies to the materials used in this device.
- Do not use if the patient's distal rectum cannot accommodate the inflated volume of the retention cuff, or if the distal rectum/anal canal is severely strictured (e.g., secondary to tumor, inflammatory condition, radiation injury, scarring).
- Do not use on patients having impacted stool.
- Do not use on patients with a recent (less than 6 weeks old) rectal anastomosis or a recent (less than 6 weeks old) anal or sphincter reconstruction.
- Do not use on patients with compromised rectal wall integrity (e.g., ischemic proctitis, mucosal ulcerations).

WARNINGS

(Failure to comply with the following warnings may result in patient injury).

- Do not use if package has been opened or damaged.
 - Do not use improper amount or type of fluids for irrigation or retention cuff inflation. NEVER use hot liquids.
 - Do not over inflate retention cuff.
 - Use only gravity or slow manual irrigation. Do not connect mechanical pumping devices to catheter irrigation connector (IRRIG).
 - Perform irrigations via the CLEAR connector (IRRIG) AND NOT via the BLUE connector (CUFF (☒)= 35-40 mL H₂O).
 - Blood per rectum should be investigated to ensure there is no evidence of pressure necrosis from the device. Discontinue use of the device if evident.
 - Abdominal distention that occurs while using the device should be investigated.
 - Excessive prolonged traction on the catheter, resulting in the retention cuff migration into the anal canal, could result in temporary or permanent clinical sphincter dysfunction or catheter expulsion.
- PRECAUTIONS**
- Do not sterilize.
 - The InstaFlo Catheter is not intended for use longer than 29 days.
 - The InstaFlo Catheter is not recommended for pediatric use.
 - To avoid damage to retention cuff, DO NOT contact with ANY sharp edge.
 - The InstaFlo Catheter may not be effective in individuals who have had their distal rectum significantly altered by surgical resection or reconstruction.
 - Patients with very weak sphincter function may expel the catheter under normal use, or may have increased leakage of stool compared to patients with normal sphincter function.
 - Caution should be observed in patients whose rectum may be altered by stricture due to radiation or affected with radiation proctitis.
 - Patients with severe tenesmus, or patients who experience tenesmus or severe pain after insertion of the device, may not tolerate the catheter in place.
 - Avoid inserting anything (e.g., thermometer, suppository, etc.) into the anal canal with the catheter in place to minimize the chance of patient injury or catheter damage.
 - Care should be taken when disconnecting the syringe from the CLEAR connector (IRRIG). Fluids may drain or splatter from the connector when it is disconnected.
 - Use WATER ONLY to inflate the retention cuff. Do not use saline solution, which may adversely affect valve function.
 - Do not use vigorous aspiration to remove fluid from the retention cuff. Vigorous aspiration may collapse the inflation lumen and/or pilot balloon, and may prevent deflation.
 - Do not allow ointments or lubricants having a petroleum base (e.g., Vaseline®, petroleum-based hand/body lotion) to contact the catheter. They may damage the silicone and compromise the integrity of the device.
 - Use only Hollister® branded bowel catheter collection bags with the InstaFlo Catheter.
 - Feces contain infectious material. Protect from splatter which may occur when disconnecting the collection bags or during catheter removal.
 - After use, this system is a biohazard. Handle and dispose of in accordance with institutional protocol and universal precautions for contaminated waste.

ADVERSE EVENTS

The following adverse events may be associated with the use of any rectal device:

- Perforation
- Pressure necrosis
- Loss of sphincter tone
- Obstruction
- Infection
- Excessive leakage of fecal contents

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









Figure/Figura/Figure 9

US InstaFlo™ Bowel Catheter System

Supplies needed in addition to the InstaFlo Bowel Catheter System Kit:

- Non-sterile gloves
- Water soluble lubricant
- Water (to inflate retention cuff)
- Catheter tip syringe, 60 mL

CATHETER INSERTION

1. The preferred patient position for catheter insertion is the left lateral knee-chest position (the patient's clinical situation may dictate the use of an alternate position).
2. Examine the rectum and anal canal for fecal impaction, formed and semi-formed stool, and anatomical abnormalities (masses, lesions, strictures, etc.). Do not use on patients having formed or impacted stool. If anatomical abnormalities are present, consult a physician prior to inserting the catheter.
3. Connect the end of the catheter drain tube to the collection bag by inserting the catheter twist lock fitting into the bag. Twist clockwise to lock in place. **WARNING: Prior to use of the collection bag, see Figure 3, Illustrations A-F.**
4. Using the Luer syringe provided, withdraw all air from the retention cuff via the BLUE connector (CUFF ① = 35-40 mL H₂O) then disconnect syringe and expel the air from the syringe. Confirm that the pilot balloon is collapsed (**See Figure 4**).
5. Fill the empty Luer syringe with 35-40 mL of water and connect it to the BLUE connector (CUFF ① = 35-40 mL H₂O).
6. Dispense water soluble lubricant into the kit tray or onto a flat surface.
7. Squeeze retention cuff flat and fold in half as shown in **Figures 6A and 6B**.
8. Grasping the catheter at the distal edge of the retention cuff, generously lubricate the end of the catheter. NOTE: To avoid damage to the retention cuff, DO NOT allow contact with ANY sharp edge.
9. Insert the retention cuff end of the catheter into the distal rectum as shown in **Figure 7**.
10. Fill the retention cuff by depressing the plunger on the water-filled Luer syringe, then disconnect the syringe. The pilot balloon should fill as the retention cuff is inflated (**See Figure 5**).
11. Use a gentle tug and release to seat the retention cuff (e.g., approximately 1 cm of the blue tubing should be visible **See Figure 8**).
12. Ensure that the catheter is not twisted.
13. Hang the collection bag at the foot of the bed so the tube is not twisted and the drainage bag is lower than the patient (**See Figure 1**).

ROUTINE CARE

- Frequently verify that the catheter and collection bag are positioned so that the catheter drain tube is not twisted, kinked or externally compressed.
- Frequently verify that waste is not accumulating in the drain tube. To remove accumulated waste:
 1. Hold the catheter close to the anus with one hand.
 2. Constrict the drain tube with the other hand.
 3. Strip the tubing toward the collection bag.
- Make sure that the patient is not lying on the drain tube or catheter connectors. This could cause discomfort or localized pressure.
- At least twice daily, flush the inside of the drain tube to remove any accumulated fecal matter:
 1. Using a catheter tip syringe, inject water into the flush/sampling port (**See Figure 2, Item #5**).
 2. Strip the tubing toward the collection bag.
- Small volumes of mucus or feces may leak onto the perianal region and can be managed with routine hygiene and absorbent pads.
- Replace the collection bag as needed (**See Figure 3**). Handle and dispose of in accordance with institutional protocol and universal precautions for contaminated waste.
- Verify the retention cuff volume at least every 7 days:
 1. Connect the Luer syringe to the BLUE connector (CUFF ① = 35-40 mL H₂O).
 2. Slowly aspirate all fluid from the retention cuff, then disconnect the syringe.
 3. Verify that the retention cuff is deflated by confirming that the BLUE connector (CUFF ① = 35-40 mL H₂O) pilot balloon is collapsed (**See Figure 4**).
 4. Refill the retention cuff via the BLUE connector (CUFF ① = 35-40 mL H₂O) with 35-40 mL of water, then disconnect the syringe.
 5. Verify that the BLUE connector (CUFF ① = 35-40 mL H₂O) pilot balloon indicates that the retention cuff is inflated (**See Figure 5**).

Irrigation of the Device

If stool is not flowing into the catheter, irrigate the catheter. Fill the Luer syringe with water, connect the syringe to the CLEAR connector (IRRIG) and slowly depress the plunger. **WARNING: Verify connection to correct catheter connector.** Irrigation is optional and used only to help reduce obstruction of the catheter tip with stool. If repeated irrigation with water does not return the flow of stool into the catheter, discontinue use of the InstaFlo System. **WARNING: Do not connect mechanical pumping devices to CLEAR connector (IRRIG).**

STOOL SAMPLING

1. Open the flush/sampling port cap and insert a catheter tip syringe into the flush/sampling port.
2. Draw an appropriate sample of fecal matter into the syringe.
3. Remove the syringe and close the flush/sampling port cap.

CATHETER REMOVAL

When the patient's stool consistency is no longer liquid or semi-liquid, remove the InstaFlo Catheter.

1. Deflate the retention cuff by connecting the Luer syringe to the BLUE connector (CUFF ① = 35-40 mL H₂O) and slowly aspirate all water. Then, disconnect the Luer syringe.
2. Verify that the retention cuff is deflated by confirming that the BLUE connector (CUFF ① = 35-40 mL H₂O) pilot balloon is collapsed (**See Figure 4**).
3. Apply water-soluble lubricant to the anal canal.
4. Grasp the catheter as close to the anus as possible, ask the patient to bear down (if capable), and apply steady traction to slide the catheter out of the anus. (**See Figure 9**). **The elastic nature of the device may result in abrupt exit of the catheter from the anus. Care should be taken to protect caregiver from splatter.**

NOTE: After use, handle and dispose of in accordance with institutional protocol and universal precautions for contaminated waste.

HELPFUL HINTS

Odor Management

- Flush the tubing with water using the flush/sampling port at least twice a day.
- Routinely cleanse any leakage from around the anus.

Leakage

- Make sure the retention cuff is inflated.
- Use a gentle tug and release to seat the retention cuff (e.g., approximately 1 cm of the blue tubing should be visible **See Figure 8**).

- Be sure the transsphincteric zone is not twisted, the tubing is straight and free of stool, and the drainage bag is lower than the patient.

- Excessive leakage during use may be secondary to catheter occlusion with stool. If leakage is excessive, consider discontinuing use of the InstaFlo System.

Catheter Expulsion

- Verify that no stool is present in the distal rectum.
- Rinse the catheter and reinsert, following the instructions for "Catheter Insertion".
- Verify that the catheter and collection bag are positioned so that the catheter drain tube is not twisted, kinked or externally compressed.
- If expulsion frequency is excessive, consider discontinuing use of the InstaFlo System.

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY
HOLLISTER INCORPORATED WARRANTS THAT IT EXERCISES REASONABLE CARE IN THE MANUFACTURE OF ITS PRODUCTS AND THEY ARE FREE FROM DEFECTS IN MATERIALS AND WORKMANSHIP UNDER NORMAL USE AND SERVICE. THERE IS NO EXPRESSED OR IMPLIED WARRANTY, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE HOLLISTER INCORPORATED PRODUCT(S) DESCRIBED IN THIS PUBLICATION. UNDER NO CIRCUMSTANCES SHALL HOLLISTER INCORPORATED BE RESPONSIBLE FOR MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES, OTHER THAN AS EXPRESSLY PROVIDED BY LAW, ARISING FROM THE PURCHASE, USE, REMOVAL, OR REPLACEMENT OF THE PRODUCT(S) DESCRIBED IN THIS PUBLICATION.

Description or specifications in Hollister Incorporated printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

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