

Hollister Wound Care

Hollister
Wound Care

HCPCS Codes Disclaimer

The reimbursement information provided herein is intended to provide general information concerning coding of Hollister products only. Hollister does not guarantee coverage or payment for any products. The ultimate responsibility for proper coding, satisfying reimbursement requirements, and obtaining reimbursement remains with the provider. Coding and coverage policies and guidelines are complex, can vary from one carrier or region to another, and are updated frequently. Providers should check with their local carriers or intermediaries often and should consult with counsel, a reimbursement specialist, and/or the DMECS website for any coding, coverage, reimbursement, or billing questions.

HCPCS codes shown in this Product Catalog are current as of 1-6-12.

Explanation of Third Party Code Reimbursement System (National Health Related Items Code)

To facilitate reimbursement of Hollister products, some Third Party Payment Programs require the use of an eleven-digit reimbursement number. Every Hollister manufactured product can be identified by an eleven-digit number. In general, the first five numbers identify the manufacturer's labeler code, and the last six numbers represent the product stock number. The labeler code for Hollister Incorporated is 08380. The labeler code for Hollister Wound Care LLC is 08567, and the labeler code for the Apogee product line is 08598. In order to utilize this coding system, follow these examples.

To code a four-digit stock number:

The Third Party Reimbursement Number for the Apogee 12 Fr Intermittent Catheter, Stock Number 1048 listed in the catalog on page 92, is: 08598-0010-48

08598 identifies Apogee. 0010-48 identifies the product Stock Number.

To code a five-digit stock number:

The Third Party Reimbursement Number for the New Image Drainable Pouch, Stock Number 18113 listed in the catalog on page 30, is: 08380-0181-13

08380 identifies Hollister Incorporated. 0181-13 identifies the product Stock Number.

To code a six-digit stock number:

The Third Party Reimbursement Number for the Restore Contact Layer Dressing with TRIACT Technology, Stock Number 509338 listed in the catalog on page 130, is: 08567-5093-38

08567 identifies Hollister Wound Care LLC. 5093-38 identifies the product Stock Number.

Additionally, three products in the Hollister Wound Care line require an NDC number due to their Dimethicone-based nature. The labeler code for their OTC drug products is 42828. The three NDC numbers are listed below.

NDC	Stock No.	Description
42828-9977-02	529977	Restore Cleanser & Moisturizer, 7.5 oz
42828-9978-03	529978	Restore Cleanser & Moisturizer, 11.5 oz
42828-9979-01	529979	Restore Dimethicreme Skin Protectant, 4 oz

Latex Statement

The specifications for the Hollister Continence Care products listed below include natural rubber latex as a component; our suppliers add natural rubber latex as part of the production process.

If a Hollister product(s) found in this catalog is not listed below, the specification for that product(s) and packaging does not include natural rubber latex as a component, and our component suppliers do not add natural rubber as part of their production process.

Current Stock Numbers:

1051	9107	9207	9835 §	B12RB
1067	9108	9208	B12R	B14RB
1084	9109	9209	B14R	B16RB
9106	9206	9834 §	B16R	

Discontinued Stock Numbers:

9606 *	9628 *	9656 †	9807 *	9818 †	9829 *
9607 *	9629 *	9657 †	9808 *	9824 ‡	9897 †
9608 *	9636 †	9658 †	9809 *	9825 ‡	9898 †
9609 *	9637 †	9659 †	9813 ‡	9826 *	
9626 *	9638 †	9803 ‡	9816 †	9827 *	
9627 *	9639 †	9806 *	9817 †	9828 *	

How to Use The Hollister Products Catalog

Hollister has created this simplified Product Catalog to help our customers find what they need in a quick and efficient manner.

If you...

know the stock number (or stock series number)

Go to the numeric index at the back of the catalog on pages 162–163 to find the page number where it appears.

If you...

know the product brand or series name

Go to the Table of Contents on pages 9-10, 86, 106, 114 and 120 to find the page number where it appears or reference the Product Name Index on pages 164–166.

If you...

know the HCPCS code

Go to the HCPCS Code Index at the back of the catalog on page 167.

† Product Discontinued as of July 2007

‡ Product Discontinued as of January 2009

* Product Discontinued as of September 2009

§ As of February 2009, this product has been redesigned to exclude natural rubber latex either as a component or from the production process. The packaging/insert of the new non-latex version of the product will not contain a latex statement. Please verify the latex status of the individual product prior to use by reviewing the packaging and product insert for a latex statement. A latex statement indicates the product is from the previous latex version.



Antimicrobial dressings vs. enzymatic debriders: Now you don't have to choose.

Hydrofera Blue antimicrobial dressings are highly absorbent and — unlike many other antimicrobial dressings — do not decrease the activity or effectiveness of enzymatic debriders or growth factors.^{1,2} This allows you to manage wound exudate and bioburden simultaneously with these therapies.



**The Choice is Clear:
Hydrofera Blue – When Managing BioBurden is a Priority.**

Hydrofera**BLUE**[®]

Exclusively Distributed by: **hollisterwoundcare**[®]

Complex Problems. Innovative Solutions.

www.hollisterwoundcare.com

¹ Data on file at Hydrofera, LLC

² Ramsey et al., Comparison of the Effects of Antimicrobial Wound Dressing on Cell Viability, Proliferation, and Growth Factor Activity, SAWC 2008

Caution: Federal law restricts this device for sale or on the order of a physician or licensed healthcare professional

Hydrofera Blue dressings are not indicated for third-degree burns and are not for use as an implant. If the dressing dries out, thoroughly moisten it with sterile saline or sterile water before removal. Regular evaluation and cleansing of wounds should be a common practice.

Table of Contents • Wound Care

General Information	121
HCPCS Codes Disclaimer.....	121
Explanation of Third Party Code Reimbursement System..... (National Health Related Items Code)	121
Medicare Utilization Chart	122–128
Hydrofera Blue Foam Dressings	129
Restore Wound Dressings with TRIACT Technology	130–132
Restore Contact Layer FLEX Dressing	130
Restore Contact Layer Dressing	130
Restore Contact Layer Dressing, Silver	130
Restore Adhesive Foam Dressing	131
Restore Non-Adhesive Foam Dressing	131
Restore Non-Adhesive Foam Dressing, Silver	131
Restore DUO Absorbent Dressing.....	131
Restore TRIO Absorbent Dressing, Adhesive	132
Restore Wound Dressings with TRIACT ADVANCED Technology	133
Restore Foam Dressing without Border	133
Restore Foam Border Dressing	133
Restore Heel Dressing without Border.....	133
Restore LITE Foam Dressing without Border.....	133
Restore Hydro-Shield Foam Dressings	134
Restore Hydro-Shield Foam Dressing without Border	134
Restore Hydro-Shield Foam Dressing	134
Restore Wound Dressings	135–136
Restore Calcium Alginate Dressing.....	135
Restore Calcium Alginate Dressing, Silver	135
Restore Hydrogel Dressing.....	135
Restore Hydrocolloid Dressing.....	136
Restore Extra Thin Hydrocolloid Dressing.....	136
Restore Hydrocolloid Dressing with Foam Backing	136
Restore Odor-Absorbent Dressing	136
Restore Wound Cleanser	137
Restore Wound Cleanser	137
ConnectEd Learning Management System	137
Skin Care	138
Restore Skin Cleanser	138
Restore Cleanser & Moisturizer	138
Restore DimethiCreme Skin Protectant.....	138
Restore Skin Conditioning Creme	138
Glossary • Wound Care	139–142

General Information

HCPCS Codes Disclaimer

The reimbursement information provided herein is intended to provide general information concerning coding of Hollister products only. Hollister does not guarantee coverage or payment for any products. The ultimate responsibility for proper coding, satisfying reimbursement requirements, and obtaining reimbursement remains with the provider. Coding and coverage policies and guidelines are complex, can vary from one carrier or region to another, and are updated frequently. Providers should check with their local carriers or intermediaries often and should consult with counsel, a reimbursement specialist, and/or the DMECS website for any coding, coverage, reimbursement, or billing questions.

HCPCS codes shown in this Product Catalog are current as of 1-6-12.

Explanation of Third Party Code Reimbursement System (National Health Related Items Code)

To facilitate reimbursement of Hollister Wound Care products, some Third Party Payment Programs require the use of an eleven-digit reimbursement number. Every Hollister Wound Care manufactured product can be identified by an eleven-digit number. The first five numbers identify the manufacturer's labeler code and the last six numbers represent the product stock number. The labeler code for Hollister Wound Care LLC is 08567. In order to utilize this coding system, simply combine the five digit labeler code with a product's six-digit stock number. You can follow this example.

To code a six-digit stock number:

The Third Party Reimbursement Number for Restore Contact Layer Silver Dressing with TRIACT Technology, Stock Number 509341, listed in the catalog on page 130 is:

08567-5093-41

08567 identifies Hollister Wound Care LLC, and 5093-41 identifies the product Stock Number.

Additionally, three products in the Hollister Wound Care line require an NDC number due to their Dimethicone-based nature. The labeler code for Hollister Wound Care LLC's OTC drug products is 42828. The three NDC numbers are listed below.

NDC	Stock No.	Description
42828-9977-02	529977	Restore Cleanser & Moisturizer, 7.5 oz
42828-9978-03	529978	Restore Cleanser & Moisturizer, 11.5 oz
42828-9979-01	529979	Restore Dimethicreme Skin Protectant, 4 oz

Medicare Utilization Chart

Medicare Usual Dressing Change Frequency for Wound Care Supplies

The quantity and type of dressings dispensed at any one time must take into account the current status of the wound(s), the likelihood of change, and the recent use of dressings. Dressing needs may change frequently (e.g., weekly) in the early phases of wound treatment and/or with heavily draining wounds. Surgical dressings must be tailored to the specific needs of an individual patient. The table below lists some coverage guidelines and usual dressing change frequencies. The medical necessity for more frequent dressing changes must be documented in the patient's medical record and submitted with the claim.

Per LCD for Surgical Dressings (L11471), NHIC effective Feb. 4, 2011.

HCPSC	Long Description	Frequency of Change	Wound Condition
A4450	Tape, non-waterproof, per 18 square inches	based on drsg size	
A4452	Tape, waterproof, per 18 square inches	based on drsg size	
A4461	Surgical dressing holder, non-reusable, each	*	
A4463	Surgical dressing holder, reusable, each	*	
A4465	Non-elastic binder for extremity	*	
A4490	Surgical stockings above knee length, each	*	
A4495	Surgical stockings thigh length, each	*	
A4500	Surgical stockings below knee length, each	*	
A4510	Surgical stockings full length, each	*	
A4649	Surgical supply; miscellaneous	*	
A6010	Collagen-based wound filler, dry form, sterile, per gram of collagen	*	
A6011	Collagen-based wound filler, gel/paste, per gram of collagen	*	
A6021	Collagen dressing, sterile, pad size 16 sq in or less, each	*	
A6022	Collagen dressing, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, each	*	
A6023	Collagen dressing, sterile, pad size more than 48 sq in, each	*	
A6024	Collagen dressing wound filler, per 6 inches	*	
A6025	Gel Sheet for dermal or epidermal application, (e.g., silicone, hydrogel, other), each	*	
A6154	Wound pouch, each	3/wk	
A6196	Alginate or other fiber gelling dressing, wound cover, sterile, pad size 16 sq in or less, each dressing	1/day	mod-high exudate
A6197	Alginate or other fiber gelling dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, each dressing	1/day	mod-high exudate
A6198	Alginate or other fiber gelling dressing, wound cover, sterile, pad size more than 48 sq in, each dressing	1/day	mod-high exudate

* No established usual maximum, refer to Medicare Utilization Guidelines.

HCPCS	Long Description	Frequency of Change	Wound Condition
A6199	Alginate or other fiber gelling dressing, wound filler, sterile, per 6 inches	2 units (12 in)/day	mod-high exudate
A6203	Composite dressing, sterile, pad size 16 sq in or less, with any size adhesive border, each dressing	3/wk	
A6204	Composite dressing, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, with any size adhesive border, each dressing	3/wk	
A6205	Composite dressing, sterile, pad size more than 48 sq in, with any size adhesive border, each dressing	3/wk	
A6206	Contact layer, sterile, 16 sq in or less, each dressing	1/wk	line the wound
A6207	Contact layer, sterile, more than 16 sq in but less than or equal to 48 sq in, each dressing	1/wk	line the wound
A6208	Contact layer, sterile, more than 48 sq in, each dressing	1/wk	line the wound
A6209	Foam dressing, wound cover, sterile, pad size 16 sq in or less, without adhesive border, each dressing	3/wk	mod-heavy exudate
A6210	Foam dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, without adhesive border, each dressing	3/wk	mod-heavy exudate
A6211	Foam dressing, wound cover, sterile, pad size more than 48 sq in, without adhesive border, each dressing	3/wk	mod-heavy exudate
A6212	Foam dressing, wound cover, sterile, pad size 16 sq in or less, with any size adhesive border, each dressing	3/wk	mod-heavy exudate
A6213	Foam dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, with any size adhesive border, each dressing	3/wk	mod-heavy exudate
A6214	Foam dressing, wound cover, sterile, pad size more than 48 sq in, with any size adhesive border, each dressing	3/wk	mod-heavy exudate
A6215	Foam dressing, wound filler, sterile, per gram	1/day	mod-heavy exudate
A6216	Gauze, non-impregnated, non sterile, pad size 16 sq in or less, without adhesive border, each dressing	3/day	
A6217	Gauze, non-impregnated, non sterile, pad size more than 16 sq in but less than or equal to 48 sq in, without adhesive border, each dressing	3/day	
A6218	Gauze, non-impregnated, non sterile, pad size more than 48 sq in, without adhesive border, each dressing	3/day	
A6219	Gauze, non-impregnated, sterile, pad size 16 sq in or less, with any size adhesive border, each dressing	1/day	
A6220	Gauze, non-impregnated, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, with any size adhesive border, each dressing	1/day	
A6221	Gauze, non-impregnated, sterile, pad size more than 48 sq in, with any size adhesive border, each dressing	1/day	

* No established usual maximum, refer to Medicare Utilization Guidelines.

Medicare Utilization Chart

HCPGS	Long Description	Frequency of Change	Wound Condition
A6222	Gauze, impregnated with other than water, normal saline, or hydrogel, sterile pad size 16 sq in or less, without adhesive border, each dressing	1/day	
A6223	Gauze, impregnated with other than water, normal saline, or hydrogel, sterile pad size more than 16 sq in but less than or equal to 48 sq in, without adhesive border, each dressing	1/day	
A6224	Gauze, impregnated with other than water, normal saline, or hydrogel, sterile pad size more than 48 sq in, without adhesive border, each dressing	1/day	
A6228	Gauze, impregnated with water or normal saline, sterile, pad size 16 sq in or less, without adhesive border, each dressing	*	
A6229	Gauze, impregnated with water or normal saline, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, without adhesive border, each dressing	*	
A6230	Gauze, impregnated with water or normal saline, sterile, pad size more than 48 sq in, without adhesive border, each dressing	*	
A6231	Gauze, impregnated, hydrogel, for direct wound contact, sterile, pad size 16 sq in or less, each dressing	1/day	min-no exudate
A6232	Gauze, impregnated, hydrogel, for direct wound contact, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, each dressing	1/day	min-no exudate
A6233	Gauze, impregnated, hydrogel, for direct wound contact, sterile, pad size more than 48 sq in, each dressing	1/day	min-no exudate
A6234	Hydrocolloid dressing, wound cover, sterile, pad size 16 sq in or less, without adhesive border, each dressing	3/wk	light-mod exudate
A6235	Hydrocolloid dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, without adhesive border, each dressing	3/wk	light-mod exudate
A6236	Hydrocolloid dressing, wound cover, sterile, pad size more than 48 sq in, without adhesive border, each dressing	3/wk	light-mod exudate
A6237	Hydrocolloid dressing, wound cover, sterile, pad size 16 sq in or less, with any size adhesive border, each dressing	3/wk	light-mod exudate
A6238	Hydrocolloid dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, with any size adhesive border, each dressing	3/wk	light-mod exudate
A6239	Hydrocolloid dressing, wound cover, sterile, pad size more than 48 sq in, with any size adhesive border, each dressing	3/wk	light-mod exudate
A6240	Hydrocolloid dressing, wound filler, paste, sterile, per ounce	3/wk	light-mod exudate
A6241	Hydrocolloid dressing, wound filler, dry form, sterile, per gram	3/wk	light-mod exudate

* No established usual maximum, refer to Medicare Utilization Guidelines.

HCPCS	Long Description	Frequency of Change	Wound Condition
A6242	Hydrogel dressing, wound cover, sterile, pad size 16 sq in or less, without adhesive border, each dressing	1/day	min-no exudate
A6243	Hydrogel dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, without adhesive border, each dressing	1/day	min-no exudate
A6244	Hydrogel dressing, wound cover, sterile, pad size more than 48 sq in, without adhesive border, each dressing	1/day	min-no exudate
A6245	Hydrogel dressing, wound cover, sterile, pad size 16 sq in or less, with any size adhesive border, each dressing	3/wk	min-no exudate
A6246	Hydrogel dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, with any size adhesive border, each dressing	3/wk	min-no exudate
A6247	Hydrogel dressing, wound cover, sterile, pad size more than 48 sq in, with any size adhesive border, each dressing	3/wk	min-no exudate
A6248	Hydrogel dressing, wound filler, gel, per fluid ounce	1/day up to 3 oz/mo	min-no exudate
A6250	Skin sealants, protectants, moisturizers, ointments, any type, any size	*	
A6251	Specialty absorptive dressing, wound cover, sterile, pad size 16 sq in or less without adhesive border, each dressing	1/day	mod-high exudate
A6252	Specialty absorptive dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, without adhesive border, each dressing	1/day	mod-high exudate
A6253	Specialty absorptive dressing, wound cover, sterile, more than 48 sq in, without adhesive border, each dressing	1/day	mod-high exudate
A6254	Specialty absorptive dressing, wound cover, sterile, pad size 16 sq in or less with any size adhesive border, each dressing	every other day	mod-high exudate
A6255	Specialty absorptive dressing, wound cover, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, with any size adhesive border, each dressing	every other day	mod-high exudate
A6256	Specialty absorptive dressing, wound cover, sterile, more than 48 sq in, with any size adhesive border, each dressing	every other day	mod-high exudate
A6257	Transparent film, sterile, 16 sq in or less, each dressing	3/wk	min exudate
A6258	Transparent film, sterile, more than 16 sq in but less than or equal to 48 sq in, each dressing	3/wk	min exudate
A6259	Transparent film, sterile, more than 48 sq in, each dressing	3/wk	min exudate
A6260	Wound cleanser, any type, any size	*	
A6261	Wound filler, gel/paste, per fluid ounce, not otherwise classified	1/day	
A6262	Wound filler, dry form, per gram, not otherwise classified	1/day	

* No established usual maximum, refer to Medicare Utilization Guidelines.

Medicare Utilization Chart

HCPGS	Long Description	Frequency of Change	Wound Condition
A6266	Gauze, impregnated, other than water, normal saline, or zinc paste, sterile, any width, per linear yard	1/day	
A6402	Gauze, non-impregnated, sterile, pad size 16 sq in or less, without adhesive border, each dressing	3/day	
A6403	Gauze, non-impregnated, sterile, pad size more than 16 sq in but less than or equal to 48 sq in, without adhesive border, each dressing	3/day	
A6404	Gauze, non-impregnated, sterile, pad size more than 48 sq in, without adhesive border, each dressing	3/day	
A6407	Packing strips, non-impregnated, sterile up to 2 inches in. width, per linear yard	3/day	
A6410	Eye pad, sterile, each	*	
A6411	Eye pad, non sterile, each	*	
A6412	Eye patch, occlusive, each	*	
A6413	Adhesive bandage, first-aid type, any size each	*	
A6441	Padding bandage, non-elastic, non-woven/non-knitted, width greater than or equal to three inches and less than five inches, per yard	same as ds chg freq	
A6442	Conforming bandage, non-elastic, knitted/woven, non sterile, width less than three inches, per yard	same as ds chg freq	
A6443	Conforming bandage, non-elastic, knitted/woven, non sterile, width greater than or equal to three inches and less than five inches, per yard	same as ds chg freq	
A6444	Conforming bandage, non-elastic, knitted/woven, non sterile, width greater than or equal to five inches, per yard	same as ds chg freq	
A6445	Conforming bandage, non-elastic, knitted/woven, sterile, width less than three inches, per yard	same as ds chg freq	
A6446	Conforming bandage, non-elastic, knitted/woven, sterile, width greater than or equal to three inches and less than five inches, per yard	same as ds chg freq	
A6447	Conforming bandage, non-elastic, knitted/woven, sterile, width greater than or equal to five inches, per yard	same as ds chg freq	
A6448	Light compression bandage, elastic, knitted/woven, width less than three inches, per yard	1/wk	
A6449	Light compression bandage, elastic, knitted/woven, width greater than or equal to three inches and less than five inches, per yard	1/wk	
A6450	Light compression bandage, elastic, knitted/woven, width greater than or equal to five inches, per yard	1/wk	

* No established usual maximum, refer to Medicare Utilization Guidelines.

HCPCS	Long Description	Frequency of Change	Wound Condition
A6451	Moderate-compression bandage, elastic, knitted/woven, load resistance of 1.25 to 1.34 foot pounds at 50% maximum stretch, width greater than or equal to three inches and less than five inches, per yard	1/wk	
A6452	High-compression bandage, elastic, knitted/woven, load resistance greater than or equal to 1.35 foot pounds at 50% maximum stretch, width greater than or equal to three inches and less than five inches, per yard	1/wk	
A6453	Self-adherent bandage, elastic, non-knitted/non-woven, width less than three inches, per yard	1/wk	
A6454	Self-adherent bandage, elastic, non-knitted/non-woven, width greater than or equal to three inches and less than five inches, per yard	1/wk	
A6455	Self-adherent bandage, elastic, non-knitted/non-woven, width greater than or equal to five inches, per yard	1/wk	
A6456	Zinc paste-impregnated bandage, non-elastic, knitted/woven, width greater than or equal to three inches and less than five inches, per yard	1/wk	
A6457	Tubular dressing with or without elastic, any width, per linear yard	*	
A6501	Compression burn garment, bodysuit (head to foot), custom fabricated	*	
A6502	Compression burn garment, chin strap, custom fabricated	*	
A6503	Compression burn garment, facial hood, custom fabricated	*	
A6504	Compression burn garment, glove to wrist, custom fabricated	*	
A6505	Compression burn garment, glove to elbow, custom fabricated	*	
A6506	Compression burn garment, glove to axilla, custom fabricated	*	
A6507	Compression burn garment, foot to knee length, custom fabricated	*	
A6508	Compression burn garment, foot to thigh length, custom fabricated	*	
A6509	Compression burn garment, upper trunk to waist including arm openings (vest), custom fabricated	*	
A6510	Compression burn garment, trunk, including arms down to leg openings (leotard), custom fabricated	*	
A6511	Compression burn garment, lower trunk including leg openings (panty), custom fabricated	*	
A6512	Compression burn garment, not otherwise classified	*	
A6513	Compression burn mask, face and/or neck, plastic or equal, custom fabricated	*	

* No established usual maximum, refer to Medicare Utilization Guidelines.

Medicare Utilization Chart

HCPGS	Long Description	Frequency of Change	Wound Condition
A6530	Gradient compression stocking, below knee, 18-30 mm HG, each	*	open venous ulcer
A6531	Gradient compression stocking, below knee, 30-40 mm HG, each	*	open venous ulcer
A6532	Gradient compression stocking, below knee, 40-50 mm HG, each	*	open venous ulcer
A6533	Gradient compression stocking, thigh length, 18-30 mm HG, each	*	open venous ulcer
A6534	Gradient compression stocking, thigh length, 30-40 mm HG, each	*	open venous ulcer
A6535	Gradient compression stocking, thigh length, 40-50 mm HG, each	*	open venous ulcer
A6536	Gradient compression stocking, full length/chap style, 18-30 mm HG, each	*	open venous ulcer
A6537	Gradient compression stocking, full length/chap style, 30-40 mm HG, each	*	open venous ulcer
A6538	Gradient compression stocking, full length/chap style, 40-50 mm HG, each		open venous ulcer
A6539	Gradient compression stocking, waist length, 18-30 mm HG, each	*	open venous ulcer
A6540	Gradient compression stocking, waist length, 30-40 mm HG, each	*	open venous ulcer
A6541	Gradient compression stocking, waist length, 40-50 mm HG, each	*	open venous ulcer
A6544	Gradient compression stocking, garter belt	*	open venous ulcer
A6545	Gradient compression wrap, non-elastic, below knee, 30-50 mm HG, each	1/6 mo	open venous ulcer
A6549	Gradient compression stocking/sleeve, not otherwise specified	*	open venous ulcer
A9720	Non-covered item or service		

* No established usual maximum, refer to Medicare Utilization Guidelines.

Hydrofera Blue Foam Dressings



**NOW
AVAILABLE**
exclusively from
**Hollister
Wound Care**

- ▶ Highly absorptive bacteriostatic foam dressing
- ▶ Broad spectrum activity against microorganisms commonly found in wounds, including MRSA and VRE
- ▶ Compatible for use with enzymatic debriders
- ▶ Does not harm growth factors

Standard Dressings:

Stock No	HCPCS	Description	Quantity/Box
HB2214	A6209	2" x 2"	10
HB4414	A6209	4" x 4"	10
HB6614	A6210	6" x 6"	10
HBT0906	A6215	9 mm Tunneling Dressing	10 pair

Heavy Drainage Dressings

Stock No	HCPCS	Description	Quantity/Box
HBHD4450	A6209	4" x 4" x 1/2" thick	5
HBHD6675	A6210	6" x 6" x 3/4" thick	5

Dressings with Moisture-Retentive Film Backing

Stock No	HCPCS	Description	Quantity/Box
HBRF4414	A6209	4" x 4"	10
HBRF2258	A6209	2-1/4" x 8"	10
HBRF4475	A6212	4" x 4-3/4" Island Dressing (pad size 2" x 2-3/4")	10
HBRF2650	A6209	2-1/2" diameter Ostomy Dressing	10

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or licensed healthcare professional.

Restore Wound Dressings with TRIACT Technology

TRIACT Technology is an innovative formulation in *select* Restore wound dressings. The unique composition of TRIACT Technology provides a moist wound interface and allows virtually pain-free removal. It also reduces trauma to delicate tissue upon removal. Dressings with TRIACT Technology are designated with this symbol on the packaging.



Restore Contact Layer FLEX Dressing with TRIACT Technology: Sterile



- ▶ **Flexible** – Comfortable for patient during application and while in place
- ▶ **Highly Conformable** – Can be shaped to fit contours of deep or irregular wounds, and in difficult anatomical locations
- ▶ **Non-Adherent** – Virtually pain-free dressing removal for enhanced patient comfort
- ▶ **Maintains a Moist Wound Interface** – Forms a gel in the presence of exudate
- ▶ **Minimizes Trauma to the Wound Bed** – The dressing's gentle tack helps keep it in place over the wound
- ▶ **Versatile** – Can be used on a wide variety of wound types, and in conjunction with NPWT

Stock No	HCPCS	Description	Quantity/Box
506487	A6206	2" x 2"	10
506488	A6207	4" x 5"	10
506489	A6207	6" x 8"	10

Restore Contact Layer Dressing with TRIACT Technology: Sterile



- ▶ Non-Adhesive, Non-Occlusive Wound Contact Dressing
- ▶ Indicated For Low to Moderate Exuding Partial- and Full-Thickness Wounds
- ▶ Ideal to Use on Wounds with Fragile Surrounding Skin
- ▶ Apply Directly to the Wound and Cover with a Secondary Dressing

Stock No	HCPCS	Description	Quantity/Box
509378	A6206	2" x 2"	10
509338	A6207	4" x 5"	10
509339	A6207	6" x 8"	10
509379	A6208	8" x 12"	10

Restore Contact Layer Dressing with Silver and TRIACT Technology: Sterile



- ▶ Non-Adhesive, Non-Occlusive Antimicrobial Wound Contact Dressing
- ▶ May Help to Reduce Infection in Low to Moderately Exuding Partial- and Full-Thickness Wounds
- ▶ Shown to be Effective Against Bacteria Most Frequently Associated with Wound Infections
- ▶ Apply Directly to the Wound and Cover with a Secondary Dressing
- ▶ Contains Silver Sulfate

Stock No	HCPCS	Description	Quantity/Box
509341	A6207	4" x 5"	10
509342	A6207	6" x 8"	10

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or licensed healthcare professional.

Restore Adhesive Foam Dressing with TRIACT Technology: Sterile



- ▶ Semi-Permeable, Absorbent Dressing with Adhesive Border
- ▶ Manages Moderate to Heavily Exuding Chronic and Acute Wounds
- ▶ Super-Absorbent Foam Pad Ensures Drainage of Exudates and Helps Protect the Skin Around the Lesion from Maceration
- ▶ Soft, Semi-Permeable Backing

Stock No	HCPCS	Description	Quantity/Box
509347	A6212	4" x 4" (1.5" x 1.5" Foam Pad)	10
509348	A6212	6" x 6" (4" x 4" Foam Pad)	10
509349	A6213	6" x 8" (3.5" x 5.5" Foam Pad)	10

Restore Non-Adhesive Foam Dressing with TRIACT Technology: Sterile



- ▶ Non-Adhesive Absorbent Dressing
- ▶ Indicated for the Treatment of All Types of Moderate to Heavily Exuding Chronic and Acute Wounds
- ▶ Super-Absorbent Foam Pad Ensures Drainage of Exudates and Helps Protect the Skin Around the Lesion from Maceration
- ▶ Soft, Semi-Permeable Backing

Stock No	HCPCS	Description	Quantity/Box
509343	A6209	4" x 4"	10
509344	A6210	6" x 8"	10

Restore Non-Adhesive Foam Dressing with Silver and TRIACT Technology: Sterile*



- ▶ Non-Adhesive, Non-Occlusive Antimicrobial Absorbent Dressing
- ▶ May Help to Reduce Infection in Moderately to High Exuding Partial- and Full-Thickness Wounds
- ▶ Super-Absorbent Foam Pad Ensures Drainage of Exudates and Helps Protect the Skin Around the Lesion from Maceration
- ▶ Soft, Semi-Permeable Backing
- ▶ Contains Silver Sulfate



Stock No	HCPCS	Description	Quantity/Box
509345	A6209	4" x 4"	10
509346	A6210	6" x 8"	10

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or licensed healthcare professional.

Restore DUO Absorbent Dressing with TRIACT Technology: Sterile



- ▶ Absorbent, Non-Adhesive Dressing that Incorporates TRIACT Non-Adherent Layer
- ▶ Flexible, Comfortable, and Conformable for Versatile Use
- ▶ Indicated for Low Exuding Acute and Chronic Wounds Including Skin Tears
- ▶ Apply Directly to the Wound and Cover with Secondary Dressing

Stock No	HCPCS	Description	Quantity/Box
509375	A6234	2" x 4"	10
509376	A6235	4" x 5"	10
509377	A6235	6" x 8"	10

Restore Wound Dressings with TRIACT Technology

Restore TRIO Absorbent Dressing, Adhesive, with TRIACT Technology: Sterile




- ▶ All-in-one Absorbent Dressing with Adhesive Border that Incorporates TRIACT Non-Adherent Layer
- ▶ Flexible, Comfortable, and Conformable for Versatile Use
- ▶ Indicated for Low Exuding Acute and Chronic Wounds Including Dermabrasions and Minor Cuts and Burns
- ▶ Apply Directly to The Wound — No Need for a Secondary Dressing

Stock No	HCPCS	Description	Quantity/Box
509372	A6237	3" x 3" with pad size 1.8" x 1.8"	10
509373	A6237	4" x 5" with pad size 2.3" x 3.2"	10
509374	A6238	6" x 8" with pad size 3.4" x 6"	10

Did you know...

A Quick Guide to Restore Silver Packaging



Restore Contact Layer and Foam Dressings with Silver have two-sided foil packaging. All other **Restore** products with TRIACT Technology have one clear side and one foil side.

Restore Wound Dressings with TRIACT ADVANCED Technology

Hollister Wound Care features a line of foam dressings featuring TRIACT ADVANCED Technology. These products offer all the advantages of TRIACT Technology — including virtually pain-free removal and maintaining a moist wound interface — with these additional benefits:



- ▶ Stays in Place Due to Micro-Adherent Technology Facilitating Easier Application and Reduced Risk of Maceration
- ▶ Excellent Conformability
- ▶ Available in Several Types of Dressings

Restore Foam Dressing without Border Featuring TRIACT ADVANCED Technology



- ▶ Non-Adhesive Absorbent Dressing
- ▶ Manages Moderate to Moderately Heavily Exuding Acute and Chronic Wounds

Stock No	HCPCS	Description	Quantity/Box
509380	A6209	2.5" x 2.5"	10
509381	A6209	4" x 4"	10
509382	A6210	6" x 6"	10



509384	A6210	Heel 4.7" x 7.5"	10
--------	-------	------------------	----

Restore Foam Border Dressing Featuring TRIACT ADVANCED Technology



- ▶ Absorbent Dressing with Adhesive Border
- ▶ Manages Moderate to Moderately Heavily Exuding Acute and Chronic Wounds

Stock No	HCPCS	Description	Quantity/Box
509387	A6212	4" x 4" with pad size 2" x 2"	10
509388	A6212	6" x 6" with pad size 4" x 4"	10
509389	A6213	6" x 8" with pad size 3.3" x 5.8"	10
509390	A6213	8" x 8" with pad size 5" x 5"	10

Restore LITE Foam Dressing without Border Featuring TRIACT ADVANCED Technology



- ▶ Non-Adhesive Absorbent Dressing is Thinner Than the Standard Foam Dressings
- ▶ Manages Moderate to Light to Moderately Exuding Acute and Chronic Wounds

Stock No	HCPCS	Description	Quantity/Box
509393	A6209	2.5" x 2.5"	10
509394	A6210	4" x 5"	10
509395	A6210	6" x 6"	10
509396	A6210	6" x 8"	10

Restore Hydro-Shield Foam Dressings

Restore Hydro-Shield Foam Dressing without Border



- ▶ Waterproof and Bacterial Barrier Protection
- ▶ Excellent Absorption and Retention

Stock No	HCPCS	Description	Quantity/Box
529300	A6209	3" x 3"	10
529301	A6209	4" x 4"	10
529302	A6210	6" x 6"	10
529303	A6211	8" x 8"	10

Restore Hydro-Shield Foam Dressing with Border



- ▶ Waterproof and Bacterial Barrier Protection
- ▶ Excellent Absorption and Retention

Stock No	HCPCS	Description	Quantity/Box
529304	A6212	4" x 4" with pad size 2" x 2"	10
529305	A6213	6" x 6" with pad size 4.5" x 4.5"	10
529306	A6213	8" x 8" with pad size 5.4" x 5.4"	10

Restore Wound Dressings

Restore Calcium Alginate Dressing: Sterile



- ▶ Highly Absorptive
- ▶ Made From Calcium Sodium Alginate, Creating Conformable Protective Gel When in Contact with Exudate
- ▶ Forms Moist Wound Interface
- ▶ Remains Intact When Saturated so it is Removed in One Piece

Stock No	HCPCS	Description	Quantity/Box
529938	A6196	2" x 2"	10
529937	A6196	4" x 4"	10
529939	A6197	4" x 8"	5
529940	A6199	12" Rope	5

Restore Calcium Alginate Dressing Silver: Sterile



- ▶ Releases Silver Ions in the Presence of Wound Exudate
- ▶ Effective Antimicrobial Agent Against a Broad Spectrum of Microorganisms
- ▶ Indicated for Moderate to Heavily Exuding Partial- to Full-Thickness Wounds
- ▶ Assists in Maintaining a Moist Environment and Allows Intact Removal
- ▶ Contains Silver Sodium Hydrogen Zirconium Phosphate

Stock No	HCPCS	Description	Quantity/Box
529967	A6196	2" x 2"	10
529968	A6197	4" x 4.75"	10
529969	A6199	1" x 12" Rope	5

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or licensed healthcare professional.

Restore Hydrogel Dressing



- ▶ Maintains Moist Healing Environment in Partial- and Full-Thickness Wounds
- ▶ Clear, Viscous Gel
- ▶ Up to 72-hour Wear Time
- ▶ Available in Three Forms

Stock No	HCPCS	Description	Quantity/Box
529974	A6248	Amorphous Gel in Tube 3 oz — Non Sterile	12
529972	A6231	Impregnated Gauze Sponge 4" x 4" — Sterile	15
529973	A6266	Impregnated Gauze Strip 2" x 3.5 yd — Sterile	12

Restore Wound Dressings (continued)

Restore Hydrocolloid Dressing: Sterile



- ▶ Conformable Semi-occlusive Dressing with Flexible Outer Layer
- ▶ Maintains Moist Wound Environment
- ▶ Manages Minimal to Moderate Exudate in Partial- and Full-Thickness Wounds
- ▶ Available With or Without Tapered Edges, Specialty Shapes
- ▶ Includes Disposable Wound Measuring Guide

Stock No	HCPCS	Description	Quantity/Box
519953	A6234	Without Tapered Edges 4" x 4"	5
519954	A6235	Without Tapered Edges 6" x 8"	3
519955	A6236	Without Tapered Edges 8" x 8"	3
519956	A6234	With Tapered Edges 4" x 4"	5
519963	A6235	With Tapered Edges 6" x 6"	5
519957	A6235	With Tapered Edges 6" x 8"	3
519958	A6236	With Tapered Edges 8" x 8"	3
519959	A6235	With Tapered Edges/Triangle Shaped 17 sq in	5
519965	A6235	With Tapered Edges/Triangle Shaped 26.5 sq in	5

Restore Extra Thin Hydrocolloid Dressing: Sterile



- ▶ Thinner Semi-occlusive Dressing with Flexible Backing
- ▶ Maintains Moist Wound Environment
- ▶ Manages Minimal Exudate in Shallow Wounds
- ▶ Protects Skin from Friction Injury
- ▶ Includes Disposable Wound Measuring Guide

Stock No	HCPCS	Description	Quantity/Box
519921	A6234	4" x 4"	5
519923	A6235	6" x 8"	3
519925	A6236	8" x 8"	3

Restore Hydrocolloid Dressing with Foam Backing: Sterile



- ▶ Occlusive Dressing with White Foam Backing
- ▶ For Low Friction Areas with Moderate Exudate
- ▶ Includes Disposable Wound Measuring Guide

Stock No	HCPCS	Description	Quantity/Box
519930	A6234	4" x 4"	5
519932	A6235	6" x 8"	3
519935	A6236	8" x 8"	3

Restore Odor-Absorbent Dressing: Sterile



- ▶ Foam Matrix Dressing Covered with Non-Woven Fabric for Patients with Malodorous Wounds Resulting from Necrosis, Infection, or Bacterial Contamination
- ▶ Activated Carbon Helps to Neutralize Offensive Odors
- ▶ Use as Primary or Secondary Dressing

Stock No	HCPCS	Description	Quantity/Box
519904	A6209	4" x 4"	10
519906	A6211	6" x 10"	10

Did you know...



Disposable Wound Guides

Disposable Wound Measuring Guides are included with all **Restore** Hydrocolloid Dressings. To order additional Wound Measuring Guides, please see page 158.

Restore Wound Cleanser

Restore Wound Cleanser



- ▶ Helps Remove Particulate Matter and Other Debris From Wound Bed
- ▶ No-Rinse Alternative to Saline and Harsh Antiseptic Cleansers

Stock No	HCPCS	Description	Quantity/Box
529975	A6260	8 oz pump spray	12
529976	A6260	12 oz pump spray	12



Wound Care Education

- Web-Based Learning Management System
- Customizable Tracking and Reporting
- In-Service Videos
- Online Tools

Education Modules include:

- ▶ Wound Care Foundation Courses
 - Moist Wound Healing
 - Wound Assessment
 - Wound Documentation
- ▶ Pressure Ulcer Program
 - Staging
 - Prevention
 - Knowledge Assessment
- ▶ Skin Tear Courses
- ▶ Vascular Ulcer Courses
- ▶ Pain
- ▶ Wound Bioburden

Some Modules available in Spanish

Skin Care

Restore Skin Cleanser



- ▶ For Use on Intact Skin
- ▶ Cleanses Peristomal and Perineal Skin of Urinary and Fecal Discharge
- ▶ Gentle, Non-Irritating
- ▶ Must Be Rinsed

Stock No	HCPCS	Description	Quantity/Box
517210	A9270	8 oz (237 mL) Pump Spray	12

Restore Cleanser & Moisturizer



- ▶ Effective Skin Protection for People with Urinary and/or Fecal Incontinence
- ▶ Cleanser, Moisturizer, Protective Breathable Barrier for Skin
- ▶ Dimethicone-Based
- ▶ No Rinse

Stock No	HCPCS	Description	Quantity/Box
529977	A6250	7.5 oz (222 mL) Pump Spray	12
529978	A6250	11.5 oz (340 mL) Pump Spray	12

Restore DimethiCreme Skin Protectant



- ▶ Dimethicone-Based Protective Breathable Barrier
- ▶ Use on Skin or Mucous Membrane Surfaces Exposed to Harmful Stimuli
- ▶ Moisturizes and Protects
- ▶ Not for Use Under Adhesive Products

Stock No	HCPCS	Description	Quantity/Box
529979	A6250	4 oz (118 mL) Tube	12

Restore Skin Conditioning Creme



- ▶ Moisturizes and Soothes Dry, Flaky Skin
- ▶ May be Used Under Adhesive Products

Stock No	HCPCS	Description	Quantity/Box
517220	A6250	4 oz (120 mL) Bottle	12

Glossary • Wound Care

Abrasion: Wearing away of the skin through some mechanical process (friction or trauma)

Abscess: Accumulation of pus formed in tissue as a result of infection

Alginate: A highly absorptive dressing derived from brown seaweed

Cellulitis: Inflammation of the tissues indicating a local infection; characterized by redness, edema, and tenderness

Collagen: Main supportive protein of the skin and connective tissue

Debridement: Removal of foreign material and devitalized or contaminated tissue from a wound

Dehiscence: Separation of wound edges

Denude: Removal or loss of superficial skin layers

Edema: Swelling

Epidermis: Outermost layer of the skin

Erythema: Diffuse redness of the skin

Eschar: Thick, leathery black crust; it is nonviable tissue and is colonized with bacteria

Excoriation: Linear scratches on the skin

Exudate: Wound fluid or drainage

Friction: Rubbing that causes mechanical trauma to the skin

Full-thickness: Tissue destruction extending through the dermis to involve subcutaneous level and possibly muscle, fascia, or bone

Granulation: Formation of connective tissue and many new capillaries in a full-thickness wound; typically appears as red and cobblestoned

Hydrocolloid dressing: A category of wound dressings composed of materials, such as gelatin, pectin, and carboxymethylcellulose, that provide a moist healing environment and adhere to the skin around the wound

Hydrogel: Water- or glycerin-based gels, impregnated gauzes, or sheet dressings. Hydrogels maintain a moist healing environment and absorb a minimal amount of wound exudate

Hydrophilic: Attracting moisture

Infection: Overgrowth of microorganisms in sufficient quantities to overwhelm the body's defenses

Maceration: Softening of tissue by soaking in fluids

Necrotic: Devitalized tissue; may appear yellow and moist, gray, or dark and leathery

Partial-thickness: Wounds that extend through the epidermis and may involve the dermis; these wounds heal by re-epithelialization

Peri-wound: The area immediately around the wound

Pus: Thick fluid composed of leukocytes, bacteria, and cellular debris


Shear: Sliding of skin over subcutaneous tissues and bones obstructing cutaneous capillaries, which may lead to ischemia

Pressure Ulcer Stages:

Suspected Deep Tissue Injury: Purple or maroon localized area of discolored intact skin or blood-filled blister due to damage of underlying soft tissue from pressure and/or shear. The area may be preceded by tissue that is painful, firm, mushy, boggy, warmer, or cooler as compared to adjacent tissue.



Deep tissue injury may be difficult to detect in individuals with dark skin tones. Evolution may include a thin blister over a dark wound bed. The wound may further evolve and become covered by thin eschar. Evolution may be rapid exposing additional layers of tissue even with optimal treatment.



Stage I: Intact skin with non-blanchable redness of a localized area usually over a bony prominence. Darkly pigmented skin may not have visible blanching; its color may differ from the surrounding area.

The area may be painful, firm, soft, warmer, or cooler as compared to adjacent tissue. Stage I may be difficult to detect in individuals with dark skin tones. May indicate “at risk” persons (a heralding sign of risk).

Stage II: Partial-thickness loss of dermis presenting as a shallow open ulcer with a red pink wound bed, without slough. May also present as an intact or open/ruptured serum-filled blister.



Presents as a shiny or dry shallow ulcer without slough or bruising. This stage should not be used to describe skin tears, tape burns, perineal dermatitis, maceration, or excoriation. Bruising indicates suspected deep tissue injury.



Stage III: Full-thickness tissue loss. Subcutaneous fat may be visible but bone, tendon, or muscle are not exposed. Slough may be present but does not obscure the depth of tissue loss. May include undermining and tunneling.

The depth of a stage III pressure ulcer varies by anatomical location. The bridge of the nose, ear, occiput, and malleolus do not have subcutaneous tissue and stage III ulcers can be shallow. In contrast, areas of significant adiposity can develop extremely deep stage III pressure ulcers. Bone/tendons is not visible or directly palpable.

Stage IV: Full-thickness tissue loss with exposed bone, tendon, or muscle. Slough or eschar may be present on some parts of the wound bed. Often include undermining and tunneling.



The depth of a stage IV pressure ulcer varies by anatomical location. The bridge of the nose, ear, occiput, and malleolus do not have subcutaneous tissue and these ulcers can be shallow. Stage IV ulcers can extend into muscle and/or supporting structures (e.g., fascia, tendon, or joint capsule) making osteomyelitis possible. Exposed bone/tendon is visible or directly palpable.

Pressure Ulcer Stages as revised and copyrighted by Nat'l Pressure Ulcer Advisory Panel, Feb., 2007

Pressure Ulcer Stages (continued):



Unstageable Pressure Ulcer: Full-thickness tissue loss in which the base of the ulcer is covered by slough (yellow, tan, gray, green, or brown) and/or eschar (tan, brown, or black) in the wound bed.

Further Description: Until enough slough and/or eschar is removed to expose the base of the wound, the true depth, and therefore stage, cannot be determined. Stable (dry, adherent, intact without erythema or fluctuance) eschar on the heels serves as “the body’s natural (biological) cover” and should not be removed.

Sinus Tract: A course or pathway which can extend in any direction from the wound base; results in dead space with potential for abscess formation

Slough: Stringy, necrotic tissue; usually yellow, tan, gray, green, or brown

Strip: Removal of epidermis by mechanical means, usually tape

Ulcer: Loss of epidermis/dermis or mucous membrane with definite margins

Undermine: Skin edges of a wound that have lost supporting tissue under intact skin

Wound: A break in the integrity of the skin; an injury to the body which causes a disruption of the normal continuity of the body structures

Wound Margin: Rim or border of a wound

Skin Structure:

The **epidermis**, which is the outermost layer of the skin, is characterized as follows: avascular, varies in thickness (depending on body location), a dry structure which sheds cells and replaces itself every 4 to 6 weeks; approximately the thickness of a piece of plastic wrap

The **dermis** is located directly beneath the epidermis and is characterized as follows: provides strength and structural support through a vascular network of connective tissue blood vessels, nerves, hair, and nails. Sebaceous glands and sweat glands originate from this layer, which is thicker than the epidermis

Below the dermis is the **subcutaneous tissue** which is composed of major vessels, lymphatics, fat, and connective tissue. This area provides insulation and nutritional support for the skin. Located below the subcutaneous tissue are fascia, muscles, tendons, and bone. The thickness of the dermis and subcutaneous layers vary from person to person and on different parts of the body

Types of Ulcers

Pressure Ulcers: A pressure ulcer is localized injury to the skin and/or underlying tissue usually over a bony prominence as a result of pressure, or pressure in combination with shear and/or friction. A number of contributing or confounding factors are also associated with pressure ulcers; the significance of these factors is yet to be elucidated.


Venous Ulcers: These are the most common type of lower-leg ulcers in ambulatory people. The underlying cause involves vein damage or an incompetent calf muscle pump action which leads to venous hypertension. As a result, the blood pools in the lower extremities causing edema and leakage of fibrinogen and other blood products into the tissues. Trauma to the area or increased pressure within the tissues results in ulceration.

Arterial Ulcers: Arterial ulcers result from chronic or acute arterial insufficiency to the skin and subcutaneous tissue of the lower extremities. The most common cause is a progressive disease: atherosclerosis. The precipitating event leading to ulceration is usually trauma, such as a bumped toe or tight shoes. Arterial ulcers may occur alone or in combination with diabetes, venous stasis, and numerous other conditions. Multidisciplinary management of these patients with early intervention and close monitoring is key to prevention of more serious complications.

Neuropathic Ulcers: Neuropathic ulcers may occur in individuals with diabetes, spinal cord injury, Hansen's disease, or other conditions that result in loss of sensation in the legs and feet. Diabetic foot ulcers are most commonly caused by peripheral neuropathy and peripheral vascular disease. Multidisciplinary management of these patients with early intervention and close monitoring is key to prevention of more serious complications.

Educational Resources

Educational
Resources



ConnectEd

Wound Care Education
LEARNING MANAGEMENT SYSTEM

A smarter way to learn.

INTEGRATES INTO EXISTING LMS SYSTEMS

TRACKABILITY & REPORTING

24/7 WEB-BASED ACCESS

INDIVIDUALIZED LEARNING



The future of Wound Care Education is here.

ConnectEd is a new Web-based education resource from Hollister Wound Care that is transforming the way bedside nurses learn about advanced wound care.

Unlike other systems, **ConnectEd** can be integrated into existing learning management systems and provides customizable tracking and reporting documentation.

For more information visit www.HollisterWoundCare.com/Connect-Ed or call your Hollister Wound Care representative today.

Educational Materials

Literature No

Restore Wound and Skin Care Family Brochure	911056
The Use and Efficacy of Silver Dressings	910751
Contact Layer Brochure	921818
Wound Dressing Selection Guide/Poster	911292
Foam and Specialty Dressings Brochure	911038
Alginate Dressings Brochure	911037
Create the Shape Poster (foam dressing cutting guide)	911146
Duo/Trio Brochure	910725
Skin Tear Algorithm Pocket Card	921805
Skin Tears: State of the Science: Consensus Statements for the Prevention, Prediction, Assessment, and Treatment of Skin Tears	921822
The Role of a Silver Releasing Lipido-Colloid Contact Layer in Venous Leg Ulcers Presenting Inflammatory Signs Suggesting Heavy Bacterial Colonization: Results of a Randomized Controlled Study	910548
Use of a New, Flexible Lipidocolloid Dressing on Acute and Chronic Wounds: Results of a Clinical Study	921817
Restore Wound Dressing Inservice Videos DVD	911301
Hydrofera Blue Brochure	921637
Hydrofera Blue Application Guide	921639
Hydrofera Blue Clinical Monograph	921638
Hydrofera Blue Inservice Video DVD	921763

Measuring Tools

Literature No

Restore Wound Measuring Guide, Small	874366
Restore Wound Measuring Guide, Large	874367
Wound Measuring Guide, 6" ruler	908892
Pain Assessment Scale	910132

ConnectEd Educational Resources

Literature No

ConnectEd Brochure	921784
ConnectEd Course Selection Guide	921785
ConnectEd DVD	921786

Product Inservice Tools are available through your Hollister sales representative

Algorithms, Protocols, and Other Educational Tools are available on our website, www.hollisterwoundcare.com. For a complete description, visit the Educational Resources portion of our website, or contact your Hollister sales representative.

EB Educational Resources

Literature No

Using a Novel Contact Layer for the Management of Epidermolysis Bullosa Skin Lesions	907990
EB Dressing Guide	911071

*Contact our special EB Resource hotline at **1.888.EBAID4U** with questions or needs specific to epidermolysis bullosa (EB). Additional resources including clinical posters are available on our website, www.hollisterwoundcare.com. For a complete description, visit the EB Resources portion of our website.*

Clinical Poster Reprints

Literature No

Hydrofera Blue Foam Dressings

The Effect of Various Wound Dressings on the Activity of Debriding Enzymes 921636

A Comparison of Two Antimicrobial Dressings: A Randomized Prospective Trial Comparing PVA Foam with Two Organic Pigments to a Silver Based Wound Dressing 921644

Case Study Review: Use of an Absorbent Bacteriostatic Dressing for Multiple Indications 921789

Wound Care Challenges and Opportunities with Elderly Patients 921790

Restore Foam Dressings featuring TRIACT ADVANCED Technology

Leg Ulcers and a new Micro-Adherent Absorbent Dressing 911075

Pressure Ulcer Treatment: Changing Formulary, Changing Practice and Wound Care 911208

Restore Contact Layer FLEX with TRIACT Technology

Benefits of a Flexible Non-Adherent Lipido-Colloid Contact Layer Used in Conjunction with Negative Pressure Wound Therapy (NPWT) 921788

Restore Contact Layer with TRIACT Technology

What is the Cost of Pain? A Quantitative Model 910489

Evaluating a New Non-Adherent Contact Layer in the Treatment of Skin Tears 910483

Simple Changes in Wound Trauma Management Yield Significant Benefits for Everyone 910487

Using a Novel Contact Layer for the Management of Epidermolysis Bullosa Skin Lesions 907990

Assessment of a New Contact Layer in the Local Management of Pediatric Wounds: Results of a European Clinical Study 907984

Use of Negative Pressure Wound Therapy with a New Wound Contact Layer: Results of a Clinical Evaluation 907982

Use of a Novel Contact Layer in Wound Care 907978

Clinical Evaluation of a New Contact Layer Dressing for the Local Management of Acute and Chronic Wounds 909916

The Role of "Touch" in the Assessment of New Wound Care Dressings 909918

Clinical Evaluation of Non-Adherent Contact Layer and Foam Dressings with and without Silver for the Local Management of Acute and Chronic Wounds 910644

Restore Contact Layer with TRIACT Technology, Silver

A Cost Comparison of a Silver, Pain-free Donor Site Dressing vs. a Traditional Donor Site Dressing 910472

Maximizing Reimbursement for Wound Care is Crucial to the Survival of the Agency: One Agency's Success Story 910485

Use of a Contact Layer with Silver Under a Negative Pressure Dressing 910278

Minimizing Complications During Dressing Changes of a Bilayer Skin Equivalent in the Treatment of Venous Leg Ulcers 910469

Clinical Poster Reprints (continued)

Literature No

Clinical Evaluation of Non-Adherent Contact Layer and Foam Dressings with and without Silver for the Local Management of Acute and Chronic Wounds 910644

Hospital-Acquired Surgical Site Infections: Considering a Simple, Effective, Versatile Silver Technology to Improve Patient Care 910654

Restore Foam Dressings with TRIACT Technology

Leg Ulcers and a New Foam Dressing: Results of Two French Multicenter Clinical Studies 907986

Use of a Novel Foam Dressing with Non-Adherent Contact Layer in Leg Ulcer Management 907980

Clinical Evaluation of a New Foam Dressing with and without Adhesive for the Local Management of Acute and Chronic Wounds 909915

Clinical Trial of a New Foam Dressing for the Local Management of 10 Donor Sites 909917

Clinical Evaluation of Non-Adherent Foam Dressings with and without Silver for the Local Management of Split-Thickness Donor Sites 910689

Clinical Evaluation of Non-Adherent Contact Layer and Foam Dressings with and without Silver for the Local Management of Acute and Chronic Wounds 910644

Restore Foam Dressings Silver, Non-Adhesive with TRIACT Technology

Critical Colonization and a New Absorbent Dressing Impregnated with Silver Salts: Results of a Clinical Study in the Management of Leg Ulcers 907992

Application of Silver Foam Dressings to Donor Sites in Immunocompromised Patients Can Enhance Wound Healing 910275

Use of a Non-Adhesive Silver Absorbent Foam Dressing Underneath an Unna Boot 910317

Clinical Evaluation of Non-Adherent Foam Dressings with and without Silver for the Local Management of Split-Thickness Donor Sites 910689

Clinical Evaluation of Non-Adherent Contact Layer and Foam Dressings with and without Silver for the Local Management of Acute and Chronic Wounds 910644

Restore Calcium Alginate with Silver

Visualizing the Death of Wound Microorganisms Within a Silver Alginate Dressing 921762

Easy Does It: Evaluation of a Silver Alginate vs. Silver Fiber Dressing in Cost Effectiveness and Ease of Use 910655

Restore Hydro-Shield Foam Dressings

The Role of Ready-To-Use, Absorbent Waterproof Dressings as Primary Wound Care Dressing 921663

Restore DUO & TRIO Absorbent Dressings with TRIACT Technology

Emergency Department Treatment of Low-Exuding Traumatic Wounds, sutured or otherwise, with a Ready-to-Apply Lipido-Colloid Dressing 911074

Clinical Poster Reprints (continued)

Literature No

Restore DUO & TRIO Absorbent Dressings with TRIACT Technology (continued)

The Role of Ready-to-Use, Lightly Absorbent Dressings with Non-Adherent Contact Layer as Primary Wound Care Dressing	911078
Clinical Evaluation of a New Absorbent, Non-Adhering Dressing with and without Adhesive for Local Management of Acute and Chronic Wounds	910982
Practice Innovation: Skin Tear Treatment that Changes Outcomes	921664

ConnectEd Learning Management System

A Novel Approach to Improving Wound Outcomes Through an Evidence-Based Learning Management System	911077
---	--------

Other Topics

WOCN Workgroup: Engaging the Power of Standardization in Wound Care	911207
---	--------

Specialty Organizations

DeBRA of America, Inc.

www.debra.org

1.866.332.7276

The Association for the Advancement of Wound Care (AAWC)

www.aawconline.org

1.866.229.9999

Wound, Ostomy and Continence Nurses Society (WOCN)

www.wocn.org

1.888.224.9626

American Academy of Wound Management (AAWM)

www.aawm.org

1.202.521.0368

American Burn Association

www.ameriburn.org

1.312.642.9260

Wound Healing Society

www.woundheal.org

1.763.765.2377

National Pressure Ulcer Advisory Panel

www.npuap.org

1.202.521.6789

American Diabetes Association

www.diabetes.org

1.800.DIABETES (342.2383)

National Spinal Cord Injury Association

www.spinalcord.org

1.800.962.9629

U.S. Department of Health and Human Services

Centers for Medicare & Medicaid Services

www.cms.hhs.gov

National Guideline Clearinghouse

www.guideline.gov

European Wound Management Association

www.ewma.org

+45 70 20 03 05

Index • Stock Number

Stock No	Page(s)	Stock No	Page(s)	Stock No	Page(s)
10xx	92	1915x	35	509349	131
1410x	27	1920x	34	509372	132
1420x	25	1925x	35	509373	132
143x	78	31003	109	509374	132
1430x	26	31005	109	509375	131
1440x	28	31008	109	509376	131
1450x	29	31009	109	509377	131
1460x	25	314x	79	509378	130
1470x	26	318x	79	509379	130
1480x	28	32004	108	509380	133
149xx	29	32005	108	509381	133
1520x	25	32006	108	509382	133
1540x	28	32007	108	509384	133
1550x	29	322x	77	509387	133
1560x	25	323x	77	509388	133
1580x	28	327x	78	509389	133
1590x	29	33004	108	509390	133
1610x	26	332x	78	509393	133
1640x	26	334x	75	509394	133
1796	51	340x	76	509395	133
1800x	30	360x	77	509396	133
1801x	32	366x	77	517210	69, 138
1810x	31	370x	73	517220	70, 138
1811x	30	372x	72	519904	136
1812x	31	373x	74	519906	136
1813x	30	374x	73	519921	136
1814x	31	376x	72	519923	136
1815x	33	3761	60	519925	136
1816x	31	3774	61	519930	136
1817x	31	3777	59	519932	136
1818x	30	3778	59	519935	136
1819x	30	3795	59	519953	136
1820x	31	3796	59	519954	136
1828x	30	3797	60	519955	136
18282	60	3799	60	519956	136
1829x	30	380x	74	519957	136
18292	60	381x	74	519958	136
1832x	32	382x	76	519959	136
1833x	33	383x	74	519963	136
1835x	33	384x	75	519965	136
1836x	32	387x	75	529300	134
1837x	32	389x	74	529301	134
1838x	32	390x	76	529302	134
18382	61	506487	130	529303	134
1839x	32	506488	130	529304	134
18392	61	506489	130	529305	134
1840x	33	509338	130	529306	134
18402	61	509339	130	529937	135
1841x	33	509341	130	529938	135
18412	61	509342	130	529939	135
1842x	33	509343	131	529940	135
1873x	33	509344	131	529967	135
1875x	33	509345	131	529968	135
1900x	34	509346	131	529969	135
1905x	35	509347	131	529972	135
1910x	34	509348	131	529973	135

<u>Stock No</u>	<u>Page(s)</u>	<u>Stock No</u>	<u>Page(s)</u>	<u>Stock No</u>	<u>Page(s)</u>
529974	135	834x	75	9343	97
529975	137	843x	52	9345	97
529976	137	844x	51	9346	97
529977	138	8440	51	9348	96
529978	138	8450	51	9349	96
529979	138	8460	51	940xx	94
716x	78	846x	51	941xx	94
7299	68	847x	52	951xx	94
7300	68	848x	52	960xx	94
7331	69	849x	52	961xx	94
7700	71	850x	46	9645	96
771x	69	851x	47	9655	96
7721	70	852x	46	970x	49
7723	70	853x	44	971xx	94
7724	70	8531	41	974xx	95
7728	70	855x	44	975xx	95
773x	69	8551	41	976xx	95
7730	68	858x	46	977x	49
7731	68	859x	47	977xx	95
7736	69	861x	47	9779	49
7740	70	862x	47	978xx	95
7760	68	8631	41	9780	116
7766	71	863x	44	9781	115
7767	70	8641	41	9782	115
780x	67	864x	45	9785	116
7800	71	8651	41	9786	116
78500	67	872x	72	9788	115
78501	67	873x	73	979xx	95
7905	71	874x	73	9799	115
7906	67	875x	72	9805	96
7910	71	877x	73	9811	98
7917	68	8770	69	9814	96
79300	67	8800	68	9821	110
79301	67	8801	68	9822	110
79400	67	881x	74	9823	110
795xx	67	883xx	45	9839	71, 98
796xx	67	88300	42	9840	98
80xx	45	88302	43	9873	98
8031	42	884xx	45	9880	110
8081	42	88400	42	Bxxx	93
81xx	45	88401	43	BxxxB	93
8131	42	88402	43	HB2214	129
8181	42	88501	44	HB4414	129
80110	48	887xx	46	HB6614	129
80070	48	88700	43	HBT0906	129
81070	48	88800	43	HBHD4450	129
81110	48	890x	75	HBHD6675	129
821xx	51	89001	52	HBRF4414	129
82100	50	89002	53	HBRF2258	129
823xx	50	89003	53	HBRF4475	129
82300	50	89004	53	HBRF2650	129
82302	50	89005	53		
82400	50	910x	95		
82402	50	920x	95		
82500	50	931xx	94		
8331	41	9342	97		

Name	Page(s)
A ctiFlo Accessories.....	109
ActiFlo Catheter Kits with Disposable Odor-Barrier Collection Bags.....	108
ActiFlo Catheter Kits with Drainable Odor-Barrier Collection Bag	108
Adapt Odor Eliminator and Lubricant.....	67
Adapt Ostomy Belts.....	68
Adapt Skin Barrier Accessories	67
Adapt Skin Barrier Paste.....	67
Adapt Stoma Powder.....	67
Adhesives and Removers.....	68
Advance Plus Touch Free Intermittent Catheter	94
Advance Plus Touch Free Intermittent Catheter System	94
Anchor Fast Oral Endotracheal Tube Fastener	115
Apogee Plus System.....	93
Apogee Plus System Kit.....	93
Apogee Intermittent Catheter	92
B edside Drainage Collection System with Anti-Reflux Valve.....	71, 98
C enterPointLock Closed Pouches	75
CenterPointLock Drainable Pouches.....	74–75
CenterPointLock Flat Skin Barriers with Floating Flanges.....	72–73
CenterPointLock Stoma Cap	76
CenterPointLock Urostomy Pouches.....	75–76
Convex Skin Barriers with Floating Flanges.....	73–74
D rainable Pouch Clamp	69
Drain Tube Attachment Device.....	116
E veryday Male External Catheter	95
Extended Wear Male External Catheter.....	95
F abric Leg Bag Straps	97
Fecal Collectors.....	110
Feeding Tube Attachment Device.....	116
Female Urinary Pouch External Collection Device.....	98
H ollister Universal Catheter Access Port	49
Horizontal Tube Attachment Device.....	115
Hydrofera Blue Foam Dressings	129
I nstaFlo Catheter Kit.....	108
InView Male External Catheters	95
Irrigation Products	70
Irrigator Sleeve	76
K araya 5 Closed Pouches	78
Karaya 5 Drainable Pouches.....	77–78
Karaya 5 Urostomy Pouches	78
L atex-Free Extension Tubing.....	97
Loop Ostomy Bridge.....	70
m 9 Cleaner/Decrystallizer	69
m9 Odor Eliminators.....	68
Male Urinary Pouch External Collection Device.....	98

Name	Page(s)
N ew Image Closed Pouches.....	32–33
New Image Convex Skin Barriers with Integrated Floating Flanges	28–29
New Image Drainable Pouches	30–31
New Image Flat Skin Barriers with Integrated Floating Flanges.....	25–27
New Image High Output Drainable Pouches	32
New Image Irrigator Sleeve	33
New Image Non Sterile Kits — Drainable.....	34
New Image Non Sterile Kits — Urostomy.....	34
New Image Sterile Kits — Drainable.....	35
New Image Sterile Kits — Urostomy.....	35
New Image Urostomy Pouches.....	33
O ther Closed Pouches	79
Other Stoma Cap.....	79
P ouchkins Newborn Pouch.....	59
Pouchkins One-Piece Pouches.....	59–60
Pouchkins Pediatric Ostomy Belt.....	61
Pouchkins Premie Pouch	59
Pouchkins Two-Piece Pouches	60–61
Pouchkins Two-Piece Skin Barrier.....	60
Premier Closed Pouches.....	50–51
Premier Drainable Pouches with Convex Barriers.....	46–47
Premier Drainable Pouches with Flat Barriers	41–46
Premier High Output Drainable Pouches.....	48
Premier Non Sterile Kits — Drainable	52–53
Premier Non Sterile Kits — Urostomy	53
Premier Stoma Cap	51
Premier Urostomy Pouches with Convex Barriers.....	52
Premier Urostomy Pouches with Flat Barriers	51
Premier Wound Drainage Collectors.....	49
R eplacement Filters.....	71
Restore Adhesive Foam Dressing with TRIACT Technology	131
Restore Calcium Alginate Dressing Silver	135
Restore Calcium Alginate Dressing.....	135
Restore Cleanser & Moisturizer	138
Restore Contact Layer Dressing with Silver and TRIACT Technology.....	130
Restore Contact Layer Dressing with TRIACT Technology.....	130
Restore Contact Layer FLEX Dressing with TRIACT Technology	130
Restore DimethiCreme Skin Protectant.....	138
Restore DUO Absorbent Dressing with TRIACT Technology	131
Restore Extra Thin Hydrocolloid Dressing.....	136
Restore Foam Border Dressing Featuring TRIACT ADVANCED Technology.....	133
Restore Foam Dressing without Border	133
Featuring TRIACT ADVANCED Technology	
Restore Heel Dressing without Border	133
Restore Hydrocolloid Dressing.....	136
Restore Hydrocolloid Dressing with Foam Backing.....	136
Restore Hydrogel Dressing.....	135

Name	Page(s)
Restore Hydro-Shield Foam Dressing	134
Restore Hydro-Shield Foam Dressing without Border	134
Restore LITE Foam Dressing without Border	133
Featuring TRIACT ADVANCED Technology	
Restore Non-Adhesive Foam Dressing with Silver and	131
TRIACT Technology	
Restore Non-Adhesive Foam Dressing with TRIACT Technology	131
Restore Odor-Absorbent Dressing	136
Restore Skin Care Products	69
Restore Skin Cleanser	69, 138
Restore Skin Conditioning Creme	138
Restore TRIO Absorbent Dressing, Adhesive, with TRIACT Technology	132
Restore Wound Cleanser	137
S kin Barriers (Other)	68, 71
Skin Barrier Pastes and Powders (Other)	71
Skin Prepping Product	68
Suction Tube Attachment Device.....	116
U rinary Leg Bag.....	96
Urinary Leg Bag Combination Pack	96
Urostomy Drain Tube Adapter	69
V ented Urinary Leg Bag Combination Pack.....	96
Vertical Tube Attachment Device.....	115

Index • HCPCS Codes

<u>HCPCS</u>	<u>Page(s)</u>	<u>HCPCS</u>	<u>Page(s)</u>
A4326	98	A5055	51, 76, 79
A4328	98	A5056	41
A4330	110	A5057	46
A4331	97	A5061	42, 43, 44, 52, 53, 59
A4349	95	A5062	59
A4351	92	A5063	30, 31, 34, 35, 60, 74, 75
A4352	92	A5071	60
A4353	93, 94	A5114	97
A4357	71, 98	A5120	68
A4358	96	A5122	68
A4362	59, 71	A5131	69
A4363	69	A5200	115, 116
A4364	68	A6154	48, 49
A4367	61, 68	A6196	135
A4368	71, 74, 77	A6197	135
A4371	67, 71	A6199	135
A4373	28, 29, 74	A6206	130
A4385	67, 68	A6207	130
A4387	78	A6208	130
A4388	41, 44, 45, 53	A6209	129, 131, 133, 134, 136
A4389	77, 78	A6210	129, 131, 133, 134
A4390	46, 47	A6211	134, 136
A4394	67, 69	A6212	129, 131, 133, 134
A4397	33, 70, 76	A6213	131, 133, 134
A4398	70	A6215	129
A4399	70	A6231	135
A4402	70	A6234	131, 136
A4405	71	A6235	131, 136
A4406	67	A6236	136
A4407	28, 29, 73	A6237	132
A4408	28	A6238	132
A4409	25, 26, 27, 34, 35, 72, 73	A6248	135
A4410	25, 72	A6250	70, 138
A4411	67	A6260	137
A4412	32, 75	A6266	135
A4414	25, 26, 34, 35, 60, 72, 73	A9270	69, 96, 138
A4415	25, 34, 35, 72, 73	A9900	115
A4416	50, 51	A9999	115, 116
A4417	78	B9998	116
A4418	79		
A4419	32, 61, 75		
A4421	69, 70		
A4424	42, 43, 45, 46, 59		
A4425	30, 31, 60, 74		
A4428	51, 53		
A4429	78		
A4430	52		
A4432	33, 34, 35, 61, 75, 76		
A4455	68		
A4456	68		
A5052	59		
A5054	33		

Hollister and logo, Hollister, ActiFlo, Adapt, Advance Plus, AF300, Anchor Fast, Apogee, Apogee Plus, CenterPointLock, ComfortWear, Extended Wear, Everyday, Flexextend, FlexWear, FormaFlex, HollisHesive, InstaFlo, InView, Karaya 5, Lock 'n Roll, m9, New Image, Pouchkins, Premier, Premium, QuietWear, and SoftFlex are trademarks of Hollister Incorporated.

C3Life, Secure Start, and "Hollister Ostomy. Details Matter." are service marks of Hollister Incorporated.

Hollisterwoundcare and wave logo are trademarks of Hollister Incorporated. Restore, Hydro-Shield and graphic, Restore DUO, Restore TRIO, TRIACT, TRIACT ADVANCED and graphic are trademarks of Hollister Wound Care LLC.

ConnectEd is a service mark of Hollister Wound Care LLC.

Hydrofera Blue is a trademark of Hydrofera, LLC.

Shadow Buddies is a trademark of Shadow Buddies.

Teflon is a trademark of E.I. du Pont de Nemours and Company.

Coloplast and SenSura are trademarks of Coloplast A/S.

ConvaTec and SUR-FIT Natura are trademarks of ConvaTec, LLC.

Zassi is a trademark of Zassi Medical Evolutions, Inc.

© 2012 Hollister Incorporated. Printed in USA.

Hollister Consumer Programs

For consumer and healthcare professional product questions and sampling needs:

1.888.740.8999

For dealers, distributors, institutions, and pharmacy questions or orders:

1.800.323.4060

fax 1.847.680.1017

For helping new patients adjust to life with an ostomy or bladder management:

Secure Start Services

1.888.808.7456

For additional information on Hollister Wound Care:

www.hollisterwoundcare.com

For helping patients and families dealing with Epidermolysis Bullosa (EB) contact our EB Resources Hotline:

1.888.EB AID4U

(1.888.322.4348)



Hollister Incorporated
2000 Hollister Drive
Libertyville, Illinois 60048 USA
www.hollister.com