



AnchorFast Guard Oral Endotracheal Tube Fastener

Introduction

This is a skills checklist for the application and routine care of the AnchorFast Guard oral endotracheal tube fastener. The AnchorFast Guard tube fastener secures oral endotracheal tubes ranging in size from 5 to 10 mm diameter. The suitability of the oral endotracheal tube fastener must be assessed for each patient. Refer to entire Instructions for Use for detailed information.

Date Completed	Name	Employee Number	Dept.
Criteria (<u>Refer to Instructions for Use for detailed information</u>)			
Prepare the Skin and Apply the Fastener			
1. Skin should be clean, dry, and free of oily residue			
2. Do not use Skin Gel Wipes or other brands of skin preps with this device			
3. Remove the release liners from the two skin barriers			
4. Ensure the tube is positioned adjacent to the opening of the tube protection sleeve			
5. Center the device on the upper lip, so the lip stabilizer touches the skin (do not apply pressure). Press barriers on the skin and hold in place for 30 seconds.			
Apply the Neck Band and Secure the Tube			
1. Secure the neck band using the plastic loop closures. Do not over tighten. Allow two fingers width between the neck band and the back of the patient's head.			
2. Squeeze the tabs on the sides of the shuttle and move the clamp until it is over the ET tube			
3. Carefully slide the tube into the tube protection sleeve and place the inflation lumen within the channel of the tube protection sleeve			
4. Make sure the ET tube is dry, and remove the release liner from the strap and wrap it tightly around the tube. Secure the strap by snapping the clamp shut (an audible click will be heard)			

(Continued)

AnchorFast Guard Oral Endotracheal Tube Fastener

Routine Care		
1. Reposition the tube by squeezing the tabs on the outer edges and moving in either direction along the track. Reposition every two hours ^{1,2} or more often as patient condition dictates, to minimize the risk of injury to the skin and/or lips from unrelieved pressure		
2. Inspect the patient's lips and skin every two hours or more frequently, and discontinue use of the device if redness or skin irritation occurs		
3. Check the patient frequently to ensure that both the AnchorFast Guard device and the ET tube are secure and correctly positioned. REMINDER: Take precaution that repeated adjustment of the ET tube in a distal or proximal direction may affect the performance of the adhesive strap.		
Removal		
1. Release the clamp holding the flexible strap in place		
2. Unwrap the adhesive strap from around the tube then remove the tube and inflation lumen from the tube protection sleeve		
3. Release the neck band by unfastening the hook and loop closures		
4. Remove the skin barriers by gently peeling them away from the patient's skin		

* Evaluation Method: (L) return demos live in class or skill lab, (O) observation in worksetting, (CR) chart review, (V) verbalizes information, (R) read only – not performed – for information, (NA) not applicable.

Reviewed by (signature) _____.

References

1. Pieper, Barbara. "Mechanical Forces: Pressure, Shear, and Friction." Chapter 12 in *Acute & Chronic Wounds: Nursing Management*, 3rd ed., ed. Ruth A. Bryant & Denise P. Nix. (St. Louis: Mosby, Inc., 2007).
2. Panel for the Prediction and Prevention of Pressure Ulcers in Adults. "Pressure Ulcers in Adults: Prediction and Prevention." *Clinical Practice Guidelines*, Number 3. AHCPR Publication No. 92-0047. (Rockville, MD: Agency for Health Care Policy and Research, Public Health Service, U.S. Department of Health and Human Services. May 1992.)

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician. Prior to use of the AnchorFast Guard Oral Endotracheal Tube Fastener, be sure to read the entire product instructions for use packet insert that accompanies the product.

The Hollister logo and AnchorFast Guard are trademarks of Hollister Incorporated. © 2017 Hollister Incorporated.



Hollister Incorporated
 2000 Hollister Drive
 Libertyville, Illinois 60048
 1.800.323.4060

www.hollister.com

Hollister Limited
 95 Mary Street
 Aurora, Ontario L4G 1G3
 1.800.263.7400