

Using Convexity

Definition

Convexity is the outward curving of the skin barrier and is designed to interface with the immediate peristomal skin. Convex skin barriers promote a good fit between the ostomy pouching system and the peristomal skin where flat barriers would be unsuccessful. This is achieved as the convexity barrier pushes on the surrounding skin, opens or flattens skin folds or helps the stoma protrude more.



Stoma that is below skin level



Convex skin barrier in place

General Information

- A secure and reliable seal between the peristomal skin and ostomy skin barrier is necessary to protect healthy peristomal skin from stoma discharge
- Skin barriers are either convex or flat. Convex simply describes the shape of the barrier
- Convexity may be recommended to increase the wear time of the ostomy pouching system

Assessing for Convexity


- A simple patient assessment methodology can assist in determining the need for convexity. Refer to the Hollister Fit Indicator Tool
- Assessing when convexity is required should be done in consultation with an ostomy care nurse
- Ideally the patient should be assessed in the sitting position
- The pouch should be removed prior to assessment
- A distal stoma assessment is recommended for all loop stomas

Convexity May be Considered When:

- Stoma Assessment
 - Liquid/Loose stoma output
 - The opening of the stoma is at or below the skin surface
 - The opening of the stoma is off-center
 - The stoma protrudes <math>< 3/4''</math> (20mm) or is retracted below the skin surface
 - The stoma becomes flush/retracted during peristalsis (telescopes)
- Peristomal Assessment
 - The peristomal topography is uneven or irregular
 - The abdominal tone is flaccid



Off Centered Stoma Opening



Retracted Stoma



Flush Stoma



Uneven/Irregular



Flaccid Abdomen

Using Convexity

Convex Products

- Convexity can either be built into the pouching system (integrated) or added with an accessory to the skin barrier
- Convexity is available in one-piece and two-piece options
- Convexity is available in pre-sized or cut-to-fit options
 - Note: When using cut-to-fit options for oval or irregular shape stomas, choose a barrier that is closest to the largest measurement so the convex curvature of the barrier is as close to the stoma as possible
- **Soft** convexity can be integrated into the pouching system or created by adding a convex barrier ring to a flat one-piece pouching system (Refer to Ostomy Care Tip on Adapt convex barrier rings)
- **Firm** convexity can be integrated into the pouching system or created by adding a convex barrier ring to a convex barrier to deepen the curvature when necessary




Soft integrated convex barrier



Soft convexity created by adding a convex barrier ring to flat skin barrier

Convexity Considerations

- Selecting the right convex product may require regular reassessment of products used, condition of the peristomal skin and the pouching system wear time
- Some ostomy care nurses recommend an ostomy belt to enhance the effect of the convexity. (Refer to the Ostomy Care Tip for Adapt ostomy belts)
- The following peristomal conditions do not preclude the use of convex products, however, caution and frequent reassessment is recommended in patients with:
 - Pyoderma Gangrenosum
 - Crohn's ulcers
 - Caput medusa (peristomal varices)
 - Parastomal hernia
 - Pressure injury
- Patients should be encouraged to contact an ostomy care nurse if there are stomal/peristomal changes, pouching system leakage or weight gain/loss



Firm integrated convex barrier



Firm integrated convex barrier with convex barrier ring to deepen the curvature

References:

- "Patient Assessment Guidelines for Convexity" – Hollister Incorporated © 2013
- "Developing Patient Assessment Guidelines for Convexity" – Drolshagen et al. 2014 Poster Presentation, WOCN Conference, Nashville TN USA
- Best Practice Document: Convex Products – WOCN 2012
- "Convex Pouching Systems: Best Practice Document for Patients - WOCN 2012
- Use of Convexity in Pouching – A Comprehensive Review – Hoefflok J, Kittscha J, Purnell P. JWOCN 2012:40 (5):506-512
- "Adapt Convex Barrier Rings" – Ostomy Care Tip – Hollister Incorporated © 2012
- Lyon, C. & Smith, A. 2010 – Abdominal Stomas and Their Skin Disorders

For product questions, sampling needs, or detailed clinical questions concerning our products in the US, call **1.888.740.8999**. In Canada call **1.800.263.7400**.

Routine follow-up with your healthcare professional is recommended.

Prior to using any ostomy products/accessories be sure to read all product inserts and labels for any indications, contraindications, hazards, warnings, cautions and directions for use.



Hollister Incorporated
2000 Hollister Drive
Libertyville, Illinois 60048
1.800.323.4060

www.hollister.com

Hollister Limited
95 Mary Street
Aurora, Ontario L4G 1G3
1.800.263.7400