

Adapt Medical Adhesive Remover Spray/Adapt Universal Remover Wipe

General Information

- Helps remove adhesive skin barriers and/or tape
- Removes adhesive residue from the skin
- Skin Friendly Formulation
 - Silicone-based
 - Alcohol free
 - No sting



Adapt Medical Adhesive Remover Spray/Adapt Universal Remover Wipe
Stock No 7737/7760

Using Adhesive Remover Spray

1. Peel back one edge of the skin barrier or tape border.
2. Hold the can so that the arrow faces the base of the skin barrier or tape, and spray until the adhesive releases. The 360° spray technology minimizes waste, and will allow you to spray from any angle.
3. Ensure your skin is dry and apply your new ostomy pouching system as usual.



Spray - Spray until the skin barrier releases.

Using Universal Remover Wipes

1. Tear open the package and take out the wipe.
2. Press or rub the wipe against your skin as you peel off your skin barrier and/or tape.
3. Press or rub the wipe over any areas of the skin where tape residue remains.
4. Ensure your skin is dry and apply your new ostomy pouching system as usual.



Wipes - Rub the wipe over skin as you peel off skin barrier.

FREQUENTLY ASKED QUESTIONS	
Is an adhesive remover always needed?	No, this product is optional. Use it to remove adhesive residue from your skin, or to make it easier to remove your skin barrier or tape.
Will this be uncomfortable to use?	This silicone-based product is formulated to be a no-sting alternative. When used correctly, the adhesive remover will help to make the removal of your pouching system easier and help reduce the potential for skin stripping.
Do I need to wash the adhesive remover off of my skin?	The adhesive remover evaporates quickly and dries without residue, so the next barrier can attach securely. If you wash your skin after use be sure to rinse and dry thoroughly before the next pouch application.



Wipes - Clean off remaining tape residue.

