

Extended Wear Drainage Bag

Key Features

- 2 litre capacity
- sterile fluid pathway
- integral bed hanger and sheet clip for added security
- easy action drain tube clamp
- soft, non-kink inlet line with self sealing sample port
- non-latex

Ext. Wear Drainage Bag	Code No.	Capacity	Tubing	Box Quantity
	9839	2000ml	120cm	1

Extended Wear Leg Bag Kits

Key Features

- discreet 540ml or 900ml capacity
- sterile fluid pathway
- non-latex tough polyvinyl construction
- supplied with adjustable extension tube
- secure, positive action outlet valve
- Available as complete kits (bags, straps and tubing) or as individual items

Ext Wear Leg Bags and Kits	Code No.	Description
	9814	540ml - Leg Bag
	9805	900ml - Leg Bag
	9342	38cm - Straps
	9343	58cm - Straps
	G9000	Silicone Tubing - 42cm
	G9006	540ml Leg Bag Kit
	G9007	900ml Leg Bag Kit

Leg Bag Kits Contain

- 1 x leg bag (540ml or 900ml)
- 1 x pair leg bag straps
- 1 x sterile extension tubing
(for connection of leg bag to sheath)



Australia

6/345 Ingles Street
Port Melbourne, Victoria 3207

FreeCall: 1800 335 911

New Zealand

58 Richard Pearse Drive
Airport Oaks, Auckland

FreeCall 0800 167 866

www.ostomyworld.com.au

©2004 Hollister Incorporated. Hollister and logo, Flexend, InView are registered trademarks of Hollister Incorporated.

HICRFPB0DLV12-04547

Continence Products

Cost Effective Solutions

Helping you maintain dignity





Urinary incontinence is a major quality of life problem affecting many people.

Those who experience this problem may avoid social activities due to the feelings of shame, embarrassment and fear of exposure.

Hollister is leading the way with successful treatment and management options to decrease or correct incontinence.

InView sheaths are a new range of all silicone products, which compliment the growing Hollister latex-free family. This formulation is non-irritating and is perfect for people with known latex sensitivity.


InView Sheaths

The secure and sensitive solution

For incontinent men suffering from urge or reflex emptying of the bladder, the InView Sheath offers unrivalled security, comfort and flexibility

Key Features

- **Two lengths to choose from**
 InView Standard - unrivalled ease of use, reliability and comfort
 InView Special - its shorter length and strategically placed adhesive provide a better fit for patients with a shorter or retracted penis who find standard sheath lengths difficult to manage - supplied with plastic disposable cap to maintain adhesive integrity
- **Skin-friendly self-adhesive band**
 adheres gently and securely to the skin
 an easy method of attachment
- **Silicone, latex-free material**
 adapts to fluctuating penis size
 ideal for caregivers and patients with latex sensitivity
 breathable silicone is more comfortable to wear than sheaths made from some other materials
- **Clear material**
 a discreet solution for the patient
 permits routine inspection of the skin without removing the sheath

InView Standard	Colour Code	Size	Code No.	Box Quantity
	●	25mm	97525	30
	●	29mm	97529	30
	●	32mm	97532	30
	●	36mm	97536	30
	●	41mm	97541	30

InView Special	Colour Code	Size	Code No.	Box Quantity
	●	25mm	97425	30
	●	29mm	97429	30
	●	32mm	97432	30
	●	36mm	97436	30
	●	41mm	97441	30

Leg Bag (Sterile) Drainage Bags

Leg Bags

Hollister Leg Bags offer flexibility to suit a wide range of patient needs

This discreet collection system complements an active lifestyle and blends user-friendly features with the optimum standards of hygiene


- 500ml and 800ml capacities
- tubing lengths of 10cm or 50cm (cut to fit)
- waterproof comfort backing material that inhibits unwanted odour
- anti-reflux valve protects against leakage and urinary tract infections caused by urine backflow
- outlet valve features easy to handle lever tap
- protected sampling port to minimise needle stick injuries


Leg Bag Straps

- comfortable and cut to fit
- anti-slip silicone bands
- attachment to thigh or calf
- code no. 9625 - pack size 5 pairs


Drainage Bags (Sterile)

- 120cm inlet tubing - allows freedom of movement
- 2000ml capacity
- anti-reflux valve protects against leakage and infection caused by urine backflow
- easy to handle lever tap and sampling port

Leg Bag Kit (Sterile)	Code No.	Capacity	Tubing	Box Quantity
	G9621K	500ml	10cm	10
	G9624K	500ml	50cm	10
	G9631K	800ml	10cm	10
	G9632K	800ml	50cm	10

Drainage Bag	Code No.	Capacity	Tubing	Box Quantity
	9650	2000ml	120cm	10

More choice,
more comfort



Easy-to-handle twist-action lever tap