

Hollister

Continence Management Systems


Extended Wear Drainage Bags



Extended Wear Drainage Bag

Key Features

- 2 litre capacity
- sterile fluid pathway
- integral bed hanger and sheet clip for added security
- easy action drain tube clamp
- soft, non-kink inlet line with self sealing sample port
- non-latex

Ext. Wear Drainage Bag	Code No.	Capacity	Tubing	Box Quantity
	9839	2000ml	120cm	1

Extended Wear Leg Bag Kits

Key Features

- discreet 540ml or 900ml capacity
- sterile fluid pathway
- non-latex tough polyvinyl construction
- supplied with adjustable extension tube
- secure, positive action outlet valve
- Available as complete kits (bags, straps and tubing) or as individual items

Ext Wear Leg Bags and Kits	Code No.	Description	Code No.	Description
	9814	540ml - Leg Bag	G9006	540ml Leg Bag Kit
	9805	900ml - Leg Bag	G9007	900ml Leg Bag Kit
	9342	38cm - Straps	Leg Bag Kits Contain 1 x leg bag (540ml or 900ml) 1 x pair leg bag straps 1 x sterile extension tubing (for connection of leg bag to sheath)	
	9343	58cm - Straps		
	G9000	Silicone Tubing - 42cm		



Hollister Australia
6/345 Ingles Street
Port Melbourne, Victoria
Australia 3207

Hollister New Zealand
58 Richard Pearse Drive
Airport Oaks
Auckland, New Zealand

For more information, Freecall
Australia 1800 335 911
New Zealand 0800 167 866