

Hollister

Continence Management Systems

Leg and Drainage Bags



Leg Bag (Sterile) Drainage Bags

More choice,
more comfort



Leg Bags

Hollister Leg Bags offer flexibility to suit a wide range of patient needs

This discreet collection system complements an active lifestyle and blends user-friendly features with the optimum standards of hygiene


- 500ml and 800ml capacities
- tubing lengths of 10cm or 50cm (cut to fit)
- waterproof comfort backing material that inhibits unwanted odour
- anti-reflux valve protects against leakage and urinary tract infections caused by urine backflow
- outlet valve features easy to handle lever tap
- protected sampling port to minimise needle stick injuries

Leg Bag Straps

- comfortable and cut to fit
- anti-slip silicone bands
- attachment to thigh or calf
- code no. 9625 - pack size 5 pairs

Drainage Bags (Sterile)

- 120cm inlet tubing - allows freedom of movement
- 2000ml capacity
- anti-reflux valve protects against leakage and infection caused by urine backflow
- easy to handle lever tap and sampling port




Easy-to-handle twist-action lever tap




'Christmas tree' adapter is designed to fit smoothly, while providing maximum security

Leg Bag Kit (Sterile)

	Code No.	Capacity	Tubing	Box Quantity
Includes 1 pr Leg Bag Straps 	G9621K	500ml	10cm	10
	G9624K	500ml	50cm	10
	G9631K	800ml	10cm	10
	G9632K	800ml	50cm	10

Drainage Bag

	Code No.	Capacity	Tubing	Box Quantity
	9650	2000ml	120cm	10



Hollister Australia
6/345 Ingles Street
Port Melbourne, Victoria
Australia 3207

Hollister New Zealand
58 Richard Pearse Drive
Airport Oaks
Auckland, New Zealand

For more information, Freecall
Australia 1800 335 911
New Zealand 0800 167 866