


Hollister

Continence Management Systems

InView Sheaths



InView Sheaths



The secure and sensitive solution

For incontinent men suffering from urge or reflex emptying of the bladder, the InView Sheath offers unrivalled security, comfort and flexibility

Key Features

- **Two lengths to choose from**

InView Standard - unrivalled ease of use, reliability and comfort

InView Special - its shorter length and strategically placed adhesive provide a better fit for patients with a shorter or retracted penis who find standard sheath lengths difficult to manage - supplied with plastic disposable cap to maintain adhesive integrity

- **Skin-friendly self-adhesive band**

adheres gently and securely to the skin

an easy method of attachment

- **Silicone, latex-free material**

adapts to fluctuating penis size


ideal for caregivers and patients with latex sensitivity

breathable silicone is more comfortable to wear than sheaths made from some other materials

- **Clear material**

a discreet solution for the patient

permits routine inspection of the skin without removing the sheath

InView Standard	Colour Code	Size	Code No.	Box Quantity
	●	25mm	97525	30
	●	29mm	97529	30
	●	32mm	97532	30
	●	36mm	97536	30
	●	41mm	97541	30

InView Special	Colour Code	Size	Code No.	Box Quantity
	●	25mm	97425	30
	●	29mm	97429	30
	●	32mm	97432	30
	●	36mm	97436	30
	●	41mm	97441	30



Hollister Australia
6/345 Ingles Street
Port Melbourne, Victoria
Australia 3207

Hollister New Zealand
58 Richard Pearse Drive
Airport Oaks
Auckland, New Zealand

For more information, Freecall
Australia 1800 335 911
New Zealand 0800 167 866