

Moderma Flex


Continued skin protection
with each pouch application and removal

Now all colostomates can benefit from Moderma Flex with *air space technology*. A full range of closed pouches in pre-cut and cut-to-fit sizes means that Moderma Flex is the pouch of choice for all your needs.

The repeated application and removal of ostomy appliances can put the peristomal skin (the skin surrounding your stoma) at risk of damage, thus reducing the skin's natural barrier properties against the corrosive action of stomal discharge.

The Moderma Flex skin barrier with *air space technology* has numerous indentations, that allow the peristomal skin to rest and recover with each application without compromising adhesion. The skin barrier is designed to vary the position of in-contact-skin with every application.

*Moderma Flex with
air space technology*



Australia

6/345 Ingles Street
Port Melbourne, Victoria 3207

FreeCall: 1800 335 911

New Zealand

58 Richard Pearse Drive
Airport Oaks, Auckland

FreeCall 0800 167 866

www.ostomyworld.com.au


©2004 Hollister Incorporated.
Hollister and logo, CenterPointLock, Lock 'n Roll, Moderma Flex,
New Image are registered trademarks of Hollister Incorporated.

HICRPPBDL/05-04378

Colostomy Products

One and Two Piece Solutions

Better products with you in mind




Moderma Flex

Highly flexible, giving you the ultimate in confidence

- The numerous indentations that create *air space technology* allow the skin to rest and recover with each application without compromising security.
- This highly flexible skin barrier with *air space technology* provides unsurpassed comfort and freedom of movement.
- The Moderma Flex skin barrier is shaped to ensure that maximum adhesion is achieved at the point of greatest stress on the barrier.
- Designed to conform to your body shape.
- The anatomical shape provides capacity where it is required.
- Soft backing for improved comfort in everyday wear.
- Quiet film pouch material is rustle free, for greater discretion.

Increased adhesion at the point of greatest stress




A high tech protected filter - both one and two piece

- Offers superior odour control for up to 72 hours.
- Helps prevent 'ballooning' and pancaking.
- Copes well with diarrhoea from a colostomate.



Ultra-efficient, crescent shaped filter for improved flatus disposal



New Image

Highly secure and discreet, the ultimate in security

- Exceptionally low profile due to an extra thin flange and slim tapered skin barrier.
- New Image products are designed to conform to your body shape and are discreet under clothing.
- Soft backing improves your comfort, especially on hot days.
- Odour barrier, rustle-free quiet pouch film increases confidence that no odour is escaping from your pouch and provides greater discretion.
- Barrier flange is locked on both sides by the pouch flange profile, providing absolute security against accidental separation.
- An audible 'click of confidence' lets you know the system is securely locked.
- Three types of adhesive to choose from. A selection of flat and convex flanges are also available.

