


Moderma Flex

Continued skin protection
with each pouch application and removal

Now all urostomates can benefit from Moderma Flex with *air space technology*. The repeated application and removal of ostomy appliances can put the peristomal skin (the skin surrounding your stoma) at risk of damage, thus reducing the skin's natural barrier properties against the corrosive action of urine. The numerous indentations allow the peristomal skin to rest and recover.

*Moderma Flex with
air space technology*



Hollister has the **largest** selection of one and two piece urostomy pouches available on the market

These include:

- Convex skin barriers for people requiring extra security.
- A range of four skin barriers to choose from, allowing optimum wear time and choice.
- The widest range of pre-cut sizes to suit your needs.



Australia

6/345 Ingles Street
Port Melbourne, Victoria 3207

FreeCall: 1800 335 911

New Zealand

33 Vestey Drive
Mount Wellington, Auckland

FreeCall 0800 167 866

www.hollister.com.au


©2004 Hollister Incorporated.
Hollister and logo, CenterPointLock, Moderma Flex, New Image
are registered trademarks of Hollister Incorporated.

HU9RPPBDLV05-04380

Urostomy Products

One and Two Piece Solutions

Better products with you in mind




Moderma Flex

Highly flexible, giving you the ultimate in confidence

- The numerous indentations that create *air space technology* allow the skin to rest and recover with each application without compromising security.
- This highly flexible skin barrier with *air space technology* provides unsurpassed comfort and freedom of movement.
- The Moderma Flex skin barrier is shaped to ensure that maximum adhesion is achieved at the point of greatest stress on the barrier.
- Designed to conform to your body shape.
- Soft backing for improved comfort in everyday wear.
- Quiet film pouch material is rustle free, for greater discretion.

Increased adhesion at the point of greatest stress



Hygienic non-reflux valve

- Helps prevent urine backwash.
- Reduces risk of infection.
- Reduces crystalline formation on the stoma.

Easy to use drainage tap

- The simple twist action wide-bore drainage tap regulates flow and allows easy drainage of mucus and urine.
- The tube adaptor (1 per box) clips directly to the tap for easy connection to the Hollister night drainage bag and other night drainage collectors.
- A golden teardrop clearly identifies when the tap is open or closed to prevent accidental leakage.




A golden teardrop clearly identifies when the tap is open or closed to prevent accidental leakage



New Image

Highly secure and discreet, the ultimate in security

- Exceptionally low profile due to an extra thin flange and slim tapered skin barrier.
- New Image products are designed to conform to your body shape and are discreet under clothing.
- Soft backing improves your comfort, especially on hot days.
- Odour barrier, rustle-free quiet pouch film increases confidence that no odour is escaping from your pouch and provides greater discretion.
- Barrier flange is locked on both sides by the pouch flange profile, providing absolute security against accidental separation.
- An audible 'click of confidence' lets you know the system is securely locked.
- Three types of adhesive to choose from. A selection of flat and convex flanges are also available.



CenterPointLock closure

