

Wound Care Products

- 02** Moist Wound Healing Dressings
- 03** Drainage Collectors
- 03** Tube Attachment Devices



Restore Extra Thin - Sterile

Thin occlusive dressing with flexible film backing provides a moist environment and is designed to manage and protect skin and superficial wounds in difficult to manage sites, especially those at risk from friction injury. Cross-linked network of barrier molecules leads to extended wear time.

Code:	10cm x 10cm	9921
Code:	15cm x 20cm	9923
Code:	20cm x 20cm	9925



Restore Plus - Sterile

Occlusive dressing with a conformable tapered edge designed for the management of dermal ulcers and partial to full thickness wounds.

Code:	Tapered Edge - 10cm x 10cm	9956
Code:	Sacral - 110cm	9959
Code:	Tapered Edge - 15cm x 20cm	9957
Code:	Tapered Edge - 20cm x 20cm	9958



Restore CalciCare - Sterile

An absorptive dressing made from calcium sodium alginate for use in management of moderate to heavily exuding, partial to full thickness wounds.

Code:	5cm x 5cm	9938
Code:	10cm x 10cm	9937
Code:	10cm x 20cm	9939
Code:	30cm Rope - 2g	9940



Wound Drainage Collector - Sterile

Code: For wounds up to 7.6cm 9701

Code: For wounds up to 9.5cm 9702



Wound Drainage Collector - Non Sterile

Code: For wounds up to 7.6cm 9773

Code: For wounds up to 9.5cm 9776

Code: For wounds up to 10cm x 20cm 9778



Universal Catheter Access Port - UCAP - Sterile

Code: 9779



Horizontal Tube Attachment Device - HTAD - Sterile

Code: 9781



Vertical Tube Attachment Device - VTAD - Sterile

Code: 9782



Feeding Tube Attachment Device - FTAD

Code: 9786



Oral Endotracheal Tube Attachment Device - ETAD

Code: 9789



Oral Endotracheal Tube Attachment Device Neck Strap

Code: 9788

