

Tube Attachment Devices

Tube Attachment Devices are designed to meet clinicians' needs by eliminating problems inherent when sutures or tape are used to secure tubes, drains and catheters.

Vertical Tube Attachment Device (V-TAD)

- Adjustable clamp locks serrated strap in place; serrated strap wraps around and secures the drain, tube or catheter vertical to the skin.
- Exclusive FlexWear Hydrocolloid skin barrier with ElastaFlex backing protects the skin.
- Microporous tape provides additional security.



Horizontal Tube Attachment Device (H-TAD)

- Adjustable clamp locks serrated strap in place; serrated strap wraps around and secures the drain, tube or catheter horizontal to the skin.
- Exclusive FlexWear hydrocolloid skin barrier with ElastaFlex backing protects the skin.
- Microporous tape provides additional security.



Feeding Tube Attachment Device (F-TAD)

- SoftFlex hydrocolloid skin barrier gently adheres to the nose, protecting the skin and providing improved patient comfort.
- Uniquely designed clamp secures the feeding tube without kinking or occluding the tube.
- Clamp swivels to allow insertion into right or left nasal passage.



Device	Code
Horizontal Tube Attachment Device (HTAD)	9781
Vertical Tube Attachment Device (VTAD)	9782
Feeding Tube Attachment Device (FTAD)	9786
Oral Endotracheal Tube Attachment Device (ETAD) <i>(not shown)</i>	9789
Oral Endotracheal Tube Attachment Device Neck Strap <i>(not shown)</i>	9788

Restore Plus	Size	Code
	10x10cm	9956
	15x20cm	9957
	20x20cm	9958
	Sacral 110cm	9959

Restore Extra Thin	Size	Code
	10x10cm	9921
	15x20cm	9923
	20x20cm	9925

Wound Drainage Collectors	Size	Code
Sterile	up to 7.6cm	9701
	up to 9.5cm	9702
	up to 10x20cm	9703

Catheter Access Port	Code
	9779

Pouchkins	Size	Code
One Piece	cuts to 35x22mm	3778
One Piece	cuts to 16mm	3777

Restore CalciCare	Size	Code
	5x5cm	9938
	10x10cm	9937
	10x20cm	9939
	30cm Rope	9940

Wound Drainage Collectors	Size	Code
Non-Sterile	up to 7.6cm	9773
	up to 9.5cm	9776
	up to 10x20cm	9778

Accessories	Size	Code
Tube		3779



Australia

6/345 Ingles Street
Port Melbourne, Victoria 3207

FreeCall: 1800 335 911

New Zealand

58 Richard Pearse Drive
Airport Oaks, Auckland

FreeCall 0800 167 866

www.ostomyworld.com.au

©2004 Hollister Incorporated. Hollister and logo, ElastaFlex, FlexTend, FlexWear, Pouchkins, SoftFlex, Restore and CalciCare are registered trademarks of Hollister Incorporated.
Kaltostat and Aquacel are registered trademarks of ConvaTec, a Bristol-Myers Squibb Company.
Sorbsan is a trademark of Dow B Hickam Pharmaceuticals Inc. Curasorb is a trademark of Kendall Healthcare Products Co.

HWGPEPBB/DLV/03-04408



Wound Care Products

Innovations in Wound Management

The edge in wound healing





Patients with limited mobility should not have to endure the added health risks associated with confinement to a bed. This basic need led Hollister to work with healthcare professionals to develop a quality array of wound care products backed, as all Hollister products are, with our policy of unconditional customer satisfaction.

Our product innovations in wound care are extensive, ranging from wound dressings to fistula management products.

All Hollister wound care and therapy products are designed to help healthcare professionals provide quality patient care, while increasing patient comfort.

Restore Plus and Restore Extra Thin

Restore Plus and Restore Extra Thin Hydrocolloid Wound Care Dressings create a moist wound healing environment. Extended wear time, up to 7 days, results in fewer dressing changes, less skin irritation, reduced care time and increased value for money. The dressings are easily conformable to the body contours and ideal for difficult-to-manage wound sites.



- Excellent total fluid handling capacity
- Highly conformable
- Improved adhesion
- Facilitates moist wound healing and autolytic debridement
- Indicated for partial and full thickness wounds
- Easy to cut and customise to wound shape
- Easy to apply and remove
- Transparency permits visual assessment
- Occlusive and waterproof
- Lightweight dressing stretches; patient has added freedom of movement
- Protects skin from friction injury
- Ideally suited for bony prominences



Transparency permits visual Assessment of wound site without dressing removal

Restore CalciCare Alginate


If you think all alginates are alike, try Restore CalciCare... *only from Hollister.*

These dressings are significantly more absorbent. That means less trauma for the patient and reduced nursing time and costs. The exclusive manufacturing technology results in advanced fibres that maintain a moist wound interface, absorb exudate uniformly, help avoid lateral wicking and protect the periwound skin.



Superior Absorption

- Restore CalciCare Dressing is the most absorbent of the leading alginates and hydro-fibre dressings. Up to 96% more absorbent than other alginates¹



- The propriety technology used to manufacture Restore CalciCare Dressings allow a unique combination of superior absorbency without sacrificing tensile strength. As a result, you can remove the dressing in one piece.
- Restore CalciCare Dressing leaves less residue than older calcium alginate technology, avoiding the need for irrigation that might damage fragile new tissue.

1. Based on laboratory tests of Restore CalciCare versus other leading alginates and hydro-fibre dressings. Test data on file. *Using the test method of the British Pharmacopoeia, dressings were soaked for 30 minutes and then removed.



Wound Drainage Collectors

Wound Drainage Collectors designed to effectively collect and contain exudate from draining wounds and fistulae. Made of transparent Odour-Barrier film, each collector features an attached Access Cap for easy inspection and cleaning of the wound; Convenience Drain allows fluids in the collector to be emptied and measured; an attached FlexWear hydrocolloid barrier provides long-wearing skin protection and security to prevent leakage.



Universal Catheter Access Port allows the passage of a drainage catheter through the collector and provides a secure seal between the catheter and the collector.

Pouchkins

Variety of small pouches, ideal for managing small draining wounds and fistulae.

