

Advance Hydro Soft

Hydrophilic Catheter

Choices
that empower life



Advance Hydro Soft

Hydrophilic Catheter

Pure

The self-contained sterile water packet in each single-use **Advance Hydro Soft** hydrophilic catheter eliminates the worries about discreet catheter preparation and access to clean water sources.

The latex-free, silicone construction of Advance Hydro Soft catheter is gentle and designed to follow the natural contours of the urethra.

Simple

The all-in-one design and peel-back packaging make the Advance Hydro Soft catheter quick to open and use. With a simple Snap and Swish, the catheter is easy to prepare and requires no external source of clean water.

Comfortable

The catheter is made of soft silicone that is highly flexible and designed to minimize the amount of pressure required during insertion. Its hydrophilic coating provides virtually friction-free insertion and removal.

Dependable

The goal of **Hollister** is simple: to improve the quality of life for people who use our products. We tailor dependable products that help make life healthier and easier – instilling confidence and control.

Snap 'n Swish



Simply *Snap* the enclosed packet of sterile water by folding it, and *Swish* the water with a gentle seesaw motion to activate the slippery hydrophilic coating for even lubrication.



Caring for the details to help give peace of mind.

It takes a community of caregivers to provide families with hope, instill confidence and even help them regain some control. Hollister is honored to be a part of this circle of care – sharing advances through our products and services and developing trusted bonds with the people we help.

The role of Hollister is to support you by focusing on the details. After all, we understand that sometimes the smallest things can make the biggest difference. You'll find our approach reflected in the design of all our products, our services and everything we do.

For more than 80 years, Hollister has been making a difference where it matters, playing a role in making life rewarding and dignified. As we continue to develop products and services, our goal remains the same – to promote independence and quality of life.



Advance Hydro Soft

Stock No.	Tip	Length	Box Qty.
80103	Straight	10Fr/6in (15cm)	30
80123	Straight	12Fr/6in (15cm)	30
80143	Straight	14Fr/6in (15cm)	30
80163	Straight	16Fr/6in (15cm)	30
80085	Straight	08Fr/10in (25cm)	30
80105	Straight	10Fr/10in (25cm)	30
80108	Straight	10Fr/16in (40cm)	30
80128	Straight	12Fr/16in (40cm)	30
80148	Straight	14Fr/16in (40cm)	30
80168	Straight	16Fr/16in (40cm)	30
80188	Straight	18Fr/16in (40cm)	30



Hollister

Hollister Incorporated
2000 Hollister Drive
Libertyville, Illinois 60048 USA

Distributed in Canada by
Hollister Limited
95 Mary Street
Aurora, Ontario L4G 1G3 CANADA

For a free sample or more information, call Hollister at 1-800-323-4060. In Canada, call 1-800-263-7400.

www.hollister.com