

Legacy Product Migrations

Key Considerations

- The following Migration Charts are for reference only. At this point, Hollister does NOT have plans to discontinue the legacy products, such as **CenterPointLock** Two-Piece Pouching System and **Karaya** One-Piece Pouching System.
- The following product brands are the focus products from Hollister and feature our most advanced technology:
 - **New Image** Two-Piece Pouching System
 - **Premier** One-Piece Pouching System
 - **Adapt** Accessories
- When attempting to migrate from a legacy product to a focus product, it is important to realize that such a transition often cannot be made with a one-to-one match. Thus, more than one option should be considered, and all product features must be reviewed, example: suggested product may have a different skin barrier, or some pre-sized conversions may result in a cut-to-fit product.
- With Two-Piece product conversions, first select an appropriate skin barrier and then choose an appropriate pouch.
 - Flange sizes may change in some instances
 - Flange sizes must match with the barrier and pouch selected, (e.g., green to green)
 - New Image pouches must be used with New Image skin barriers. They cannot be used with the **CenterPointLock** System
- Regarding drainable pouches, the focus is on drainable pouches that feature the **Lock 'n Roll** Closure. Drainable pouches with a traditional pouch clamp may, however, be more appropriate for certain people, depending on their individual product needs.

CenterPointLock to New Image Two-Piece Pouching System

CenterPointLock Skin Barrier with Floating Flange → New Image Skin Barrier with Floating Flange

FlexTend Stock No	Extended Flange Size	Wear, with Tape, Barrier Opening	Cut-to-Fit
8727	1-1/2" (38mm)	up to 1" (25mm)	
8722	1-3/4" (44mm)	up to 1-1/4" (32mm)	
8723	2-1/4" (57mm)	up to 1-3/4" (44mm)	
8724	2-3/4" (70mm)	up to 2-1/4" (57mm)	

FlexTend Stock No	Extended Flange Size	Wear, without Tape, Barrier Opening	Cut-to-Fit
8752	1-3/4" (44mm)	up to 1-1/4" (32mm)	
8753	2-1/4" (57mm)	up to 1-3/4" (44mm)	
8754	2-3/4" (70mm)	up to 2-1/4" (57mm)	

FlexWear Stock No	Standard Flange Size	Wear, with Tape, Barrier Opening	Cut-to-Fit
3727	1-1/2" (38mm)	up to 1" (25mm)	
3722	1-3/4" (44mm)	up to 1-1/4" (32mm)	
3723	2-1/4" (57mm)	up to 1-3/4" (44mm)	
3724	2-3/4" (70mm)	up to 2-1/4" (57mm)	
3726	4" (102mm)	up to 3-1/2" (89mm)	

SoftFlex Stock No	Standard Flange Size	Wear, without Tape, Barrier Opening	Cut-to-Fit
3767	1-1/2" (38mm)	up to 1" (25mm)	
3762	1-3/4" (44mm)	up to 1-1/4" (32mm)	
3763	2-1/4" (57mm)	up to 1-3/4" (44mm)	
3764	2-3/4" (70mm)	up to 2-1/4" (57mm)	
3766	4" (102mm)	up to 3-1/2" (89mm)	

HolliHesive Stock No	Standard Flange Size	Wear, with Tape, Barrier Opening	Cut-to-Fit
3707	1-1/2" (38mm)	up to 3/4" (19mm)	
3702	1-3/4" (44mm)	up to 1" (25mm)	
3703	2-1/4" (57mm)	up to 1-1/2" (38mm)	
3704	2-3/4" (70mm)	up to 2" (51mm)	
3706	4" (102mm)	up to 3-1/2" (89mm)	

FlexTend Stock No	Extended Flange Size	Wear, with Tape, Barrier Opening	Pre-Sized
8743	1-1/2" (38mm)	7/8" (22mm)	
8744	1-1/2" (38mm)	1" (25mm)	
8745	1-3/4" (44mm)	1-1/8" (29mm)	
8746	1-3/4" (44mm)	1-1/4" (32mm)	
8747	2-1/4" (57mm)	1-3/8" (35mm)	
8748	2-1/4" (57mm)	1-1/2" (38mm)	

FlexWear Stock No	Standard Flange Size	Wear, with Tape, Barrier Opening	Pre-Sized
3747	1-1/2" (38mm)	3/4" (19mm)	
3742	1-1/2" (38mm)	1" (25mm)	
3748	1-3/4" (44mm)	1-1/4" (32mm)	
3743	2-1/4" (57mm)	1-1/2" (38mm)	
3749	2-1/4" (57mm)	1-3/4" (44mm)	

FlexTend Stock No	Extended Flange Size	Wear, with Tape, Barrier Opening	Cut-to-Fit
14602	1-3/4" (44mm)	up to 1-1/4" (32mm)	
14602	1-3/4" (44mm)	up to 1-1/4" (32mm)	
14603	2-1/4" (57mm)	up to 1-3/4" (44mm)	
14604	2-3/4" (70mm)	up to 2-1/4" (57mm)	

FlexTend Stock No	Extended Flange Size	Wear, without Tape, Barrier Opening	Cut-to-Fit
15602	1-3/4" (44mm)	up to 1-1/4" (32mm)	
15603	2-1/4" (57mm)	up to 1-3/4" (44mm)	
15604	2-3/4" (70mm)	up to 2-1/4" (57mm)	

FlexWear Stock No	Standard Flange Size	Wear, with Tape, Barrier Opening	Cut-to-Fit
14202	1-3/4" (44mm)	up to 1-1/4" (32mm)	
14202	1-3/4" (44mm)	up to 1-1/4" (32mm)	
14203	2-1/4" (57mm)	up to 1-3/4" (44mm)	
14204	2-3/4" (70mm)	up to 2-1/4" (57mm)	
14206	4" (102mm)	up to 3-1/2" (89mm)	

FlexWear Stock No	Standard Flange Size	Wear, without Tape, Barrier Opening	Cut-to-Fit
15202	1-3/4" (44mm)	up to 1-1/4" (32mm)	
15202	1-3/4" (44mm)	up to 1-1/4" (32mm)	
15203	2-1/4" (57mm)	up to 1-3/4" (44mm)	
15204	2-3/4" (70mm)	up to 2-1/4" (57mm)	
14206	4" (102mm)	up to 3-1/2" (89mm)	

FlexWear Stock No	Standard Flange Size	Wear, with Tape, Barrier Opening	Cut-to-Fit
14202	1-3/4" (44mm)	up to 1-1/4" (32mm)	
14202	1-3/4" (44mm)	up to 1-1/4" (32mm)	
14203	2-1/4" (57mm)	up to 1-3/4" (44mm)	
14204	2-3/4" (70mm)	up to 2-1/4" (57mm)	
14206	4" (102mm)	up to 3-1/2" (89mm)	

FlexTend Stock No	Extended Flange Size	Wear, with Tape, Barrier Opening	Pre-Sized
14703	1-3/4" (44mm)	7/8" (22mm)	
14704	1-3/4" (44mm)	1" (25mm)	
14705	1-3/4" (44mm)	1-1/8" (29mm)	
14706	1-3/4" (44mm)	1-1/4" (32mm)	
14707	2-1/4" (57mm)	1-3/8" (35mm)	
14708	2-1/4" (57mm)	1-1/2" (38mm)	

FlexWear Stock No	Standard Flange Size	Wear, with Tape, Barrier Opening	Pre-Sized
14302	1-3/4" (44mm)	3/4" (19mm)	
14304	1-3/4" (44mm)	1" (25mm)	
14306	1-3/4" (44mm)	1-1/4" (32mm)	
14308	2-1/4" (57mm)	1-1/2" (38mm)	
14309	2-1/4" (57mm)	1-3/4" (44mm)	

FlexTend Stock No	Extended Flange Size	Wear, without Tape, Barrier Opening	Cut-to-Fit
15602	1-3/4" (44mm)	up to 1-1/4" (32mm)	
15602	1-3/4" (44mm)	up to 1-1/4" (32mm)	
15603	2-1/4" (57mm)	up to 1-3/4" (44mm)	
15604	2-3/4" (70mm)	up to 2-1/4" (57mm)	

FlexTend Stock No	Extended Flange Size	Wear, with Tape, Barrier Opening	Cut-to-Fit
14602	1-3/4" (44mm)	up to 1-1/4" (32mm)	
14603	2-1/4" (57mm)	up to 1-3/4" (44mm)	
14604	2-3/4" (70mm)	up to 2-1/4" (57mm)	

FlexWear Stock No	Standard Flange Size	Wear, without Tape, Barrier Opening	Cut-to-Fit
15202	1-3/4" (44mm)	up to 1-1/4" (32mm)	
15202	1-3/4" (44mm)	up to 1-1/4" (32mm)	
15203	2-1/4" (57mm)	up to 1-3/4" (44mm)	
15204	2-3/4" (70mm)	up to 2-1/4" (57mm)	

FlexWear Stock No	Standard Flange Size	Wear, with Tape, Barrier Opening	Cut-to-Fit
14202	1-3/4" (44mm)	up to 1-1/4" (32mm)	
14202	1-3/4" (44mm)	up to 1-1/4" (32mm)	
14203	2-1/4" (57mm)	up to 1-3/4" (44mm)	
14204	2-3/4" (70mm)	up to 2-1/4" (57mm)	

FlexWear Stock No	Standard Flange Size	Wear, without Tape, Barrier Opening	Cut-to-Fit
15202	1-3/4" (44mm)	up to 1-1/4" (32mm)	
15202	1-3/4" (44mm)	up to 1-1/4" (32mm)	
15203	2-1/4" (57mm)	up to 1-3/4" (44mm)	
15204	2-3/4" (70mm)	up to 2-1/4" (57mm)	

FlexTend Stock No	Extended Flange Size	Wear, without Tape, Barrier Opening	Pre-Sized
16103	1-3/4" (44mm)	7/8" (22mm)	
16104	1-3/4" (44mm)	1" (25mm)	
16105	1-3/4" (44mm)	1-1/8" (29mm)	
16106	1-3/4" (44mm)	1-1/4" (32mm)	
16107	2-1/4" (57mm)	1-3/8" (35mm)	
16108	2-1/4" (57mm)	1-1/2" (38mm)	

FlexTend Stock No	Extended Flange Size	Wear, with Tape, Barrier Opening	Pre-Sized
14702	1-3/4" (44mm)	3/4" (19mm)	
14704	1-3/4" (44mm)	1" (25mm)	
14706	1-3/4" (44mm)	1-1/4" (32mm)	
14708	2-1/4" (57mm)	1-1/2" (38mm)	
14709	2-1/4" (57mm)	1-3/4" (44mm)	

Comments

Change in flange size

Change in flange size

Change in flange size

Change in flange size

Change in flange size

Change in flange size

Change in flange size

Change in flange size

Change to cut-to-fit; FlexWear standard wear only

CenterPointLock to New Image Two-Piece Pouching System (continued)

CenterPointLock Convex Skin Barrier with Floating Flange

FlexTend Stock No	Extended Wear, with Tape, Cut-to-Fit Flange Size	Barrier Opening
8774	1-3/4" (44mm)	up to 1" (25mm)
8778	2-1/4" (57mm)	up to 1-1/2" (38mm)
87711	2-3/4" (70mm)	up to 2" (51mm)

FlexTend Stock No	Extended Wear, with Tape, Pre-Sized Flange Size	Barrier Opening
8731	1-1/2" (38mm)	5/8" (16mm)
8732	1-1/2" (38mm)	3/4" (19mm)
8733	1-3/4" (44mm)	7/8" (22mm)
8734	1-3/4" (44mm)	1" (25mm)
8735	2-1/4" (57mm)	1-1/8" (29mm)
8736	2-1/4" (57mm)	1-1/4" (32mm)
8737	2-1/4" (57mm)	1-3/8" (35mm)
8738	2-1/4" (57mm)	1-1/2" (38mm)

FlexWear Stock No	Standard Wear, with Tape, Pre-Sized Flange Size	Barrier Opening
3730	1-1/2" (38mm)	1/2" (13mm)
3731	1-1/2" (38mm)	5/8" (16mm)
3732	1-1/2" (38mm)	3/4" (19mm)
3733	1-3/4" (44mm)	7/8" (22mm)
3734	1-3/4" (44mm)	1" (25mm)
3735	2-1/4" (57mm)	1-1/8" (29mm)
3736	2-1/4" (57mm)	1-1/4" (32mm)
3737	2-1/4" (57mm)	1-3/8" (35mm)
3738	2-1/4" (57mm)	1-1/2" (38mm)
3739	2-3/4" (70mm)	1-5/8" (41mm)
37310	2-3/4" (70mm)	1-3/4" (44mm)
37311	2-3/4" (70mm)	2" (51mm)

CenterPointLock Drainable Pouch with Replaceable Filter

Pouch Clamp Closure Stock No	Flange Size	Length	Color
3897	1-1/2" (38mm)	12"	Transparent
3892	1-3/4" (44mm)	12"	Transparent
3893	2-1/4" (57mm)	12"	Transparent
3894	2-3/4" (70mm)	12"	Transparent
3837	1-1/2" (38mm)	12"	Beige
3832	1-3/4" (44mm)	12"	Beige
3833	2-1/4" (57mm)	12"	Beige
3834	2-3/4" (70mm)	12"	Beige

CenterPointLock Drainable Pouch, without Filter

Pouch Clamp Closure Stock No	Flange Size	Length	Color
8817	1-1/2" (38mm)	12"	Beige
8812	1-3/4" (44mm)	12"	Beige
8813	2-1/4" (57mm)	12"	Beige
8814	2-3/4" (70mm)	12"	Beige

Pouch Clamp Closure Stock No	Flange Size	Length	Color
3807	1-1/2" (38mm)	11"	Transparent
3802	1-3/4" (44mm)	11"	Transparent
3803	2-1/4" (57mm)	11"	Transparent
3804	2-3/4" (70mm)	11"	Transparent
3806	4" (102mm)	11"	Transparent
3817	1-1/2" (38mm)	11"	Beige
3812	1-3/4" (44mm)	11"	Beige
3813	2-1/4" (57mm)	11"	Beige
3814	2-3/4" (70mm)	11"	Beige

Pouch Clamp Closure Stock No	Flange Size	Length	Color
3872	1-3/4" (44mm)	16"	Transparent
3873	2-1/4" (57mm)	16"	Transparent
3874	2-3/4" (70mm)	16"	Transparent

CenterPointLock Drainable Mini-Pouch, without Filter

Pouch Clamp Closure Stock No	Flange Size	Length	Color
3847	1-1/2" (38mm)	9"	Beige
3842	1-3/4" (44mm)	9"	Beige
3843	2-1/4" (57mm)	9"	Beige
3844	2-3/4" (70mm)	9"	Beige

→ New Image Convex Skin Barrier with Floating Flange

FlexTend Stock No	Extended Wear, with Tape, Cut-to-Fit Flange Size	Barrier Opening
14802	1-3/4" (44mm)	up to 1" (25mm)
14803	2-1/4" (57mm)	up to 1-1/2" (38mm)
14804	2-3/4" (70mm)	up to 2" (51mm)

FlexTend Stock No	Extended Wear, with Tape, Pre-Sized Flange Size	Barrier Opening
14901	1-3/4" (44mm)	5/8" (16mm)
14902	1-3/4" (44mm)	3/4" (19mm)
14903	1-3/4" (44mm)	7/8" (22mm)
14904	1-3/4" (44mm)	1" (25mm)
14905	2-1/4" (57mm)	1-1/8" (29mm)
14906	2-1/4" (57mm)	1-1/4" (32mm)
14907	2-1/4" (57mm)	1-3/8" (35mm)
14908	2-1/4" (57mm)	1-1/2" (38mm)

FlexWear Stock No	Standard Wear, with Tape, Pre-Sized Flange Size	Barrier Opening
14402	1-3/4" (44mm)	1/2" (13mm)
14402	1-3/4" (44mm)	up to 1" (25mm)
14402	1-3/4" (44mm)	up to 1" (25mm)
14503	1-3/4" (44mm)	7/8" (22mm)
14504	1-3/4" (44mm)	1" (25mm)
14505	2-1/4" (57mm)	1-1/8" (29mm)
14506	2-1/4" (57mm)	1-1/4" (32mm)
14404	2-3/4" (70mm)	up to 2" (51mm)
14404	2-3/4" (70mm)	up to 2" (51mm)
14404	2-3/4" (70mm)	up to 2" (51mm)
14404	2-3/4" (70mm)	up to 2" (51mm)

→ New Image Drainable Pouch, with Filter

Choice of Closure Lock 'n Roll	Clamp	Flange Size	Length	Color
18192	18162	1-3/4" (44mm)	12"	Transparent
18192	18162	1-3/4" (44mm)	12"	Transparent
18193	18163	2-1/4" (57mm)	12"	Transparent
18194	18164	2-3/4" (70mm)	12"	Transparent
18182	18142	1-3/4" (44mm)	12"	Beige
18182	18142	1-3/4" (44mm)	12"	Beige
18183	18143	2-1/4" (57mm)	12"	Beige
18184	18144	2-3/4" (70mm)	12"	Beige

→ New Image Drainable Pouch, without Filter

Choice of Closure Lock 'n Roll	Clamp	Flange Size	Length	Color
18112	18122	1-3/4" (44mm)	12"	Beige
18112	18122	1-3/4" (44mm)	12"	Beige
18113	18123	2-1/4" (57mm)	12"	Beige
18114	18124	2-3/4" (70mm)	12"	Beige

Choice of Closure Lock 'n Roll	Clamp	Flange Size	Length	Color
18132	18102	1-3/4" (44mm)	12"	Transparent
18132	18102	1-3/4" (44mm)	12"	Ultra-clear
18003	18173	2-1/4" (57mm)	12"	Ultra-clear
18004	18174	2-3/4" (70mm)	12"	Ultra-clear
18006	18176	4" (102mm)	12"	Ultra-clear
18112	18122	1-3/4" (44mm)	12"	Beige
18112	18122	1-3/4" (44mm)	12"	Beige
18113	18123	2-1/4" (57mm)	12"	Beige
18114	18124	2-3/4" (70mm)	12"	Beige

Choice of Closure Soft Tap	Lock 'n Roll	Flange Size	Length	Color
18012*	18132	1-3/4" (44mm)	12"	Transparent
18013*	18133	2-1/4" (57mm)	12"	Transparent
18014*	18134	2-3/4" (70mm)	12"	Transparent

→ New Image Drainable Mini-Pouch, without Filter

Clamp	Flange Size	Length	Color
18202	1-3/4" (44mm)	9"	Beige
18202	1-3/4" (44mm)	9"	Beige
18203	2-1/4" (57mm)	9"	Beige
18204	2-3/4" (70mm)	9"	Beige

Comments

*Change to cut-to-fit and flange size

*Change to cut-to-fit and flange size

*Change to cut-to-fit and flange size

*Change to cut-to-fit and flange size

Change in flange size

Change to cut-to-fit and flange size

Change to cut-to-fit and flange size

Change to cut-to-fit and flange size

Change to cut-to-fit and flange size

Change to cut-to-fit

Change to cut-to-fit

Change to cut-to-fit

Change to cut-to-fit

Change in flange size

Change in flange size

Change in flange size

Change in flange size

Change in flange size

Shorter length, * high output pouch

Shorter length, * high output pouch

Shorter length, * high output pouch

Change in flange size

CenterPointLock to New Image Two-Piece Pouching System (continued)

CenterPointLock Closed Pouch, with Filter				→	New Image Closed Pouch, with Filter				Comments
Stock No	Flange Size	Length	Color		Stock No	Flange Size	Length	Color	
8347	1-1/2" (38mm)	7.5"	Beige		18392	1-3/4" (44mm)	7"	Beige	Change in flange size
8342	1-3/4" (44mm)	7.5"	Beige		18392	1-3/4" (44mm)	7"	Beige	
8343	2-1/4" (57mm)	7.5"	Beige		18393	2-1/4" (57mm)	7"	Beige	
8344	2-3/4" (70mm)	7.5"	Beige		18394	2-3/4" (70mm)	7"	Beige	
3347	1-1/2" (38mm)	9"	Beige		18372	1-3/4" (44mm)	9"	Beige	Change in flange size
3342	1-3/4" (44mm)	9"	Beige		18372	1-3/4" (44mm)	9"	Beige	
3343	2-1/4" (57mm)	9"	Beige		18373	2-1/4" (57mm)	9"	Beige	
3344	2-3/4" (70mm)	9"	Beige		18374	2-3/4" (70mm)	9"	Beige	

CenterPointLock Urostomy Pouch → New Image Urostomy Pouch

Stock No	Flange Size	Length	Color		Stock No	Flange Size	Length	Color	
8902	1-3/4" (44mm)	9"	Beige		18412	1-3/4" (44mm)	9"	Beige	
8903	2-1/4" (57mm)	9"	Beige		18413	2-1/4" (57mm)	9"	Beige	
3907	1-1/2" (38mm)	9"	Transparent		18402	1-3/4" (44mm)	9"	Transparent	Change in flange size; has ComfortWear panels
3902	1-3/4" (44mm)	9"	Transparent		18402	1-3/4" (44mm)	9"	Transparent	Has ComfortWear panels
3903	2-1/4" (57mm)	9"	Transparent		18403	2-1/4" (57mm)	9"	Transparent	Has ComfortWear panels
3904	2-3/4" (70mm)	9"	Transparent		18404	2-3/4" (70mm)	9"	Transparent	Has ComfortWear panels

CenterPointLock Stoma Cap, with Filter → New Image Closed Mini-Pouch, with Filter

Stock No	Flange Size	Length	Color		Stock No	Flange Size	Length	Color	
3407	1-1/2" (38mm)	4.5"	Beige		18392	1-3/4" (44mm)	7"	Beige	Longer length, Change in flange size
3402	1-3/4" (44mm)	4.5"	Beige		18392	1-3/4" (44mm)	7"	Beige	Longer length
3403	2-1/4" (57mm)	4.5"	Beige		18393	2-1/4" (57mm)	7"	Beige	Longer length
3404	2-3/4" (70mm)	4.5"	Beige		18394	2-3/4" (70mm)	7"	Beige	Longer length

CenterPointLock Irrigator Sleeve → New Image Irrigator Sleeve

Stock No	Flange Size	Length	Color		Stock No	Flange Size	Length	Color	
3827	1-1/2" (38mm)	30"	Transparent		18152	1-3/4" (44mm)	30"	Transparent	Change in flange size
3822	1-3/4" (44mm)	30"	Transparent		18152	1-3/4" (44mm)	30"	Transparent	
3823	2-1/4" (57mm)	30"	Transparent		18153	2-1/4" (57mm)	30"	Transparent	
3824	2-3/4" (70mm)	30"	Transparent		18154	2-3/4" (70mm)	30"	Transparent	
3826	4" (102mm)	30"	Transparent		18154	2-3/4" (70mm)	30"	Transparent	Change in flange size

Karaya 5 to Premier One-Piece Pouching System

Karaya 5 One-Piece Pre-Sized Drainable Pouch → Premier One-Piece Convex Pre-Sized Drainable Pouch

Karaya 5 Standard Wear, with Replaceable Filter, with Tape, Pouch Clamp Closure				Flextend Extended Wear, without Filter, with Tape					
Stock No	Barrier Opening	Length	Color	Lock 'n Roll	Pouch Clamp	Barrier Opening	Length	Color	
3662	5/8" (15mm)	12"	Transparent	8524	8624	up to 1" (25mm)	12"	Transparent	Change to extended wear barrier
3668	7/8" (22mm)	12"	Transparent	8511	8611	7/8" (22mm)	12"	Transparent	Change to extended wear barrier
3663	1-1/8" (29mm)	12"	Transparent	8513	8613	1-1/8" (29mm)	12"	Transparent	Change to extended wear barrier
3669	1-3/8" (35mm)	12"	Transparent	8515	8615	1-3/8" (35mm)	12"	Transparent	Change to extended wear barrier
3664	1-1/2" (38mm)	12"	Transparent	8528*	8616	1-1/2" (38mm)	12"	Transparent	Change to extended wear barrier, * Lock 'n Roll closure change to cut-to-fit up to 1-1/2" (38mm)
3665	2" (51mm)	12"	Transparent	85211*	8619	2" (51mm)	12"	Transparent	Change to extended wear barrier, * Lock 'n Roll closure change to cut-to-fit up to 2" (51mm)

Karaya 5 Standard Wear, without Filter, with Tape, → Flextend Extended Wear, without Filter, with Tape

Karaya 5 Standard Wear, without Filter, with Tape, Pouch Clamp Closure				Flextend Extended Wear, without Filter, with Tape					
Stock No	Barrier Opening	Length	Color	Lock 'n Roll	Pouch Clamp	Barrier Opening	Length	Color	
3608	7/8" (22mm)	12"	Transparent	8511	8611	7/8" (22mm)	12"	Transparent	Change to extended wear barrier
3603	1-1/8" (29mm)	12"	Transparent	8513	8613	1-1/8" (29mm)	12"	Transparent	Change to extended wear barrier
3609	1-3/8" (35mm)	12"	Transparent	8515	8615	1-3/8" (35mm)	12"	Transparent	Change to extended wear barrier
3604	1-1/2" (38mm)	12"	Transparent	8528*	8616	1-1/2" (38mm)	12"	Transparent	Change to extended wear barrier, *Lock 'n Roll closure change to cut-to-fit up to 1-1/2" (38mm)
3605	2" (51mm)	12"	Transparent	85211*	8619	2" (51mm)	12"	Transparent	Change to extended wear barrier, *Lock 'n Roll closure change to cut-to-fit up to 2" (51mm)
3606	2-1/2" (64mm)	12"	Transparent	8531* (Flat)	8635 (Flat)	2-1/2" (64mm)	12"	Transparent	Change to extended wear barrier, *Lock 'n Roll closure change to cut-to-fit up to 2-1/2" (64mm)

Karaya 5 One-Piece Pre-Sized Drainable Mini-Pouch → Premier One-Piece cut-to-fit Drainable Mini-Pouch

Karaya 5 Standard Wear, without Filter, with Tape, Pouch Clamp Closure				SoftFlex Standard Wear, with Filter, without Tape					
Stock No	Barrier Opening	Length	Color	Lock 'n Roll	Barrier Opening	Length	Color		
3232	5/8" (15mm)	9"	Transparent	88800	up to 2-1/8" (55mm)	7"	Transparent	Change to cut-to-fit, with filter, with Lock 'n Roll closure	
3238	7/8" (22mm)	9"	Transparent	88800	up to 2-1/8" (55mm)	7"	Transparent	Change to cut-to-fit, with filter, with Lock 'n Roll closure	
3233	1-1/8" (29mm)	9"	Transparent	88800	up to 2-1/8" (55mm)	7"	Transparent	Change to cut-to-fit, with filter, with Lock 'n Roll closure	
3239	1-3/8" (35mm)	9"	Transparent	88800	up to 2-1/8" (55mm)	7"	Transparent	Change to cut-to-fit, with filter, with Lock 'n Roll closure	
3234	1-1/2" (38mm)	9"	Transparent	88800	up to 2-1/8" (55mm)	7"	Transparent	Change to cut-to-fit, with filter, with Lock 'n Roll closure	

Karaya 5 to Premier One-Piece Pouching System (continued)

Karaya 5 One-Piece Pre-Sized Drainable Pouch

Karaya 5 Standard Wear, without Filter, with Tape, Pouch Clamp Closure

Stock No	Barrier Opening	Length	Color
3222	5/8" (15mm)	12"	Transparent
3228	7/8" (22mm)	12"	Transparent
3223	1-1/8" (29mm)	12"	Transparent
3229	1-3/8" (35mm)	12"	Transparent
3224	1-1/2" (38mm)	12"	Transparent
3225	2" (51mm)	12"	Transparent
3226	2-1/2" (64mm)	12"	Transparent
3278	7/8" (22mm)	16"	Transparent
3273	1-1/8" (29mm)	16"	Transparent
3279	1-3/8" (35mm)	16"	Transparent
3274	1-1/2" (38mm)	16"	Transparent
3275	2" (51mm)	16"	Transparent
3276	2-1/2" (64mm)	16"	Transparent

→ Premier One-Piece Convex Pre-Sized Drainable Pouch

Flextend Extended Wear, with Filter, with Tape

Lock 'n Roll	Pouch Clamp	Barrier Opening	Length	Color
8524	8624	up to 1" (25mm)	12"	Transparent
8511	8611	7/8" (22mm)	12"	Transparent
8513	8613	1-1/8" (29mm)	12"	Transparent
8515	8615	1-3/8" (35mm)	12"	Transparent
8528*	8616	1-1/2" (38mm)	12"	Transparent
85211*	8619	2" (51mm)	12"	Transparent
8531* (Flat)	8635 (Flat)	2-1/2" (64mm)	12"	Transparent
8511	8611	7/8" (22mm)	12"	Transparent
8513	8613	1-1/8" (29mm)	12"	Transparent
8515	8615	1-3/8" (35mm)	12"	Transparent
8528*	8616	1-1/2" (38mm)	12"	Transparent
85211*	8619	2" (51mm)	12"	Transparent
8531* (Flat)	8635 (Flat)	2-1/2" (64mm)	12"	Transparent

Comments

Change to extended wear barrier and cut-to-fit
 Change to extended wear barrier
 Change to extended wear barrier
 Change to extended wear barrier
 Change to extended wear barrier, * Lock 'n Roll closure change to cut-to-fit up to 1-1/2" (38mm)
 Change to extended wear barrier, * Lock 'n Roll closure change to cut-to-fit up to 2" (51mm)
 Change to extended wear barrier, * Lock 'n Roll closure change to cut-to-fit up to 2-1/2" (64mm)
 Change to extended wear barrier and shorter length
 Change to extended wear barrier and shorter length
 Change to extended wear barrier and shorter length
 Change to extended wear barrier, * Lock 'n Roll closure change to cut-to-fit up to 1-1/2" (38mm) and shorter length
 Change to extended wear barrier, * Lock 'n Roll closure change to cut-to-fit up to 2" (51mm) and shorter length
 Change to extended wear barrier, * Lock 'n Roll closure change to cut-to-fit up to 2-1/2" (64mm) and shorter length

Karaya 5 One-Piece Pre-Sized Closed Pouch

Karaya 5 Standard Wear, with Filter, with Tape

Stock No	Barrier Opening	Length	Color
3328	7/8" (22mm)	9"	Transparent
3323	1-1/8" (29mm)	9"	Transparent
3329	1-3/8" (35mm)	9"	Transparent
3324	1-1/2" (38mm)	9"	Transparent
3325	2" (51mm)	9"	Transparent
3326	2-1/2" (64mm)	9"	Transparent

Karaya 5 Standard Wear, without Filter, with Tape

Stock No	Barrier Opening	Length	Color
7168	7/8" (22mm)	9"	Transparent
7163	1-1/8" (29mm)	9"	Transparent
7169	1-3/8" (35mm)	9"	Transparent
7164	1-1/2" (38mm)	9"	Transparent
7165	2" (51mm)	9"	Transparent
7166	2-1/2" (64mm)	9"	Transparent

→ Premier One-Piece Cut-to-Fit Closed Pouch

SoftFlex Standard Wear, with Filter, without Tape

Transparent	Beige	Barrier Opening	Length	Color
82400	82300	up to 2-1/8" (55mm)	9"	Transparent
82400	82300	up to 2-1/8" (55mm)	9"	Transparent
82400	82300	up to 2-1/8" (55mm)	9"	Transparent
82400	82300	up to 2-1/8" (55mm)	9"	Transparent
82400	82300	up to 2-1/8" (55mm)	9"	Transparent
82400	82300	up to 2-1/8" (55mm)	9"	Transparent

SoftFlex Standard Wear, without Filter, without Tape

Transparent	Beige	Barrier Opening	Length	Color
82400	82300	up to 2-1/8" (55mm)	9"	Transparent
82400	82300	up to 2-1/8" (55mm)	9"	Transparent
82400	82300	up to 2-1/8" (55mm)	9"	Transparent
82400	82300	up to 2-1/8" (55mm)	9"	Transparent
82400	82300	up to 2-1/8" (55mm)	9"	Transparent
82400	82300	up to 2-1/8" (55mm)	9"	Transparent

Change to cut-to-fit
 Change to cut-to-fit
 Change to cut-to-fit
 Change to cut-to-fit
 Change to cut-to-fit
 Change to cut-to-fit

Change to cut-to-fit with filter
 Change to cut-to-fit with filter
 Change to cut-to-fit with filter
 Change to cut-to-fit with filter
 Change to cut-to-fit with filter
 Change to cut-to-fit with filter

Karaya 5 One-Piece Pre-Sized Urostomy Pouch

Karaya 5 Standard Wear, with Tape

Stock No	Barrier Opening	Length	Color
1432	5/8" (15mm)	9"	Transparent
1438	7/8" (22mm)	9"	Transparent
1433	1-1/8" (29mm)	9"	Transparent
1439	1-3/8" (35mm)	9"	Transparent
1434	1-1/2" (38mm)	9"	Transparent

→ Premier One-Piece Pre-Sized Urostomy Pouch

Flextend Extended Wear, with Tape

Flat	Convex	Barrier Opening	Length	Color
8460*	8481	5/8" (15mm)	9"	Transparent
8460*	8483	7/8" (22mm)	9"	Transparent
8460*	8485	1-1/8" (29mm)	9"	Transparent
8460*	8487	1-3/8" (35mm)	9"	Transparent
8463	8488	1-1/2" (38mm)	9"	Transparent

Change to extended wear barrier, *Flat is cut-to-fit up to 2-1/2" (64mm)
 Change to extended wear barrier, *Flat is cut-to-fit up to 2-1/2" (64mm)
 Change to extended wear barrier, *Flat is cut-to-fit up to 2-1/2" (64mm)
 Change to extended wear barrier, *Flat is cut-to-fit up to 2-1/2" (64mm)
 Change to extended wear barrier