

# Leg Bags, Drainage Bag and Accessories




# Leg Bags, Drainage Bag and Accessories

More choices,  
more comfort



The vented system with a unique vent and vacuum relief valve allows air into the system, and is designed to eliminate the “vacuum effect”.



Caring for the details makes all the difference

As a healthcare professional, you develop a bond of trust with those you help and their families – throughout their lives. The role of Hollister is to support you in these vital relationships by focusing on the details. Experience has taught us that sometimes the smallest things make the greatest difference. You’ll find this approach reflected in the design of all our products, our services and everything we do.


For more than 80 years, Hollister has been making a difference where it matters, playing a role in making life more rewarding and dignified. As we continue to develop products and services, our goal remains unchanged – to promote independence and quality of life.

**Hollister Urinary Leg Bags and Bedside Drainage Collection Systems suit a wide range of needs with discretion, security and comfort.**

## Leg Bags

The soft, pleated sides of our vented and non-vented leg bags rest discreetly against the leg, maintaining a low profile as the bag expands. The anti-reflux valve prevents urine back-flow, which helps reduce the risk of bacterial growth and urinary tract infection. The wide fabric leg straps hold the bag securely and comfortably in place, without constricting. In addition, the vacuum relief valve of our vented leg bag system reduces the chance for vacuum formation.

- 1 Unique Vacuum Relief Valve**
  - Allows urine to drain freely and helps prevent the “vacuum effect”
- 2 Vent**
  - Helps prevent air build-up in the leg bag providing a lower profile
- 3 Extension Tubing**
  - Our latex-free oval extension tubing has an anti-kink feature that helps keep the fluid path open. Also available in round latex. Sterile and non-sterile available.
- 4 Fabric Straps**
  - Wide, latex-free straps hold the bag comfortably and securely
- 5 Convenience Drain**
  - Easy to operate with one hand, even for those with limited dexterity
- 6 Buttons**
  - Low profile buttons help secure the bag
- 7 Anti-Reflux Valve**
  - Controls urine back-flow up the tubing
- 8 Soft Vinyl Film**
  - Comfortable and limits unwanted odor



The collector is easy to empty and the drain outlet tube stores conveniently in the attached holder.

### Bedside Drainage Collection System

This collector – which holds volumes of up to 2000mL – meets your needs for dependable functionality and cost-effectiveness.

- ① **Sterile Fluid Path**
  - Provides patient protection
- ② **Sample Port**
  - Self-sealing latex-free sample port
- ③ **Graduations**
  - Enable visual assessment of urine output

- ④ **Integrated Drain Tube Holder**
  - Allows for easy access and storage of outlet tube
- ⑤ **Tapered Connector**
  - Compatible with external catheters and other collection devices
- ⑥ **Sheet Clip and Bed Hanger**
  - Provides additional security



Convenient self-sealing sample port.



Cleaner/Decrystallizer helps reduce unwanted odor.

## Urinary Leg Bags

### Vented Latex-Free Urinary Leg Bag Combination Pack

Includes: (1) Vented Urinary Leg Bag with sterile fluid pathway – medium or large, (1) Latex-Free, kink-resistant extension tubing with vacuum relief valve, (1) Pair fabric leg bag straps – medium (23")

Stock No	Size	Box Qty
9645	Medium (18 oz/540mL)	Box of 10
9655	Large (30 oz/900mL)	Box of 10

### Vented Latex-Free Urinary Leg Bag System

Includes: (1) Vented Urinary Leg Bag with sterile fluid pathway – medium or large, (1) Latex-Free, kink-resistant extension tubing with vacuum relief valve

Stock No	Size	Box Qty
9845	Medium (18 oz/540mL)	Box of 10
9855	Large (30 oz/900mL)	Box of 10

### Latex-Free Urinary Leg Bag Combination Pack

Includes: (1) Urinary Leg Bag – medium or large, (1) Latex-Free extension tubing (18") with connector, (1) Pair fabric leg bag straps – medium (23")

Stock No	Size	Box Qty
9348	Medium (18 oz/540mL)	Box of 12
9349	Large (30 oz/900mL)	Box of 12

### Latex-Free Urinary Leg Bag and Tubing System

Includes: (1) Urinary Leg Bag – medium or large, (1) Latex-Free extension tubing (18") with connector

Stock No	Size	Box Qty
9356	Medium (18 oz/540mL)	Box of 12
9357	Large (30 oz/900mL)	Box of 12

### Latex-Free Urinary Leg Bag and Strap System

Includes: (1) Urinary Leg Bag – medium or large, (1) Pair fabric leg bag straps – medium (23")

Stock No	Size	Box Qty
9354	Medium (18 oz./540mL)	Box of 12
9355	Large (30 oz/900mL)	Box of 12

### Latex-Free Urinary Leg Bag

Includes: (1) Urinary Leg Bag – medium or large

Stock No	Size	Box Qty
9814	Medium (18 oz/540mL)	Box of 10
9805	Large (30 oz/900mL)	Box of 10

### Urinary Leg Bag 4 Pack

Includes: (4) Urinary Leg Bags – medium or large, (1) Latex extension tubing with connector, (1) Pair fabric leg bag straps – medium (23")

Stock No	Size	Box Qty
9834	Medium (18 oz/540mL)	Single Unit
9835	Large (30 oz/900mL)	Single Unit

### Urinary Leg Bag Combination Pack

Includes: (1) Urinary Leg Bag – medium or large, (1) Latex extension tubing with connector, (1) Pair fabric leg bag straps – medium (23")

Stock No	Size	Box Qty
9824	Medium (18 oz/540mL)	Box of 12
9825	Large (30 oz/900mL)	Box of 12

Hollister and logo, Hollister and m9 are trademarks of Hollister Incorporated.  
©2005 Hollister Incorporated. 907267-405

## Latex-Free Urinary Drainage Bags

### Bedside Drainage Collection System with Anti-Reflux Valve

Includes: 68 oz/2000mL capacity, sample port, sterile fluid pathway

Stock No	Size	Box Qty
9839	68 oz/2000mL	Box of 20

## Accessories

### Fabric Leg Bag Straps

Stock No	Size	Box Qty
9342	Small 15 in/38cm	10 Pair
9343	Medium 23 in/58cm	10 Pair
9344	Large 30 in/76cm	10 Pair

### Latex-Free Extension Tubing

Stock No	Description	Box Qty
9345	Non-sterile	Box of 10
9346	Sterile	Box of 10

### Latex Extension Tubing

Stock No	Description	Box Qty
9803	Non-sterile	Box of 10
9813	Sterile	Box of 10

### m9 Odor Eliminators and Cleaner/Decrystallizer

Odor Eliminators – pump spray bottle

Stock No	Size/Description	Box Qty
7732	2 oz/60mL – unscented	12
7733	8 oz/240mL – unscented	6
7734*	2 oz/60mL – scented	12
7735	8 oz/240mL – scented	6

Cleaner/Decrystallizer – free wash bottle included

Stock No	Size/Description	Box Qty
7736	16 oz/480mL bottle	12

\*Not available in Canada



**Hollister Incorporated**  
2000 Hollister Drive  
Libertyville, Illinois 60048

Distributed in Canada by  
**Hollister Limited**  
95 Mary Street  
Aurora, Ontario L4G 1G3

**For a free sample or more information, call Hollister at 1-800-323-4060. In Canada, call 1-800-263-7400.**

[www.hollister.com](http://www.hollister.com)