

Fecal Collector



Fecal collectors from Hollister are designed to protect perianal skin and contain stool. They provide an effective fecal containment alternative to absorbent products for acutely ill patients.

For the Patient:

- Helps protect perianal skin from irritation and maceration caused by fecal discharge
- Helps contain unpleasant odors
- Improves patient dignity

For the Clinician:

- Easy-to-apply cut-to-fit skin barrier
- Reduces nursing time compared to frequent changing of absorbent products
- Attaches easily to a collection device for continuous drainage



Attention to Detail. Attention to Life.

Fecal Collectors



9880

<i>Description</i>	<i>Stock No</i>	<i>Quantity/Box</i>	<i>Pouch Length</i>
Large	9880	10	12" (30 cm)

- Cut-to-Fit FlexTend Tapered Skin Barrier
- Odor-Barrier Pouch Film
- Convenience Drain

HCPCS: A4330



9822/9821

<i>Description</i>	<i>Stock No</i>	<i>Quantity/Box</i>	<i>Pouch Length</i>
Medium	9822	10	10" (25 cm)
Large	9821	10	12" (30 cm)

- Cut-to-Fit SoftFlex Skin Barrier
- Odor-Barrier Pouch Film
- Thermometer Port
- Convenience Drain

HCPCS: A4330



9823

<i>Description</i>	<i>Stock No</i>	<i>Quantity/Box</i>	<i>Pouch Length</i>
Medium	9823	10	10" (25 cm)

- Cut-to-Fit FlexTend Tapered Skin Barrier
- Ultra-Clear Pouch Film
- Convenience Drain

HCPCS: A4330

Visit http://www.hollister.com/us/bm/fc/how_to_use/ to view a complete instructional video



Hollister Incorporated
 2000 Hollister Drive
 Libertyville, Illinois 60048 USA
 1.888.740.8999

Distributed in Canada by
Hollister Limited
 95 Mary Street
 Aurora, Ontario L4G 1G3
 1.800.263.7400

www.hollister.com