

# Tube Attachment Devices



Tube Attachment Devices from Hollister help protect the skin and provide an easy-to-use alternative to tape for securing tubes, drains, and catheters.

## Each device is designed to:

- Protect your patients' skin with a skin-friendly barrier that adheres securely, conforms easily, and eliminates the need for tape and sutures
- Allow better access for care with adjustable clamping mechanisms and easy tube positioning
- Hold drains, tubes, and catheters securely to help minimize tube migration and occlusion



Attention to Detail. Attention to Life.

## Tube Attachment Devices



### Horizontal Tube Attachment Device (HTAD)

<i>Description</i>	<i>Stock No</i>	<i>Quantity/Box</i>
	9781	5

- Secures drains, tubes, or catheters (5-40 Fr) parallel to patient's skin
- Secures tubes at exit site or distal to exit site
- FlexWear Skin Barrier
- Sterile, latex-free



### Vertical Tube Attachment Device (VTAD)

<i>Description</i>	<i>Stock No</i>	<i>Quantity/Box</i>
	9782	5

- Secures drains, tubes, or catheters (5-40 Fr) vertical to patient's skin
- Secures tubes at exit site
- FlexWear Skin Barrier
- Sterile, latex-free



### Drain Tube Attachment Device (DTAD)

<i>Description</i>	<i>Stock No</i>	<i>Quantity/Box</i>
	9780	5

- Secures a variety of drains or tubes (10-24 Fr) vertical to patient's skin
- FlexWear Skin Barrier
- Sterile, Latex-free



### Feeding Tube Attachment Device (FTAD)

<i>Description</i>	<i>Stock No</i>	<i>Quantity/Box</i>
	9786	20

- Secures nasogastric/nasointestinal tubes (5-18 Fr)
- Butterfly-shape FlexWear Skin Barrier conforms to patient's nose
- Unique swivel clamp eliminates pressure on the nares
- Latex-free



### Suction Tube Attachment (STAD)

<i>Description</i>	<i>Stock No</i>	<i>Quantity/Box</i>
	9785	20

- Secures nasogastric suction tubes (12-18 Fr)
- Butterfly-shape FlexWear Skin Barrier conforms to patient's nose
- Unique swivel clamp eliminates pressure on the nares
- Latex-free

Visit <http://www.hollister.com/us/tube/learning/theatre.asp> to view a complete instructional video



**Hollister Incorporated**  
 2000 Hollister Drive  
 Libertyville, Illinois 60048 USA  
 1.888.740.8999

Distributed in Canada by  
**Hollister Limited**  
 95 Mary Street  
 Aurora, Ontario L4G 1G3  
 1.800.263.7400

[www.hollister.com](http://www.hollister.com)