


Tube attachment devices from Hollister provide an easy-to-use alternative to tape for securing tubes, drains, and catheters. Each device features a skin-friendly FlexWear skin barrier to help protect patient skin and an adjustable clamping mechanism to allow for easy tube positioning and better access for site care.



Stock No	Description	Quantity/Box
9781	<b>Horizontal Tube Attachment Device (HTAD)*</b> <ul style="list-style-type: none"> <li>• Secures drains, tubes, or catheters (5-40 Fr) parallel to patient's skin</li> <li>• Secures tubes at exit site or distal to exit site</li> <li>• Sterile, latex-free</li> </ul>	5



Stock No	Description	Quantity/Box
9782	<b>Vertical Tube Attachment Device (VTAD)*</b> <ul style="list-style-type: none"> <li>• Secures drains, tubes, or catheters (5-40 Fr) vertical to patient's skin</li> <li>• Secures tubes at exit site</li> <li>• Sterile, latex-free</li> </ul>	5



Stock No	Description	Quantity/Box
9780	<b>Drain Tube Attachment Device (DTAD)</b> <ul style="list-style-type: none"> <li>• Secures a variety of drains or tubes (10-24 Fr) vertical to patient's skin</li> <li>• Sterile, Latex-free</li> </ul>	5



Stock No	Description	Quantity/Box
9786	<b>Feeding Tube Attachment Device (FTAD)</b> <ul style="list-style-type: none"> <li>• Secures nasogastric/nasointestinal tubes (5-18 Fr)</li> <li>• Butterfly-shaped skin barrier conforms to patient's nose</li> <li>• Unique swivel clamp eliminates pressure on the nares</li> <li>• Latex-free</li> </ul>	20



Stock No	Description	Quantity/Box
9785	<b>Suction Tube Attachment (STAD)</b> <ul style="list-style-type: none"> <li>• Secures nasogastric suction tubes (12-18 Fr)</li> <li>• Butterfly-shaped skin barrier conforms to patient's nose</li> <li>• Unique swivel clamp eliminates pressure on the nares</li> <li>• Latex-free</li> </ul>	20

\*Visit <http://www.hollister.com/us/tube/learning/theatre.asp> to view video instructions