



 **Obstetrics**
Product Portfolio



Attention to Detail. Attention to Life.

Double-Grip Umbilical Cord Clamp



AmniHook Amniotic Membrane Perforator



The AmniHook safely “snags” the membrane, eliminating the need to poke or puncture.

Double-Grip Umbilical Cord Clamp

Since 1964, the **Double-Grip** Umbilical Cord Clamp has given over 100 million babies a healthy, safe start in life! This precision-engineered cord clamp has earned its reputation for perfection. Our innovative laser vision inspection system verifies that only clamps meeting 100% of our criteria will make their way to the Labor and Delivery room.

- Safe, disposable, easy-to-apply.
- Dual gripping surface prevents slipping during application.
- Maintains tight seal and applies constant pressure on any size cord as it shrinks.
- Hinge guard confines cord between serrated jaws to eliminate movement.
- No rough edges to snag clothing or get caught on bed linen.

Description	Unit of Sale	Stock No
Sterile Cord Clamp	125 each/box	9423
Non-Sterile Cord Clamp	250 each/box	9411

AmniHook Amniotic Membrane Perforator

Our **AmniHook** Amniotic Membrane Perforator provides physicians with a safe, comfortable method to rupture the amniotic membrane the first time, virtually every time. Each AmniHook must pass through our rigorous, laser vision inspection system to ensure the first AmniHook will be the only one required for the procedure.

- Rounded tip and protected point designed to rupture the amniotic membrane without injury to mother or fetus.
- Disposable, pre-sterilized, individually packaged for convenience.
- 10" long, with molded sides ensure a non-slip grip and precise control.

Description	Unit of Sale	Stock No
Sterile Membrane Perforator	100 each/box	9601

Cord Clamp Clipper

Designed to work with our **Double-Grip** Umbilical Cord Clamp, our Clipper safely cuts through the specially designed clamp hinge for easy removal from baby.

- Concealed cutting surface eliminates risk of sharp or pointed objects near baby.
- Single patient use only to prevent risk of cross-contamination.

Description	Unit of Sale	Stock No
Single Clipper	1 each/box	9441
Bulk Clippers	25 each/box	9445



Cord Clamp Clipper





PlastiBell Circumcision Device

In today's fast paced medical environment, a circumcision procedure requiring no special post operative care or dressing means a significant savings in both time and money.

The disposable **PlastiBell** Circumcision Device eliminates the need for repetitive and costly sterilization required of stainless steel clamps. No chance of lost or mismatched parts that often render other clamps unsafe or unusable.

- Clear, plastic construction provides visual access at all times.
- A clean line of excision; helps to promote rapid healing.
- No special dressings required.
- Sterile, individually packaged in six sizes and assortment packs.

Description	Unit of Sale	Stock No
Hospital Assortment Pack	125 each/box	9210
1.1 cm (25 each)		
1.2 cm (35 each)		
1.3 cm (45 each)		
1.4 cm (10 each)		
1.5 cm (5 each)		
1.7 cm (5 each)		

Description	Unit of Sale	Stock No
Physician Office Assortment Pack	25 each/box	9230

- | | | |
|-----------------|--|--|
| 1.1 cm (5 each) | | |
| 1.2 cm (7 each) | | |
| 1.3 cm (8 each) | | |
| 1.4 cm (3 each) | | |
| 1.5 cm (1 each) | | |
| 1.7 cm (1 each) | | |

Hospital (Bulk) Single Sizes

1.1 cm	125 each/box	9211
1.2 cm	125 each/box	9212
1.3 cm	125 each/box	9213
1.4 cm	125 each/box	9214
1.5 cm	125 each/box	9215
1.7 cm	125 each/box	9217

Physician Office Single Sizes

1.1 cm	25 each/box	9231
1.2 cm	25 each/box	9232
1.3 cm	25 each/box	9233
1.4 cm	25 each/box	9234
1.5 cm	25 each/box	9235
1.7 cm	25 each/box	9237

Parent/Physician Education

Parent Education Pamphlets

Shows parents what to expect as healing occurs.

Physician Education Video

Demonstrates infant circumcision procedure with the PlastiBell Circumcision Device.

Parent/Physician Education

Description	Unit of Sale	Stock No
Parent Education Pamphlets	25 each/pack	905087
Physician Education Video (VHS)	1 each	905407



U-Bag Urine Specimen Collectors: (above left) Newborn; (above center) Premie; (right) 24-hour Style



Patented “bump” simplifies placement on the infant and helps reduce incidence of fecal contamination.

(right) First-Time Urine Specimen Collectors



U-Bag Urine Specimen Collectors

The **U-Bag** Urine Specimen Collectors are the first choice in convenience, protection and care for an infant's tender skin. This complete system of urine collectors features products to facilitate any infant urine collection need.

- Transparent butterfly patch makes application easy and accurate.
- Compatible and comfortable for male and female anatomy.
- Double chambered construction with no-flow back valve keeps urine away from skin to help prevent irritation.
- Adhesive patch can be folded in half to send sealed specimen to the lab.
- Premie size ideal for NICU use with low birth weight infants.

Description	Unit of Sale	Sterile Stock No	Non-Sterile Stock No
U-Bag with Porous Cloth Adhesive			
Single Specimen Style			
Premie (10 mL)	100 each/box	7540	NA
Newborn (100 mL)	100 each/box	7535	7525
Pediatric (200 mL)	100 each/box	7531	7521
24-hour Style			
Newborn (100 mL)	50 each/box	7537	7527
Pediatric (200 mL)	50 each/box	7533	7523

U-Bag with Standard Adhesive

Single Specimen Style			
Newborn (100 mL)	100 each/box	7515	7505
Pediatric (200 mL)	100 each/box	7511	7501
24-hour Style			
Newborn (100 mL)	100 each/box	7517	7507
Pediatric (200 mL)	100 each/box	7513	7503

First-Time Urine Specimen Collectors


The **First-Time** Urine Specimen Collectors were designed with extensive input from nursery and pediatric nurse clinicians. Use to obtain adequate urine specimens whenever a collector is applied.

- Features patented “bump” design for even simpler placement and reduction of fecal contamination.
- Easy to use drain spout with attached, replaceable cap at base of collector.
- Accommodates syringes.
- Compatible and comfortable with male and female anatomy.
- Special foam adhesive is gentle on baby's sensitive skin.
- Double chambered construction with no-flow back valve traps urine away from skin to help prevent irritation.

Description	Unit of Sale	Stock No
Sterile Single Specimen Style		
Newborn (100 mL)	25 each/box	7545
Pediatric (200 mL)	25 each/box	7541
Sterile 24-hour Style		
Newborn (100 mL)	10 each/box	7547
Pediatric (200 mL)	10 each/box	7543



(left) Birth Certificates;
(below) Disposable Footprinter;
(right) Newborn Identification Forms



Disposable FootPrinter

A baby's footprints can provide a positive means of newborn identification, but only if the prints are clear and readable. Our disposable **Hollister FootPrinters** are designed to produce perfect prints every time. Over-inking is virtually impossible with our unique **Ready-Rolled** ink process that allows for just enough ink for two baby footprints and one mother's fingerprint. These footprinters require such minimal space that 1,250 plates can be stored in a single cubic foot.

- Disposable plastic plate fits easily in the palm of one hand.
- 4-1/2" x 2-3/8", less than 1/8" thick with stackable edges for easy, compact storage.
- Non-toxic, black ink – easy to see and wash off of baby.
- Great for official identification forms, **Hollister Birth Certificates** and hospital made mementos.

Description	Unit of Sale	Stock No
(25) Plastic Footprinters/box	4 boxes/case	5402

Hollister Birth Certificates

For over 80 years **Hollister Birth Certificates** have been welcoming babies, preserving memories and bringing joy to families. It's customary to celebrate the birth of a baby by presenting a gift to the parents. What better than one that captures the details of the very event!

- Choose from three styles and customize to your own specifications.
- All certificates are personalized with the hospital name, city, and state.
- Gold seals and pink or blue satin ribbons also available.

Call our **Customer Service Center** at **1.800.323.4060** for additional information on options available for your hospital.



HOLLISTER INCORPORATED **OBSTETRICS**

Newborn Identification Forms

The paper used to record a footprint has a great deal to do with the clarity and overall quality of the prints. Our Hollister **Newborn Identification Forms** are pre-printed on high contrast, archival quality paper, perfect for both receiving and preserving infant footprints. When used with the disposable **Hollister FootPrinter**, the Newborn Identification Forms assure the hospital of dependable, readable newborn identification to link mother and baby.

- Receives and retains every detail of a baby's tiny skin print.
- Space for recording pertinent birth information and verifying signatures.
- Enough space for two infant footprints and a mother's fingerprint.
- Some forms allow for the attachment of infant's ID band.

Description	Unit of Sale	Stock No
Minimum Order is 5 pads		
8" x 5" Pertinent Birth Information, Infant Footprints, Mother's Fingerprint	50 sheets/pad	5821
8" x 11" Same as 5821, plus signatures for ID band verification	50 sheets/pad	5823
8" x 11" Same as 5823, plus space to attach infant ID band upon discharge	50 sheets/pad	5825
8" x 11" Same as 5825, for use with mechanical imprinters	50 sheets/pad	5826



Hollister Incorporated is committed to improving the quality of life and adding value to healthcare worldwide. Through every line of our company's products including MotherBaby Care, Ostomy, Urologicals and Wound Care – the overriding goal and commitment is to make a difference – to add value and improve the quality of life for the people who use our products. Our founder, John Dickinson Schneider, established his company in Chicago in 1921 as a printing shop later called JDS Printer Craftsman, Inc. The company was built on a solid reputation for quality and service. His underlying principle was “Only First Class is Good Enough” and that principle carried over into a philosophy of quality products, quality service and quality people.

The Hollister Traditional Obstetrics Products have repeatedly assured medical professionals of the best quality of care for the families they serve. All of Hollister's products must go through world class quality inspection which virtually guarantees a perfect product every time. Hollister has been in the Labor and Delivery suite for over 80 years and has assisted in the birth of well over 100 million babies.

Hollister's commitment to safety and quality continues to be trusted each and every day by medical professionals worldwide and is unrivaled by other brands of traditional obstetric products. So we ask, “Why take a risk?” With Hollister Traditional Obstetric Products, you can count on the quality we have stood behind for years.

At Hollister, innovation links past, present and future. Hollister has set the highest standards with a continued commitment to build, serve and grow.

Obstetrics Products

- Double-Grip Umbilical Cord Clamp
- AmniHook Amniotic Membrane Perforator
- U-Bag Urine Collection Bags
- PlastiBell Circumcision Device
- Disposable FootPrinter
- Newborn Identification Forms
- Hollister Birth Certificates

Availability Notice

This product catalog represents products that may not be available in all countries. Check with your local distributor for availability.



Hollister Incorporated
2000 Hollister Drive
Libertyville, Illinois 60048 USA
1.800.323.4060

Distributed in Canada by
Hollister Limited
95 Mary Street
Aurora, Ontario L4G 1G3
1.800.263.7400

www.hollister.com