

# Continence Care

Continence  
Care

### HCPCS Codes Disclaimer

The reimbursement information provided herein is intended to provide general information concerning coding of Hollister products only. Hollister does not guarantee coverage or payment for any products. The ultimate responsibility for proper coding, satisfying reimbursement requirements, and obtaining reimbursement remains with the provider. Coding and coverage policies and guidelines are complex, can vary from one carrier or region to another, and are updated frequently. Providers should check with their local carriers or intermediaries often and should consult with counsel, a reimbursement specialist, and/or the DMECS website for any coding, coverage, reimbursement, or billing questions.

*HCPCS codes shown in this Product Catalog are current as of 1-6-12.*

### Explanation of Third Party Code Reimbursement System (National Health Related Items Code)

To facilitate reimbursement of Hollister products, some Third Party Payment Programs require the use of an eleven-digit reimbursement number. Every Hollister manufactured product can be identified by an eleven-digit number. In general, the first five numbers identify the manufacturer's labeler code, and the last six numbers represent the product stock number. The labeler code for Hollister Incorporated is 08380. The labeler code for Hollister Wound Care LLC is 08567, and the labeler code for the Apogee product line is 08598. In order to utilize this coding system, follow these examples.

#### To code a four-digit stock number:

The Third Party Reimbursement Number for the Apogee 12 Fr Intermittent Catheter, Stock Number 1048 listed in the catalog on page 92, is: 08598-0010-48

08598 identifies Apogee. 0010-48 identifies the product Stock Number.

#### To code a five-digit stock number:

The Third Party Reimbursement Number for the New Image Drainable Pouch, Stock Number 18113 listed in the catalog on page 30, is: 08380-0181-13

08380 identifies Hollister Incorporated. 0181-13 identifies the product Stock Number.

#### To code a six-digit stock number:

The Third Party Reimbursement Number for the Restore Contact Layer Dressing with TRIACT Technology, Stock Number 509338 listed in the catalog on page 130, is:

08567-5093-38

08567 identifies Hollister Wound Care LLC. 5093-38 identifies the product Stock Number.

Additionally, three products in the Hollister Wound Care line require an NDC number due to their Dimethicone-based nature. The labeler code for their OTC drug products is 42828. The three NDC numbers are listed below.

NDC	Stock No.	Description
42828-9977-02	529977	Restore Cleanser & Moisturizer, 7.5 oz
42828-9978-03	529978	Restore Cleanser & Moisturizer, 11.5 oz
42828-9979-01	529979	Restore Dimethicreme Skin Protectant, 4 oz

## Latex Statement

The specifications for the Hollister Continence Care products listed below include natural rubber latex as a component; our suppliers add natural rubber latex as part of the production process.

If a Hollister product(s) found in this catalog is not listed below, the specification for that product(s) and packaging does not include natural rubber latex as a component, and our component suppliers do not add natural rubber as part of their production process.

## Current Stock Numbers:

1051	9107	9207	9835 §	B12RB
1067	9108	9208	B12R	B14RB
1084	9109	9209	B14R	B16RB
9106	9206	9834 §	B16R	

## Discontinued Stock Numbers:

9606 *	9628 *	9656 †	9807 *	9818 †	9829 *
9607 *	9629 *	9657 †	9808 *	9824 ‡	9897 †
9608 *	9636 †	9658 †	9809 *	9825 ‡	9898 †
9609 *	9637 †	9659 †	9813 ‡	9826 *	
9626 *	9638 †	9803 ‡	9816 †	9827 *	
9627 *	9639 †	9806 *	9817 †	9828 *	

## How to Use The Hollister Products Catalog

Hollister has created this simplified Product Catalog to help our customers find what they need in a quick and efficient manner.

### *If you...*

*know the stock number (or stock series number)*

Go to the numeric index at the back of the catalog on pages 162–163 to find the page number where it appears.

### *If you...*

*know the product brand or series name*

Go to the Table of Contents on pages 9-10, 86, 106, 114 and 120 to find the page number where it appears or reference the Product Name Index on pages 164–166.

### *If you...*

*know the HCPCS code*

Go to the HCPCS Code Index at the back of the catalog on page 167.

---

† Product Discontinued as of July 2007

‡ Product Discontinued as of January 2009

\* Product Discontinued as of September 2009

§ As of February 2009, this product has been redesigned to exclude natural rubber latex either as a component or from the production process. The packaging/insert of the new non-latex version of the product will not contain a latex statement. Please verify the latex status of the individual product prior to use by reviewing the packaging and product insert for a latex statement. A latex statement indicates the product is from the previous latex version.

# Hollister Continence Care



## People First.

Hollister Continence Care is committed to people and to helping empower their lives. Our products and services are testimony – first and foremost – to the assurance that quality of life needn't be compromised by managing one's continence.



Intermittent Catheters



Male External Catheters



Leg Bag Systems

[www.hollister.com](http://www.hollister.com)

# Table of Contents • Continance Care

<b>General Information</b> .....	<b>87</b>
HCPCS Codes Disclaimer.....	87
Explanation of Third Party Code Reimbursement System..... (National Health Related Items Code)	87
<b>Medicare Utilization Chart</b> .....	<b>88–89</b>
<b>Intermittent Catheters</b> .....	<b>92–94</b>
<b>Apogee</b> Intermittent Catheter.....	92
<b>Apogee Plus</b> System .....	93
<b>Apogee Plus</b> System Kit .....	93
<b>Advance Plus</b> Touch Free Intermittent Catheter.....	94
<b>Advance Plus</b> Touch Free Intermittent Catheter System.....	94
<b>Male External Catheters</b> .....	<b>95</b>
<b>InView</b> Male External Catheters.....	95
<b>Extended Wear</b> Male External Catheter.....	95
<b>Everyday</b> Male External Catheter.....	95
<b>Leg Bag Systems</b> .....	<b>96</b>
Urinary Leg Bag.....	96
Urinary Leg Bag Combination Pack .....	96
Vented Urinary Leg Bag Combination Pack.....	96
<b>Accessories</b> .....	<b>97</b>
Fabric Leg Bag Straps .....	97
Latex-Free Extension Tubing.....	97
<b>Collection Systems</b> .....	<b>98</b>
Bedside Drainage Collection System with Anti-Reflux Valve.....	98
Male Urinary Pouch External Collection Device.....	98
Female Urinary Pouch External Collection Device .....	98
<b>Glossary • Continance Care</b> .....	<b>99</b>

### HCPCS Codes Disclaimer

The reimbursement information provided herein is intended to provide general information concerning coding of Hollister products only. Hollister does not guarantee coverage or payment for any products. The ultimate responsibility for proper coding, satisfying reimbursement requirements, and obtaining reimbursement remains with the provider. Coding and coverage policies and guidelines are complex, can vary from one carrier or region to another, and are updated frequently. Providers should check with their local carriers or intermediaries often and should consult with counsel, a reimbursement specialist, and/or the DMECS website for any coding, coverage, reimbursement, or billing questions.

***HCPCS codes shown in this Product Catalog are current as of 1-6-12.***

### Explanation of Third Party Code Reimbursement System (National Health Related Items Code)

To facilitate reimbursement of Hollister products, some Third Party Payment Programs require the use of an eleven-digit reimbursement number. Every Hollister manufactured product can be identified by an eleven-digit number. In general, the first five numbers identify the manufacturer's labeler code, and the last six numbers represent the product stock number. The labeler code for Hollister Incorporated is 08380. The labeler code for Hollister Wound Care LLC is 08567, and the labeler code for the Apogee product line is 08598. In order to utilize this coding system, follow these examples.

#### To code a four-digit stock number:

The Third Party Reimbursement Number for the Apogee 12 Fr Intermittent Catheter, Stock Number 1048 listed in the catalog on page 92, is:  
08598-0010-48

08598 identifies Apogee. 0010-48 identifies the product Stock Number.

#### To code a five-digit stock number:

The Third Party Reimbursement Number for the Advance Plus Intermittent Catheter Kit 12 Fr, straight tip, Stock Number 96124 listed in the catalog on page 94 is: 08380-0961-24

08380 identifies Hollister Incorporated. 0961-24 identifies the product Stock Number.

For more information regarding how to code the Apogee Closed System products with the alpha-numeric Stock Numbers, please contact Hollister Incorporated at 1.888.740.8999.

## Medicare's Usual Maximum Quantity for Urological Supplies

The quantity of urological supplies needed by a patient is determined by many factors. There will be variation according to individual need and their needs may vary over time. The table below lists the maximum number of items/units of service that are medically necessary. The actual quantity needed for a particular patient may be more or less than the amount listed. The medical necessity for use of a greater quantity of supplies than the amounts listed must be clearly documented in the patient's medical record, and must be available upon request.

**Per LCD for Urological Supplies (L5080), NHIC effective Feb. 4, 2011**

HCPCS	Long Description	Quantity
A4217	Sterile water/saline, 500 mL	*
A4310	Insert tray without drainage bag and without catheter (accessories only)	1 tray/insertion
A4311	Insertion tray without drainage bag with indwelling catheter, foley type, two-way latex with coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.)	1 tray/insertion
A4312	Insertion tray without drainage bag with indwelling catheter, foley type, two-way all silicone	1 tray/insertion
A4313	Insertion tray without drainage bag with indwelling catheter, foley type, three-way, for continuous irrigation	1 tray/insertion
A4314	Insertion tray with drainage bag with indwelling catheter, foley type, two-way latex with coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.)	1 ea/mo
A4315	Insertion tray with drainage bag with indwelling catheter, foley type, two-way, all silicone	1 ea/mo
A4316	Insertion tray with drainage bag with indwelling catheter, foley type, three-way, for continuous irrigation	1 ea/mo
A4320	Irrigation tray with bulb or piston syringe, any purpose	*
A4321	Therapeutic agent for urinary catheter irrigation	*
A4322	Irrigation syringe, bulb, or piston, each	*
A4326	Male external catheter with integral collection chamber, any type, each	*
A4327	Female external urinary collection device; meatal cup, each	1 ea/wk
A4328	Female external urinary collection device; pouch, each	1 ea/day
A4331	Extension drainage tubing, any type, any length, with connector/adaptor, for use with urinary leg bag or urostomy pouch, each	*
A4332	Lubricant, individual sterile packet, each	200 ea/mo
A4333	Urinary catheter anchor device, adhesive skin attachment, each	3 ea/wk
A4334	Urinary catheter anchoring device, leg strap, each	1 ea/mo
A4335	Incontinence supply; miscellaneous	*
A4336	Incontinence supply, urethral insert, any type, each	*
A4338	Indwelling catheter; foley type, two-way latex with coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.), each	1 ea/mo
A4340	Indwelling catheter; specialty type, (e.g., coudé, mushroom, wing, etc.), each	1 ea/mo
A4344	Indwelling catheter, foley type, two-way, all silicone, each	1 ea/mo
A4346	Indwelling catheter; foley type, three-way for continuous irrigation, each	1 ea/mo
A4349	Male external catheter, with or without adhesive, disposable, each	35 ea/mo

\* No established usual maximum, refer to Medicare Utilization Guidelines.

HCPGS	Long Description	Quantity
A4351	Intermittent urinary catheter; straight tip, with or without coating (Teflon, silicone, silicone elastomer, or hydrophilic, etc.), each	200 ea/mo
A4352	Intermittent urinary catheter; coudé (curved) tip, with or without coating (Teflon, silicone, silicone elastomeric, or hydrophilic, etc.), each	200 ea/mo
A4353	Intermittent urinary catheter, with insertion supplies	200 ea/mo
A4354	Insertion tray with drainage bag but without catheter	1 ea/mo
A4355	Irrigation tubing set for continuous bladder irrigation through a three-way indwelling foley catheter, each	1 ea/day
A4356	External urethral clamp or compression device (not to be used for catheter clamp), each	1 ea/3 mo
A4357	Bedside drainage bag, day or night, with or without anti-reflux device, with or without tube, each	2 ea/mo
A4358	Urinary drainage bag, leg or abdomen, vinyl, with or without tube, with straps, each	2 ea/mo
A4360	Disposable external urethral clamp or compression device, with pad and/or pouch, each	NC
A4402	Lubricant, per ounce	*
A4450	Tape, non-waterproof, per 18 square inches	10 units (18 sq in)/mo
A4452	Tape, waterproof, per 18 square inches	10 units (18 sq in)/mo
A4455	Adhesive remover or solvent (for tape, cement, or other adhesive), per ounce	*
A4520	Incontinence garment, any type (e.g., brief, diaper), each	*
A4456	Adhesive remover, wipes, any type, each	NC
A4554	Disposable underpads, all sizes	*
A5102	Bedside drainage bottle with or without tubing, rigid or expandable, each	1 ea/3 mo
A5105	Urinary suspensory with leg bag, with or without tube, each	*
A5112	Urinary drainage bag, leg or abdomen, latex, with or without tube, with straps, each	1 ea/mo
A5113	Leg strap; latex, replacement only, per set	*
A5114	Leg strap, foam or fabric, replacement only, per set	*
A5131	Appliance cleaner, incontinence and ostomy appliances, per 16 ounce	1 (16 oz)/mo
A5200	Percutaneous catheter/tube anchoring device, adhesive skin attachment	*
A9270	Non-covered item or service	

\* No established usual maximum, refer to Medicare Utilization Guidelines.

Note: Intermittent catheterization using a sterile intermittent catheter kit (A4353) is covered when the patient requires catheterization and the patient meets one of the following criteria (1-5):

- 1) The patient resides in a nursing facility
- 2) The patient is immunosuppressed, for example (not all-inclusive): on a regimen of immunosuppressive drugs post-transplant, on cancer chemotherapy, has AIDS, has a drug-induced state such as chronic oral corticosteroid use
- 3) The patient has radiologically documented vesico-urethral reflux while on a program of intermittent catheterization
- 4) The patient is a spinal-cord-injured-female with neurogenic bladder who is pregnant (for duration of pregnancy only)
- 5) The patient has had distinct, recurrent urinary tract infections, while on a program of sterile intermittent catheterization with A4351/A4352 and sterile lubricant A4332, twice within the 12 months prior to the initiation of sterile intermittent catheter kits.

## Advance Plus

Touch Free Intermittent Catheter System



*Confidence. Pure and simple.*

Hollister Advance Plus Intermittent Catheter features *Touch Free* technology enabling users to catheterize with confidence – anytime, anywhere – supporting user mobility and independence.

A unique *Touch Free* protective introducer tip helps reduce the risk of carrying bacteria into the urinary system. The Advance Plus *Touch Free* collection bag shields the sterile catheter from possible environmental contamination. Ultra-smooth catheter eyelets, together with a unique gel reservoir, help ensure trouble-free insertion and withdrawal for enhanced user convenience and comfort.

[www.hollister.com](http://www.hollister.com)

**Information for Use:** Intermittent Catheters are flexible tubular devices inserted through the urethra for the purpose of bladder drainage. Self-catheterization should only be carried out under medical advice and only in accordance with the instructions provided with the product. Should a patient experience any sign of trauma or discomfort, they should immediately discontinue use and consult their clinician. **Caution:** Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional. **Date of Issuance:** January 2012.

# Intermittent Catheters

## Apogee Intermittent Catheter



### 10xx

- ▶ Ultra-Smooth Eyelets
- ▶ Latex-Free (exceptions noted below)
- ▶ Designed for Comfort

**Box of 30**

**Box of 50**

**HCPCS: A4351**

**HCPCS for 1021, 1035, 1050, 1066, 1086, 1087 & 1097:  
A4352**

15" Length Stock No	Size	Description	Quantity/Box
1007	6 Fr	Firm	30

16" Length Stock No	Size	Description	Quantity/Box
1018	8 Fr	Firm	30
1034	10 Fr	Firm	30
1048	12 Fr	Firm	30
1064	14 Fr	Firm	30
1082	16 Fr	Firm	30
1096	18 Fr	Firm	30
1036	10 Fr	Soft	30
1046	12 Fr	Soft	30
1061	14 Fr	Soft	30
1083	16 Fr	Soft	30
1049	12 Fr	Curved Pkg	50
1065	14 Fr	Curved Pkg	50
1074A	14 Fr	Soft, Curved Pkg	50
1050	12 Fr	Coudé	30
1066	14 Fr	Coudé	30
1086	16 Fr	Coudé	30
1021	8 Fr	Coudé, Olive Tip	30
1035	10 Fr	Coudé, Olive Tip	30
1097	18 Fr	Coudé, Olive Tip	30
1087	16 Fr	Coudé, Curved Pkg	50
1051	12 Fr	Red Rubber Latex*	30
1067	14 Fr	Red Rubber Latex*	30
1084	16 Fr	Red Rubber Latex*	30

10" Length Stock No	Size	Description	Quantity/Box
1020	8 Fr	Firm, Luer End	30
1040	10 Fr	Firm, Luer End	30

6" Length Stock No	Size	Description	Quantity/Box
1016	8 Fr	Firm	30
1033	10 Fr	Firm	30
1047	12 Fr	Firm	30
1062	14 Fr	Firm	30
1063	14 Fr	No Funnel	30

**\*Caution:** This product contains natural rubber latex which may cause allergic reactions.

## Apogee Plus System



### Bxxx

Case of 100

- ▶ 16" length
- ▶ Introducer Tip
- ▶ Pre-lubricated
- ▶ Ultra-Smooth Eyelets
- ▶ 1500 mL Collection Bag
- ▶ Latex-Free (exceptions noted below)

HCPCS: A4353

Stock No	Size	Description
B6F	6 Fr	Firm
B8F	8 Fr	Firm
B10F	10 Fr	Firm
B12F	12 Fr	Firm
B14F	14 Fr	Firm
B16F	16 Fr	Firm
B18F	18 Fr	Firm
B12C	12 Fr	Coudé
B14C	14 Fr	Coudé
B16C	16 Fr	Coudé
B12S	12 Fr	Soft
B14S	14 Fr	Soft
B12R	12 Fr	Red Rubber Latex*
B14R	14 Fr	Red Rubber Latex*
B16R	16 Fr	Red Rubber Latex*
B14F-Female	14 Fr, 6"	Firm

**\*Caution:** This product contains natural rubber latex which may cause allergic reactions.

## Apogee Plus System Kit



### BxxxB

Case of 100

- ▶ 16" length
- ▶ Introducer Tip
- ▶ Pre-lubricated
- ▶ Ultra-Smooth Eyelets
- ▶ 1500 mL Collection Bag
- ▶ Kit contains gloves, underpad drape, gauze, and PVP swabsticks
- ▶ Latex-Free (exceptions noted below)

HCPCS: A4353

Stock No	Size	Description
B6FB	6 Fr	Firm
B8FB	8 Fr	Firm
B10FB	10 Fr	Firm
B12FB	12 Fr	Firm
B14FB	14 Fr	Firm
B16FB	16 Fr	Firm
B18FB	18 Fr	Firm
B12CB	12 Fr	Coudé
B14CB	14 Fr	Coudé
B16CB	16 Fr	Coudé
B12SB	12 Fr	Soft
B14SB	14 Fr	Soft
B12RB	12 Fr	Red Rubber Latex*
B14RB	14 Fr	Red Rubber Latex*
B16RB	16 Fr	Red Rubber Latex*
B14FB-Female	14 Fr, 6"	Firm

**\*Caution:** This product contains natural rubber latex which may cause allergic reactions.

## Advance Plus Touch Free Intermittent Catheter



**940xx/941xx/951xx/931xx**

**Box of 100**

- ▶ 16" length
- ▶ Latex-Free
- ▶ Gel Reservoir
- ▶ Unique Protective Tip
- ▶ 1500 mL Collection Bag
- ▶ Ultra-Smooth Eyelets

**HCPCS: A4353**

Stock No	Size	Description	Length
94064	6 Fr	Straight	16"
94084	8 Fr	Straight	16"
94104	10 Fr	Straight	16"
94124	12 Fr	Straight	16"
94144	14 Fr	Straight	16"
94164	16 Fr	Straight	16"
94184	18 Fr	Straight	16"
95124	12 Fr	Coudé	16"
95144	14 Fr	Coudé	16"
95164	16 Fr	Coudé	16"
93104*	10 Fr	Straight	16"
93124*	12 Fr	Straight	16"
93144*	14 Fr	Straight	16"
93164*	16 Fr	Straight	16"

\***Advance Plus Pocket Intermittent Catheter**

## Advance Plus Touch Free Intermittent Catheter System



**960xx/961xx/971xx**

**Box of 100**

- ▶ Latex-Free
- ▶ Gel Reservoir
- ▶ Unique Protective Tip
- ▶ 1500 mL Collection Bag
- ▶ Ultra-Smooth Eyelets
- ▶ Kit contains gloves, underpad drape, and benzalkonium chloride swabsticks

**HCPCS: A4353**

Stock No	Size	Description	Length
96064	6 Fr	Straight	16"
96084	8 Fr	Straight	16"
96104	10 Fr	Straight	16"
96124	12 Fr	Straight	16"
96144	14 Fr	Straight	16"
96164	16 Fr	Straight	16"
96184	18 Fr	Straight	16"
97124	12 Fr	Coudé	16"
97144	14 Fr	Coudé	16"
97164	16 Fr	Coudé	16"

# Male External Catheters

## InView Male External Catheters



**974xx/975xx/976xx**  
**977xx/978xx/979xx**

**Box of 30**  
**Box of 100**

- ▶ Self-Adhesive
- ▶ Latex-Free
- ▶ Clear Material Enables Monitoring of Skin

**HCPCS: A4349**

### SPECIAL: shorter length sheath



Stock No Box of 30	Stock No Box of 100	Color	Size
97425	97725	Green	25 mm
97429	97729	Purple	29 mm
97432	97732	Brown	32 mm
97436	97736	Blue	36 mm
97441	97741	Orange	41 mm

### STANDARD: the daily catheter for general use



Stock No Box of 30	Stock No Box of 100	Color	Size
97525	97825	Green	25 mm
97529	97829	Purple	29 mm
97532	97832	Brown	32 mm
97536	97836	Blue	36 mm
97541	97841	Orange	41 mm

### EXTRA: 70% more adhesive than Standard for added security



Stock No Box of 30	Stock No Box of 100	Color	Size
97625	97925	Green	25 mm
97629	97929	Purple	29 mm
97632	97932	Brown	32 mm
97636	97936	Blue	36 mm
97641	97941	Orange	41 mm

## Extended Wear Male External Catheter



**920x**

**Box of 30**

- ▶ Self-Adhesive, Latex
- ▶ Unique Inner Flap
- ▶ Easy Applicator Collar

**HCPCS: A4349**

Stock No Box of 30	Color	Size
9206	Green	22 mm-25 mm
9207	Purple	26 mm-30 mm
9208	Brown	31 mm-35 mm
9209	Blue	36 mm-39 mm

**Caution:** This product contains natural rubber latex which may cause allergic reactions.

## Everyday Male External Catheter



**910x**

**Box of 30**

- ▶ Self-Adhesive, Latex
- ▶ Easy Applicator Collar

**HCPCS: A4349**

Stock No Box of 30	Color	Size
9106	Green	22 mm-25 mm
9107	Purple	26 mm-30 mm
9108	Brown	31 mm-35 mm
9109	Blue	36 mm-39 mm

**Caution:** This product contains natural rubber latex which may cause allergic reactions.

# Leg Bag Systems

## Urinary Leg Bag



**9814/9805**

**Box of 10**

- ▶ Includes:
  - (1) Urinary Leg Bag — medium or large
- ▶ Latex-Free

**HCPCS: A9270**

**HCPCS: A4358 (if straps are added)**

Stock No	Description	Length	Width
9814	Medium — 18 oz (540 mL)	10"	4-1/4"
9805	Large — 30 oz (900 mL)	11-1/8"	5"

## Urinary Leg Bag Combination Pack



**9348/9349**

**Box of 12**

- ▶ Includes:
  - (1) Urinary Leg Bag — medium or large
  - (1) Extension tubing (18" long, oval, 11/32" effective inside diameter) with connector
  - (1) Pair fabric leg bag straps — medium (23" long, 1-1/8" wide)
- ▶ Latex-Free; Contains sterile components

**HCPCS: A4358**

Stock No	Description	Length	Width
9348	Medium — 18 oz (540 mL)	10"	4-1/4"
9349	Large — 30 oz (900 mL)	11-1/8"	5"

## Vented Urinary Leg Bag Combination Pack



**9645/9655**

**Box of 10**

- ▶ Includes:
  - (1) Vented Urinary Leg Bag — medium or large
  - (1) Kink-resistant extension tubing with vacuum relief valve, 18" long, oval, 11/32" effective inside diameter
  - (1) Pair fabric leg bag straps — medium (23" long, 1-1/8" wide)
- ▶ Latex-Free; Contains sterile components

**HCPCS: A4358**

Stock No	Description	Length	Width
9645	Medium — 18 oz (540 mL)	10"	4-1/4"
9655	Large — 30 oz (900 mL)	11-1/8"	5"

## Did you know...

### Hollister Leg Bags Offer Many Benefits:

- Soft, vinyl material is comfortable
- Pleats provide a discreet profile
- Wide, woven fabric straps are comfortable and stay in place
- Easy-to-operate drain valve requires minimal dexterity
- Protected buttons on the bags may decrease the risk of pressure points
- Anti-reflux valve helps prevent urine back flow into the tubing

**Fabric Leg Bag Straps**



**9342/9343**

**10 Pair**

▶ Latex-Free

**HCPCS: A5114**

Stock No	Description	Length	Width
9342	Small	15"	1-1/8"
9343	Medium	23"	1-1/8"

**Latex-Free Extension Tubing**



**9345/9346**

**Box of 10**

▶ Latex-Free, kink-resistant extension tubing with connector, oval, 11/32" effective inside diameter

**HCPCS: A4331**

Stock No	Description	Length
9345	Non Sterile	18"
9346	Sterile	18"

## Bedside Drainage Collection System with Anti-Reflux Valve



**9839**

**Box of 20**

- ▶ Latex-Free
- ▶ Sample port
- ▶ Sterile fluid pathway

**HCPCS: A4357**

Stock No	Description	Length	Width
9839	68 oz (2000 mL)	12"	9-3/4"

## Male Urinary Pouch External Collection Device\*



**9873 (CUT-TO-FIT)**

**Box of 10**

- ▶ FlexTend Skin Barrier
- ▶ No Additional Tape Needed
- ▶ Odor-Barrier Film
- ▶ Tubing Adaptor Port
- ▶ Alternative to Indwelling Catheter for Bedridden Males
- ▶ Latex-Free

**HCPCS: A4326**

Stock No	Description	Length	Width
9873	Male Urinary Pouch	7-1/2"	4-1/2"



**9811 (CUT-TO-FIT)**

**Box of 10**

- ▶ SoftFlex Skin Barrier
- ▶ Thinner, More Flexible Barrier
- ▶ Odor-Barrier Film
- ▶ Tubing Adaptor Port
- ▶ Tape Strips and Skin Gel Wipe Included
- ▶ Alternative to Indwelling Catheter for Bedridden Males
- ▶ Latex-Free

**HCPCS: A4326**

Stock No	Description	Length	Width
9811	Male Urinary Pouch	7-1/2"	4-1/2"

## Female Urinary Pouch External Collection Device\*



**9840 (CUT-TO-FIT)**

**Box of 10**

- ▶ SoftFlex Skin Barrier
- ▶ Odor-Barrier Film
- ▶ Tubing Adaptor Port
- ▶ Alternative to Indwelling Catheter for Bedridden Females
- ▶ Adapt Paste, Tape Strips, and Skin Gel Wipe Included
- ▶ Latex-Free

**HCPCS: A4328**

Stock No	Description	Length	Width
9840	Female Urinary Pouch	7-1/2"	4-1/2"

\* Meant to be connected to leg bag or gravity drainage

**Anti-Reflux:** Refers to a tube or collapsible material within a urine collection device to help prevent urine from re-entering the tubing

---

**Applicator Collar:** Found on some Hollister male external catheters; a plastic guide with notches for the thumb and forefinger to assist proper placement against the tip of the penis

---

**Bladder:** A balloon-like muscular organ which lies in the pelvis and functions to store and expel urine

---

**Bladder Catheterization:** A procedure in which a catheter is passed through the urethra or stoma into the bladder for the purpose of draining urine and performing diagnostic tests of bladder or urethral function

---

**Benzalkonium Chloride (BZK):** An antimicrobial solution used for disinfecting the urethral opening prior to inserting an intermittent catheter. Does not stain skin or clothing. Less likely to cause sensitivity

---

**Catheter (Urinary):** A rubber or vinyl hollow tube inserted through the urethra or a stoma to the bladder to withdraw urine or instill medication

---

**Closed System:** Refers to a "touchless" catheter located within a urine collection bag. The enclosed catheter is manipulated through the bag, guided through the introducer tip as it is inserted into the urethra or stoma without being directly touched

---

**Closed System Kit:** Closed system catheter includes the catheter located within the urine collection bag, and usually includes gloves, PVP or BZK swab sticks, and underpad drape

---

**Collection Bag:** A leg bag or bedside bag designed to collect urine by gravity that is connected by tubing to an indwelling or external catheter, urinary pouch, or also to wound drainage systems

---

**Color Coded:** On Hollister intermittent catheters, the funnel ends are color and size coordinated so catheters may be consistently ordered by size or funnel color

---

**Connector:** A device that secures the extension tubing to the catheter or urinary pouch system

---

**Coudé Tip:** A slight bend manufactured in the tip of the catheter which makes insertion past the prostate easier for some men. A notch at the funnel end is a guide for alignment during insertion

---

**External Condom Catheter:** Device that is secured externally with various adhesives to the penis shaft for the purpose of urine collection. Must be connected with tubing to a urine collection bag

---

**Extension Tubing:** Tubing that provides the connection between an external or indwelling catheter and a urine collection bag or leg bag. May also serve as a connection between an ostomy pouch or wound drainage collector. Made of latex or vinyl

---

**Eyelet:** An oval-shaped hole in the insertion end of the catheter to facilitate drainage of urine from the bladder. Eyelets are usually two in number, and may be across from each other or offset. Smooth eyelets protect the urethra from insertion trauma

---

**Female Urinary Pouch:** A cut-to-fit, one piece, vinyl pouch with a flexible synthetic barrier which is attached externally to the female urethral/vaginal area. The pouch is usually connected to a bedside collection bag. Recommended for the more immobile female

---

**Firm Catheter:** Refers to lesser pliability of a urethral catheter

---

**Flexible:** Refers to greater pliability of a urethral catheter or of a skin barrier found on the external female and male urinary pouches

---

**Flextend Barrier:** A skin barrier with special additives that achieve a stronger adhesive seal and are more resistant to breakdown from fluids. Skin-prepping agents are not recommended under Flextend Skin Barriers

---

**French Size:** Abbreviated Fr, the measuring gauge for the outer diameter of a straight or indwelling catheter

---

**Funnel:** Colored non-insertion end of an intermittent catheter allows for ease of fluid control. A catheter tipped syringe may be attached at the funnel for fluid instillation into the bladder. Color corresponds to specific French size

---

**Gel Reservoir:** A small, lubricant-filled flexible device unique to Hollister built into the Advance Plus closed catheter system. The catheter is evenly lubricated as it passes through the blue gel reservoir on its way to the urethra

---

**Grip Section:** A finely ribbed area on the catheter in a closed catheter system which assists with catheter control for continuous catheter advancement

---

**Hydrophilic Catheter:** A catheter designed to be lubricated when moistened with water which eases friction on the urethra

---

**Incontinence:** Loss of control of bowel and/or bladder

---

**Indwelling Catheter:** A flexible tube that remains in the bladder continuously to drain urine. May be referred to as a Foley catheter

---

**Inner Flap:** A thin, latex extra membrane inside Hollister Extended Wear latex male external catheters. The inner flap rests against the tip of the penis to prevent urine backwash and undermining of the catheter adhesive

---

**Intermittent Catheter:** A flexible tube that is used to empty the bladder on a regular schedule; used for self-catheterization

---

**Introducer Plate:** A nickel-size disc on the introducer tip that acts as a stop mechanism just outside the urethral opening. The plate helps to stabilize the catheter against the meatus area during insertion

---

**Introducer Tip:** A specially designed system on some intermittent catheters which functions to prevent urethra bacterial contamination of the catheter. It is inserted into the urethra. The 15 mm length on the Hollister introducer tip helps decrease the risk of introducing bacteria into the urethra. The tapered tip with multiple slits is designed for comfortable insertion

---

**Latex:** A material made from natural rubber. Red rubber catheters and some male external catheters are made from latex

---

**Leg Bag Straps:** Hold the leg bag in place

---

**Luer End:** No funnel on non insertion end. A syringe may be connected to the Luer end for bladder irrigation or insertion of medication

---

**Lubricant:** A water soluble jelly applied to a catheter prior to insertion which allows for easier insertion and helps minimize urethral irritation

---

**Male External Urinary Pouch:** A cut-to-fit, one-piece vinyl pouch with a synthetic adhesive backing that may be used for the very short penis, or the flush or retracted penis. Usually attached to a bedside collection bag, and is recommended more for those with limited mobility

---

**Meatus:** The opening of the urethra in both men and women

---

**Neurogenic Bladder:** An atonic or unstable bladder associated with a neurological condition, such as diabetes, stroke, or spinal cord injury

---

**Non Latex:** A material not made from rubber; usually vinyl or silicone for catheters or plastic for collection bags. Minimizes risk of allergic reaction

---

**Olive Tip:** An oval-shaped catheter tip which makes it easier for a female to find the urethral opening when catheterizing

---

**Pre-Lubricated Catheter:** Refers to a catheter which is lubricated by passing through a gel reservoir, or a lubricated catheter in a closed system, or a catheter in packaging which contains a lubricant packet which is ruptured to lubricate the catheter prior to opening the packaging

---

**PVP:** Povidone iodine solution used as an antimicrobial disinfectant applied to the urethral opening prior to catheter insertion. May stain clothing

---

**Ring Cap:** Cap on the introducer tip protects the catheter tip. Ring assists with cap removal especially for those with limited dexterity. Oval cap design prevents roll away

---

**Sizing Guide:** Used for external catheters, these semicircular cutouts help determine the diameter of the penis in order to guide selection for correct fit

---

**Soft Catheter:** A gentle, flexible catheter with limited firmness

---

**SoftFlex Barrier:** A flexible synthetic material that will conform to rounded and some irregular surfaces. Gentle enough to be removed in 24 hours or less

---

**Straight Tip:** Refers to the tapered, rounded insertion end of a catheter

---

**Urine Crystallization:** Crystals of salts and minerals may form in alkaline urine. Alkaline urine may also allow bacteria to grow in the bladder which may result in a urinary tract infection (UTI). Adequate fluid intake will keep urine acidic which will help prevent crystal formation and possible UTI (Citric juices and fruits turn alkaline in the stomach)

---

**Urethra:** The natural tube through which urine exits the bladder

---

**Vacuum Relief Valve:** A feature on the leg bag tubing that helps prevent collapse of the tubing, as well as collapse of an external urinary pouch or catheter

---

**Vented Leg Bag:** A special vacuum relief mechanism featured on the oval kink-resistant tubing and the companion pouch which also features an air vent. These combined features help to minimize a vacuum in the leg bag as well as in the tubing. A vacuum is created as the urine cools outside the body. The vacuum will stall urine drainage through the tubing and/or prevent the leg bag from emptying completely

---

# Educational Resources

Educational  
Resources

<b>Educational Materials</b>	<b>Literature No</b>
Patient Education Anatomy Guide	903789
Apogee Intermittent Catheter Guide for Use	908984
Apogee Closed System Guide for Use	908983
Advance Plus Guide for Use (English)	907266
Advance Plus Guide for Use (Spanish)	907373
Advance Plus Patient Education CD	906574
Understanding Neurogenic Bladder (English)	911195
Understanding Neurogenic Bladder (Spanish)	921745
Living with Spinal Cord Injury (English)	907500
Living with Spinal Cord Injury (Spanish)	907640
Pediatric Bladder Booklet – Girl (English)	908923
Pediatric Bladder Booklet – Girl (Spanish)	909964
Pediatric Bowel Booklet – Girl (English)	908925
Pediatric Bowel Booklet – Girl (Spanish)	909965
Pediatric Bladder Booklet – Boy (English)	908924
Pediatric Bladder Booklet – Boy (Spanish)	909966
Pediatric Bowel Booklet – Boy (English)	908926
Pediatric Bowel Booklet – Boy (Spanish)	909967
Catheter Coloring Book – Girl	907706
Catheter Coloring Book – Boy	907705
InView Male External Catheter Care Tips	911229
Everyday Male External Catheter Care Tips	921798
Extended Wear Male External Catheter Care Tips	921799
Clinician Documentation	910598
<hr/>	
<b>Clinical and Laboratory Evidence</b>	<b>Literature No</b>
The Hygienic Benefits of a No-Touch Intermittent Catheter with a Protective Tip and/or sleeve	911162
<hr/>	
<b>Secure Start Contenance Care Program</b>	<b>Literature No</b>
Intermittent Catheter, Male External Catheter and Urine Collector Prescription Form	910703
Leg Bag Prescription Form	907523
Patient Information	910767
Clinician Information	910768
The Benefits of a Comprehensive Urological Discharge Program	910601
<hr/>	
<b>Intermittent Catheter Shadow Buddies</b>	<b>Item No</b>
Light Skin Tone Boy	907281
Light Skin Tone Girl	907283
Medium Skin Tone Boy	907324
Medium Skin Tone Girl	907325
Dark Skin Tone Boy	907282
Dark Skin Tone Girl	907086
<hr/>	
<b>Global Sizing Guide</b>	911174

### Professional Organizations

American Academy of Physical Medicine and Rehabilitation (AAPMR)  
www.aapmr.org  
1.847.737.6000

American Spinal Injury Association (ASIA)  
www.asia-spinalinjury.org  
1.404.355.9772

American Urological Association (AUA)  
www.auanet.org  
1.410.689.3700

Association of Rehabilitation Nurses (ARN)  
www.rehabnurse.org  
1.800.229.7530

International Continence Society  
www.icsoffice.org  
+44.117.9444881

Society of Urological Nurses and Associates (SUNA)  
www.sunna.org  
1.888.827.7862

Wound, Ostomy and Continence Nurses Society (WOCN)  
www.wocn.org  
1.888.224.WOCN (9626)

---

### Support Organizations

The Buoniconti Fund to Cure Paralysis  
www.thebuonicontifund.com  
1.888.STANDUP (782.6387)

The Miami Project to Cure Paralysis  
www.miamiproject.miami.edu  
1.800.STAND-UP (782.6387)

National Association for Continence (NAFC)  
www.nafc.org  
1.800.BLADDER (252.3337)

National Rehabilitation Information Center  
www.naric.com  
1.800.346.2742

National Spinal Cord Injury Association  
www.spinalcord.org  
1.800.962.9629

Paralyzed Veterans of America  
www.pva.org  
1.800.555.9140

The Simon Foundation for Continence  
www.simonfoundation.org  
1.800.23SIMON (237.4666)

Spina Bifida Association  
www.sbaa.org  
1.800.621.3141

**Support Organizations (continued)**

Spinal Cord Injury Information Network

[www.spinalcord.uab.edu](http://www.spinalcord.uab.edu)

1.205.934.3283

Us Too International, Inc.

[www.ustoo.org](http://www.ustoo.org)

1.800.80.US-TOO (1.800.808.7866)

# Index • Stock Number

Stock No	Page(s)	Stock No	Page(s)	Stock No	Page(s)
10xx	92	1915x	35	509349	131
1410x	27	1920x	34	509372	132
1420x	25	1925x	35	509373	132
143x	78	31003	109	509374	132
1430x	26	31005	109	509375	131
1440x	28	31008	109	509376	131
1450x	29	31009	109	509377	131
1460x	25	314x	79	509378	130
1470x	26	318x	79	509379	130
1480x	28	32004	108	509380	133
149xx	29	32005	108	509381	133
1520x	25	32006	108	509382	133
1540x	28	32007	108	509384	133
1550x	29	322x	77	509387	133
1560x	25	323x	77	509388	133
1580x	28	327x	78	509389	133
1590x	29	33004	108	509390	133
1610x	26	332x	78	509393	133
1640x	26	334x	75	509394	133
1796	51	340x	76	509395	133
1800x	30	360x	77	509396	133
1801x	32	366x	77	517210	69, 138
1810x	31	370x	73	517220	70, 138
1811x	30	372x	72	519904	136
1812x	31	373x	74	519906	136
1813x	30	374x	73	519921	136
1814x	31	376x	72	519923	136
1815x	33	3761	60	519925	136
1816x	31	3774	61	519930	136
1817x	31	3777	59	519932	136
1818x	30	3778	59	519935	136
1819x	30	3795	59	519953	136
1820x	31	3796	59	519954	136
1828x	30	3797	60	519955	136
18282	60	3799	60	519956	136
1829x	30	380x	74	519957	136
18292	60	381x	74	519958	136
1832x	32	382x	76	519959	136
1833x	33	383x	74	519963	136
1835x	33	384x	75	519965	136
1836x	32	387x	75	529300	134
1837x	32	389x	74	529301	134
1838x	32	390x	76	529302	134
18382	61	506487	130	529303	134
1839x	32	506488	130	529304	134
18392	61	506489	130	529305	134
1840x	33	509338	130	529306	134
18402	61	509339	130	529937	135
1841x	33	509341	130	529938	135
18412	61	509342	130	529939	135
1842x	33	509343	131	529940	135
1873x	33	509344	131	529967	135
1875x	33	509345	131	529968	135
1900x	34	509346	131	529969	135
1905x	35	509347	131	529972	135
1910x	34	509348	131	529973	135

<u>Stock No</u>	<u>Page(s)</u>	<u>Stock No</u>	<u>Page(s)</u>	<u>Stock No</u>	<u>Page(s)</u>
529974	135	834x	75	9343	97
529975	137	843x	52	9345	97
529976	137	844x	51	9346	97
529977	138	8440	51	9348	96
529978	138	8450	51	9349	96
529979	138	8460	51	940xx	94
716x	78	846x	51	941xx	94
7299	68	847x	52	951xx	94
7300	68	848x	52	960xx	94
7331	69	849x	52	961xx	94
7700	71	850x	46	9645	96
771x	69	851x	47	9655	96
7721	70	852x	46	970x	49
7723	70	853x	44	971xx	94
7724	70	8531	41	974xx	95
7728	70	855x	44	975xx	95
773x	69	8551	41	976xx	95
7730	68	858x	46	977x	49
7731	68	859x	47	977xx	95
7736	69	861x	47	9779	49
7740	70	862x	47	978xx	95
7760	68	8631	41	9780	116
7766	71	863x	44	9781	115
7767	70	8641	41	9782	115
780x	67	864x	45	9785	116
7800	71	8651	41	9786	116
78500	67	872x	72	9788	115
78501	67	873x	73	979xx	95
7905	71	874x	73	9799	115
7906	67	875x	72	9805	96
7910	71	877x	73	9811	98
7917	68	8770	69	9814	96
79300	67	8800	68	9821	110
79301	67	8801	68	9822	110
79400	67	881x	74	9823	110
795xx	67	883xx	45	9839	71, 98
796xx	67	88300	42	9840	98
80xx	45	88302	43	9873	98
8031	42	884xx	45	9880	110
8081	42	88400	42	Bxxx	93
81xx	45	88401	43	BxxxB	93
8131	42	88402	43	HB2214	129
8181	42	88501	44	HB4414	129
80110	48	887xx	46	HB6614	129
80070	48	88700	43	HBT0906	129
81070	48	88800	43	HBHD4450	129
81110	48	890x	75	HBHD6675	129
821xx	51	89001	52	HBRF4414	129
82100	50	89002	53	HBRF2258	129
823xx	50	89003	53	HBRF4475	129
82300	50	89004	53	HBRF2650	129
82302	50	89005	53		
82400	50	910x	95		
82402	50	920x	95		
82500	50	931xx	94		
8331	41	9342	97		

Name	Page(s)
<b>A</b> ctiFlo Accessories.....	109
ActiFlo Catheter Kits with Disposable Odor-Barrier Collection Bags.....	108
ActiFlo Catheter Kits with Drainable Odor-Barrier Collection Bag .....	108
Adapt Odor Eliminator and Lubricant.....	67
Adapt Ostomy Belts.....	68
Adapt Skin Barrier Accessories .....	67
Adapt Skin Barrier Paste.....	67
Adapt Stoma Powder.....	67
Adhesives and Removers.....	68
Advance Plus Touch Free Intermittent Catheter .....	94
Advance Plus Touch Free Intermittent Catheter System .....	94
Anchor Fast Oral Endotracheal Tube Fastener .....	115
Apogee Plus System.....	93
Apogee Plus System Kit.....	93
Apogee Intermittent Catheter .....	92
<b>B</b> edside Drainage Collection System with Anti-Reflux Valve.....	71, 98
<b>C</b> enterPointLock Closed Pouches .....	75
CenterPointLock Drainable Pouches.....	74–75
CenterPointLock Flat Skin Barriers with Floating Flanges.....	72–73
CenterPointLock Stoma Cap .....	76
CenterPointLock Urostomy Pouches.....	75–76
Convex Skin Barriers with Floating Flanges.....	73–74
<b>D</b> rainable Pouch Clamp .....	69
Drain Tube Attachment Device.....	116
<b>E</b> veryday Male External Catheter .....	95
Extended Wear Male External Catheter.....	95
<b>F</b> abric Leg Bag Straps .....	97
Fecal Collectors.....	110
Feeding Tube Attachment Device.....	116
Female Urinary Pouch External Collection Device.....	98
<b>H</b> ollister Universal Catheter Access Port .....	49
Horizontal Tube Attachment Device.....	115
Hydrofera Blue Foam Dressings .....	129
<b>I</b> nstaFlo Catheter Kit.....	108
InView Male External Catheters .....	95
Irrigation Products .....	70
Irrigator Sleeve .....	76
<b>K</b> araya 5 Closed Pouches .....	78
Karaya 5 Drainable Pouches.....	77–78
Karaya 5 Urostomy Pouches .....	78
<b>L</b> atex-Free Extension Tubing.....	97
Loop Ostomy Bridge.....	70
<b>m</b> 9 Cleaner/Decrystallizer .....	69
m9 Odor Eliminators.....	68
Male Urinary Pouch External Collection Device.....	98

<b>Name</b>	<b>Page(s)</b>
<b>N</b> ew Image Closed Pouches.....	32–33
New Image Convex Skin Barriers with Integrated Floating Flanges .....	28–29
New Image Drainable Pouches .....	30–31
New Image Flat Skin Barriers with Integrated Floating Flanges.....	25–27
New Image High Output Drainable Pouches .....	32
New Image Irrigator Sleeve .....	33
New Image Non Sterile Kits — Drainable.....	34
New Image Non Sterile Kits — Urostomy.....	34
New Image Sterile Kits — Drainable.....	35
New Image Sterile Kits — Urostomy.....	35
New Image Urostomy Pouches.....	33
<b>O</b> ther Closed Pouches .....	79
Other Stoma Cap.....	79
<b>P</b> ouchkins Newborn Pouch.....	59
Pouchkins One-Piece Pouches.....	59–60
Pouchkins Pediatric Ostomy Belt.....	61
Pouchkins Premie Pouch .....	59
Pouchkins Two-Piece Pouches .....	60–61
Pouchkins Two-Piece Skin Barrier.....	60
Premier Closed Pouches.....	50–51
Premier Drainable Pouches with Convex Barriers.....	46–47
Premier Drainable Pouches with Flat Barriers .....	41–46
Premier High Output Drainable Pouches.....	48
Premier Non Sterile Kits — Drainable .....	52–53
Premier Non Sterile Kits — Urostomy .....	53
Premier Stoma Cap .....	51
Premier Urostomy Pouches with Convex Barriers.....	52
Premier Urostomy Pouches with Flat Barriers .....	51
Premier Wound Drainage Collectors.....	49
<b>R</b> eplacement Filters.....	71
Restore Adhesive Foam Dressing with TRIACT Technology .....	131
Restore Calcium Alginate Dressing Silver .....	135
Restore Calcium Alginate Dressing.....	135
Restore Cleanser & Moisturizer .....	138
Restore Contact Layer Dressing with Silver and TRIACT Technology.....	130
Restore Contact Layer Dressing with TRIACT Technology.....	130
Restore Contact Layer FLEX Dressing with TRIACT Technology .....	130
Restore DimethiCreme Skin Protectant.....	138
Restore DUO Absorbent Dressing with TRIACT Technology .....	131
Restore Extra Thin Hydrocolloid Dressing.....	136
Restore Foam Border Dressing Featuring TRIACT ADVANCED Technology.....	133
Restore Foam Dressing without Border .....	133
Featuring TRIACT ADVANCED Technology	
Restore Heel Dressing without Border .....	133
Restore Hydrocolloid Dressing.....	136
Restore Hydrocolloid Dressing with Foam Backing.....	136
Restore Hydrogel Dressing.....	135

Name	Page(s)
Restore Hydro-Shield Foam Dressing .....	134
Restore Hydro-Shield Foam Dressing without Border .....	134
Restore LITE Foam Dressing without Border .....	133
Featuring TRIACT ADVANCED Technology	
Restore Non-Adhesive Foam Dressing with Silver and .....	131
TRIACT Technology	
Restore Non-Adhesive Foam Dressing with TRIACT Technology .....	131
Restore Odor-Absorbent Dressing .....	136
Restore Skin Care Products .....	69
Restore Skin Cleanser .....	69, 138
Restore Skin Conditioning Creme .....	138
Restore TRIO Absorbent Dressing, Adhesive, with TRIACT Technology .....	132
Restore Wound Cleanser .....	137
<b>S</b> kin Barriers (Other) .....	68, 71
Skin Barrier Pastes and Powders (Other) .....	71
Skin Prepping Product .....	68
Suction Tube Attachment Device.....	116
<b>U</b> rinary Leg Bag.....	96
Urinary Leg Bag Combination Pack .....	96
Urostomy Drain Tube Adapter .....	69
<b>V</b> ented Urinary Leg Bag Combination Pack.....	96
Vertical Tube Attachment Device.....	115

# Index • HCPCS Codes

HCPCS	Page(s)	HCPCS	Page(s)
A4326	98	A5055	51, 76, 79
A4328	98	A5056	41
A4330	110	A5057	46
A4331	97	A5061	42, 43, 44, 52, 53, 59
A4349	95	A5062	59
A4351	92	A5063	30, 31, 34, 35, 60, 74, 75
A4352	92	A5071	60
A4353	93, 94	A5114	97
A4357	71, 98	A5120	68
A4358	96	A5122	68
A4362	59, 71	A5131	69
A4363	69	A5200	115, 116
A4364	68	A6154	48, 49
A4367	61, 68	A6196	135
A4368	71, 74, 77	A6197	135
A4371	67, 71	A6199	135
A4373	28, 29, 74	A6206	130
A4385	67, 68	A6207	130
A4387	78	A6208	130
A4388	41, 44, 45, 53	A6209	129, 131, 133, 134, 136
A4389	77, 78	A6210	129, 131, 133, 134
A4390	46, 47	A6211	134, 136
A4394	67, 69	A6212	129, 131, 133, 134
A4397	33, 70, 76	A6213	131, 133, 134
A4398	70	A6215	129
A4399	70	A6231	135
A4402	70	A6234	131, 136
A4405	71	A6235	131, 136
A4406	67	A6236	136
A4407	28, 29, 73	A6237	132
A4408	28	A6238	132
A4409	25, 26, 27, 34, 35, 72, 73	A6248	135
A4410	25, 72	A6250	70, 138
A4411	67	A6260	137
A4412	32, 75	A6266	135
A4414	25, 26, 34, 35, 60, 72, 73	A9270	69, 96, 138
A4415	25, 34, 35, 72, 73	A9900	115
A4416	50, 51	A9999	115, 116
A4417	78	B9998	116
A4418	79		
A4419	32, 61, 75		
A4421	69, 70		
A4424	42, 43, 45, 46, 59		
A4425	30, 31, 60, 74		
A4428	51, 53		
A4429	78		
A4430	52		
A4432	33, 34, 35, 61, 75, 76		
A4455	68		
A4456	68		
A5052	59		
A5054	33		

Hollister and logo, Hollister, ActiFlo, Adapt, Advance Plus, AF300, Anchor Fast, Apogee, Apogee Plus, CenterPointLock, ComfortWear, Extended Wear, Everyday, FlexTend, FlexWear, FormaFlex, HollisHesive, InstaFlo, InView, Karaya 5, Lock 'n Roll, m9, New Image, Pouchkins, Premier, Premium, QuietWear, and SoftFlex are trademarks of Hollister Incorporated.

C3Life, Secure Start, and "Hollister Ostomy. Details Matter." are service marks of Hollister Incorporated.

Hollisterwoundcare and wave logo are trademarks of Hollister Incorporated. Restore, Hydro-Shield and graphic, Restore DUO, Restore TRIO, TRIACT, TRIACT ADVANCED and graphic are trademarks of Hollister Wound Care LLC.

ConnectEd is a service mark of Hollister Wound Care LLC.

Hydrofera Blue is a trademark of Hydrofera, LLC.

Shadow Buddies is a trademark of Shadow Buddies.

Teflon is a trademark of E.I. du Pont de Nemours and Company.

Coloplast and SenSura are trademarks of Coloplast A/S.

ConvaTec and SUR-FIT Natura are trademarks of ConvaTec, LLC.

Zassi is a trademark of Zassi Medical Evolutions, Inc.

© 2012 Hollister Incorporated. Printed in USA.

## Hollister Consumer Programs

For consumer and healthcare professional product questions and sampling needs:

**1.888.740.8999**

For dealers, distributors, institutions, and pharmacy questions or orders:

**1.800.323.4060**

**fax 1.847.680.1017**

For helping new patients adjust to life with an ostomy or bladder management:

**Secure Start Services**

**1.888.808.7456**

For additional information on Hollister Wound Care:

**[www.hollisterwoundcare.com](http://www.hollisterwoundcare.com)**

For helping patients and families dealing with Epidermolysis Bullosa (EB) contact our EB Resources Hotline:

**1.888.EB AID4U**

**(1.888.322.4348)**



**Hollister Incorporated**  
2000 Hollister Drive  
Libertyville, Illinois 60048 USA  
**[www.hollister.com](http://www.hollister.com)**