

Ostomy Care

Table of Contents • Ostomy Care

Quick Reference/Product Selection Guides	pages 11-16
New Image Two-Piece Pouching System Skin Barriers	11
New Image Two-Piece Pouching System Pouches	12-13
New Image Two-Piece Pouching System Single Use Kits	13
Premier One-Piece Pouching System Pouches	14-16
Premier One-Piece Pouching System Single-Use Kits	16
Medicare Utilization Chart	17-19
New Image Two-Piece Pouching System	pages 20-29
Skin Barriers	20-23
Drainable Pouches	23-25
High Output Drainable Pouches	25
Closed Pouches	26
Urostomy Pouches	27
Irrigator Sleeve	27
Stoma Cap	27
Non-Sterile and Sterile Kits	27-29
Premier One-Piece Pouching System	pages 30-39
Drainable Pouches	30-34
Post-Operative Drainable Pouches	35
High Output Drainable Pouches	35
Closed Pouches	36
Urostomy Pouches	37-38
Non-Sterile Kits	38-39
Wound Drainage Collectors	39
Pouchkins Pediatric Pouching Systems	pages 40-43
Premie Pouch	40
Newborn Pouch	40
One-Piece Pouches	40-41
Two-Piece Skin Barrier	41
Two-Piece Pouches	41-42
Accessories	42
Other Pouchkins Products	42-43
Adapt Accessories	pages 44-45
Skin Barrier Accessories	44
Adapt Barrier Rings	44
Adapt Convex Barrier Rings	44
Adapt Barrier Strips	44
Skin Barrier Pastes and Powders	44
Adapt Paste	44
Adapt Stoma Powder	44
Skin Barriers	44
Flextend (extended wear) Skin Barrier	44
Skin Gel Wipes	44

Table of Contents • Ostomy Care





Adapt Accessories (continued)	pages 44-45
Skin Gel Protective Dressing Wipes	44
Adhesives and Removers	44-45
Universal Remover Wipes	44
Medical Adhesive	45
Medical Adhesive Remover	45
Odor Eliminators & Lubricants	45
Adapt Lubricating Deodorant	45
Odor Eliminator Spray	45
Cleaner/Decrystallizer	45
Ostomy Belt, Drainable Pouch Clamp, and Urostomy Drain Tube Adapter	45
OTHER OSTOMY PRODUCTS	
Two-Piece Pouching Systems	pages 46-50
CenterPointLock System	46-50
One-Piece Pouching Systems	pages 51-53
Karaya 5 System	51-52
Other One-Piece Products	53
Ostomy Accessories	pages 54-55
Skin Care Products	54
Irrigation Products	54
Loop Ostomy Bridge	54
Replacement Filters	55
Bedside Drainage Collector	55
Skin Barrier Pastes and Powders	55
Karaya Paste	55
Karaya Powder	55
Skin Barriers	55
Premium (standard wear) Skin Barrier	55
HolliHesive (standard wear) Skin Barrier	55
Odor Eliminator Drops	55
Discontinued Products	pages 56-58
Glossary	pages 59-60

How to Size Hollister Ostomy Pouches




A properly fitted pouch is essential to minimize stomal and skin irritation and to obtain optimum pouch wear time. A stoma measuring card allows correct measurement of the stoma.

Measure the stoma accurately. Select the smallest barrier opening which fits comfortably around the stoma without actually touching it. Too large an opening will leave skin exposed and could lead to skin irritation or erosion.





New Image Two-Piece Flat Skin Barriers with Floating Flange (see pages 20-21)

New Image Flat Barriers		 1-3/4" (44mm)	 2-1/4" (57mm)	 2-3/4" (70mm)	 4" (102mm)	
FlexTend	Cut-to-Fit	With Tape	14602 up to 1-1/4" (32mm)	14603 up to 1-3/4" (44mm)	14604 up to 2-1/4" (57mm)	14606 up to 3-1/2" (89mm)
		Without Tape	15602 up to 1-1/4" (32mm)	15603 up to 1-3/4" (44mm)	15604 up to 2-1/4" (57mm)	—
	Pre-Sized	With Tape	1470x	14707 1-3/8" (35mm) 14708 1-1/2" (38mm) 14709 1-3/4" (44mm)	—	—
		Without Tape	1610x	16107 1-3/8" (35mm) 16108 1-1/2" (38mm)	—	—
FlexWear	Cut-to-Fit	With Tape	14202 up to 1-1/4" (32mm)	14203 up to 1-3/4" (44mm)	14204 up to 2-1/4" (57mm)	14206 up to 3-1/2" (89mm)
		Without Tape	15202 up to 1-1/4" (32mm)	15203 up to 1-3/4" (44mm)	15204 up to 2-1/4" (57mm)	—
	Pre-Sized	With Tape	1430x	14307 1-3/8" (35mm) 14308 1-1/2" (38mm) 14309 1-3/4" (44mm)	—	—
		Without Tape	1640x	16407 1-3/8" (35mm) 16408 1-1/2" (38mm)	—	—




New Image Two-Piece Convex Skin Barriers with Floating Flange (see pages 22-23)

New Image Convex Barriers		 1-3/4" (44mm)	 2-1/4" (57mm)	 2-3/4" (70mm)	
FlexTend	Cut-to-Fit	With Tape	14802 up to 1" (25mm)	14803 up to 1-1/2" (38mm)	14804 up to 2" (51mm)
		Without Tape	15802 up to 1" (25mm)	15803 up to 1-1/2" (38mm)	15804 up to 2" (51mm)
	Pre-Sized	With Tape	14901 5/8" (16mm) 14902 3/4" (19mm) 14903 7/8" (22mm) 14904 1" (25mm)	14905 1-1/8" (29mm) 14906 1-1/4" (32mm) 14907 1-3/8" (35mm) 14908 1-1/2" (38mm)	14909 1-5/8" (41mm) 14910 1-3/4" (44mm) 14911 2" (51mm)
		Without Tape	15903 7/8" (22mm) 15904 1" (25mm)	15905 1-1/8" (29mm) 15906 1-1/4" (32mm)	—
FlexWear	Cut-to-Fit	With Tape	14402 up to 1" (25mm)	14403 up to 1-1/2" (38mm)	14404 up to 2" (51mm)
		Without Tape	15402 up to 1" (25mm)	15403 up to 1-1/2" (38mm)	15404 up to 2" (51mm)
	Pre-Sized	With Tape	14503 7/8" (22mm) 14504 1" (25mm)	14505 1-1/8" (29mm) 14506 1-1/4" (32mm)	—
		Without Tape	15503 7/8" (22mm) 15504 1" (25mm)	15505 1-1/8" (29mm) 15506 1-1/4" (32mm)	—

New Image Two-Piece Drainable Pouches (see pages 23-25)

New Image Drainable Pouches			 1-3/4" (44mm)	 2-1/4" (57mm)	 2-3/4" (70mm)	 4" (102mm)
Lock 'n Roll	Filter	Transparent	18192 18292	18193 18293	18194 18294	—
		Beige	18182 18282	18183 18283	18184 18284	—
	No Filter	Transparent	18132	18133	18134	—
		Ultra-Clear	18002	18003	18004	18006
		Beige	18112	18113	18114	—
Pouch Clamp	Filter	Transparent	18162	18163	18164	—
		Beige	18142	18143	18144	—
	No Filter	Transparent	18102	18103	18104	—
		Ultra-Clear	18172	18173	18174	18176
		Beige	18122 18202	18123 18203	18124 18204	—
High Output With Soft Tap	No Filter	Ultra-Clear	18012	18013	18014	18016





New Image Two-Piece Closed Pouches and Stoma Cap (see pages 26-27)

New Image Closed Pouches			 1-3/4" (44mm)	 2-1/4" (57mm)	 2-3/4" (70mm)
9" (23cm)	Filter	Transparent	18362	18363	18364
		Beige	18372	18373	18374
	No Filter	Beige	18332 (box of 60) 18732 (box of 30)	18333 (box of 60) 18733 (box of 30)	18334 (box of 60) 18734 (box of 30)
7" (18cm)	Filter	Transparent	18382	18383	18384
		Beige	18392	18393	18394
	No Filter	Beige	18352 (box of 60) 18752 (box of 30)	18353 (box of 60) 18753 (box of 30)	18354 (box of 60) 18754 (box of 30)

New Image Stoma Caps				
4-1/2" (11cm)	Beige	18502	18503	18504

New Image Two-Piece Urostomy Pouches and Single Use Kits (see pages 27-29)

New Image Urostomy Pouches				
9" (23cm)	Transparent	18402	18403	18404
	Ultra-Clear	18422	18423	18424
	Beige	18412	18413	18414

New Image Single Use Kits				 1-3/4" (44mm)	 2-1/4" (57mm)	 2-3/4" (70mm)	 4" (102mm)
Drainable (FlexWear)	Lock 'n Roll	Non-Sterile	Ultra-Clear	19002	19003	19004	19006
		Sterile	Ultra-Clear	19052	19053	19054	19056
	Pouch Clamp	Non-Sterile	Ultra-Clear	19102	19103	19104	19106
		Sterile	Ultra-Clear	19152	19153	19154	19156
Urostomy (FlexTend)	Non-Sterile	Ultra-Clear	19202	19203	19204	—	
	Sterile	Ultra-Clear	19252	19253	19254	—	

Premier One-Piece Drainable Pouches with Flat Barriers (see pages 30-33)

Premier Pouches Drainable with Flat Barriers				Cut-To-Fit		Pre-Sized								
				up to 2-1/2" (64mm)	3/4" (19mm)	1" (25mm)	1-3/16" (30mm)	1-1/4" (32mm)	1-3/8" (35mm)	1-1/2" (38mm)	1-9/16" (40mm)	1-3/4" (44mm)	2" (51mm)	2-1/2" (64mm)
Flexextend, No Filter, With Tape	Lock 'n Roll	12" (30cm)	Transparent	—	—	8532	—	8538	—	8533	—	8539	—	—
			Ultra-Clear	8531	—	—	—	—	—	—	—	—	—	—
			Beige	8551	—	8552	—	8558	—	8553	—	8559	—	—
	Pouch Clamp	9" (23cm)	Transparent	—	—	—	—	—	—	—	—	—	—	—
			Beige	8641	8647	8642	—	8648	—	8643	—	8649	8644	—
		12" (30cm)	Transparent	8631	8637	8632	—	8638	—	8633	—	8639	8634	8635
Ultra-Clear	8651		—	—	—	—	—	—	—	—	—	—		
SoftFlex, Filter Without Tape	Lock 'n Roll	7" (18cm)	Transparent	88800*	—	88825	88830	—	88835	—	88840	—	—	—
			Beige	88700*	—	88725	88730	—	88735	—	88740	—	—	—
		12" (30cm)	Transparent	88400*	—	88425	88430	—	88435	—	88440	—	—	—
			Ultra-Clear	88401**	—	—	—	—	—	—	—	—	—	—
	Pouch Clamp	12" (30cm)	Beige	88300*	—	88325	88330	—	88335	—	88340	—	—	—
			Ultra-Clear	88501**	—	—	—	—	—	—	—	—	—	—

*Cut-to-fit only to 2-1/8" (55mm) **No filter

Premier One-Piece Drainable Pouches with Convex Barriers (see pages 33-34)

Premier Pouches Drainable with Convex Barriers			Cut-To-Fit				Pre-Sized									
			up to 1" (25mm)	up to 1-1/2" (38mm)	up to 2" (51mm)	3/4" (19mm)	7/8" (22mm)	1" (25mm)	1-1/8" (29mm)	1-1/4" (32mm)	1-3/8" (35mm)	1-1/2" (38mm)	1-5/8" (41mm)	1-3/4" (44mm)	2" (51mm)	
Flexextend, No Filter, With Tape	Lock 'n Roll	12" (30cm)	Trans- parent	8524	8528	85211	8510	8511	8512	8513	8514	8515	—	—	—	—
			Beige	8504	8508	85011	8590	8591	8592	8593	8594	8595	—	—	—	—
	Pouch Clamp	Trans- parent	8624	8628	86211	8610	8611	8612	8613	8614	8615	8616	8617	8618	8619	




Premier High-Output Drainable Pouches (see page 35)

Premier Pouches High Output Drainable with Flat Barriers				Cut-To-Fit	
				up to 2-3/4" (70mm)	up to 4-1/3" (110mm)
Flextend, No Filter, Without Tape	With Soft Tap	12" (30cm)	Non-Sterile	80070	80110
			Sterile	81070	81110

Premier One-Piece Closed Pouches (see page 36)

Premier Closed Pouches			Cut-To-Fit	Pre-Sized		
			up to 2-1/8" (55mm)	1" (25mm)	1-3/16" (30mm)	1-3/8" (35mm)
SoftFlex, Filter, Without Tape	9" (23cm)	Transparent	82400	—	—	—
		Beige	82300	82325	82330	82335
	7" (18cm)	Transparent	82500	—	—	—
		Beige	82100	82125	82130	82135

Premier One-Piece Urostomy Pouches (see pages 37-38)

Premier Urostomy Pouches		Cut-To-Fit				Pre-Sized											
		up to 1" (25mm)	up to 1-1/2" (38mm)	up to 2" (51mm)	up to 2-1/2" (64mm)	1/2" (13mm)	5/8" (15mm)	3/4" (19mm)	7/8" (22mm)	1" (25mm)	1-1/8" (29mm)	1-1/4" (32mm)	1-3/8" (35mm)	1-1/2" (38mm)	1-3/4" (44mm)	2" (51mm)	
Flexextend	Flat	Trans-parent	—	—	—	8460	—	—	8467	—	8462	—	8468	—	8463	8469	8464
		Ultra-Clear	—	—	—	8450	—	—	—	—	—	—	—	—	—	—	—
		Beige	—	—	—	8440	—	—	8447	—	8442	—	8448	—	8443	8449	8444
	Convex	Trans-parent	8474	8478	84711	—	8480	8481	8482	8483	8484	8485	8486	8487	8488	8489	—
		Beige	8434	8438	84311	—	8490	8491	8492	8493	8494	8495	8496	8497	8498	8499	—



Premier Single Use Kits (see pages 38-39)

Premier Single Use Kits				up to 2-1/8" (55mm)	up to 2-1/2" (64mm)
Drainable (SoftFlex)	Lock 'n Roll	Non-Sterile	Ultra-Clear	89001	—
	Pouch Clamp	Non-Sterile	Ultra-Clear	89002	—
Drainable (Flexextend)	Lock 'n Roll	Non-Sterile	Ultra-Clear	89001	89004
	Pouch Clamp	Non-Sterile	Ultra-Clear	89002	89003
Urostomy (Flexextend)	Non-Sterile		Ultra-Clear	—	89005

Medicare's Usual Maximum Quantity for Ostomy Supplies

The quantity of ostomy supplies needed by a patient is determined to a great extent by the type of ostomy, its location, its construction, and the condition of the skin surface surrounding the stoma. There will be variation according to individual patient need and their needs may vary over time. The table below lists the maximum number of items/units of service that are usually medically necessary. The actual quantity needed for a particular patient may be more or less than the amount listed depending on the factors that affect the frequency of barrier and pouch change. The medical necessity for use of a greater quantity of supplies than the amounts listed must be clearly documented in the patient's medical record.

*Per LCD for Ostomy Supplies (L11502) TriCenturion, Effective Jan. 1, 2006
Incorporated code changes per Medicare DME PSC Bulletin Bul20061201HCPCS2007*

HCPCS	Long Description	Quantity
A4331	Extension drainage tubing, any type, any length, with connector/adaptor, for use with urinary leg bag or urostomy pouch, each	*
A4357	Bedside drainage bag, day or night, with or without anti-reflux device, with or without tube, each	2 ea /mo
A4361	Ostomy faceplate, each	3 ea/6 mo
A4362	Skin barrier, solid, 4x4 or equivalent, each	20 ea/mo
A4363	Ostomy clamp, any type, replacement only, each	*
A4364	Adhesive, liquid or equal, any type, per oz	4 oz/mo
A4365	Adhesive remover wipes, any type, per 50	*
A4366	Adhesive, liquid or equal, any type, per oz	*
A4367	Ostomy belt, each	1 ea/mo
A4368	Ostomy filter, any type, each	*
A4369	Ostomy skin barrier, liquid (spray, brush, etc), per oz.	2 oz/mo
A4371	Ostomy skin barrier, powder, per oz.	10 oz/6 mo
A4372	Ostomy skin barrier, solid 4 x 4 or equivalent, with built-in convexity, each	*
A4373	Ostomy skin barrier, with flange (solid, flexible or accordion), with built-in convexity, any size, each	*
A4375	Ostomy pouch, drainable, with faceplate attached, plastic, each	*
A4376	Ostomy pouch, drainable, with faceplate attached, rubber, each	*
A4377	Ostomy Pouch drainable for use on faceplate, plastic, each	10 ea/mo
A4378	Ostomy pouch, drainable, for use on faceplate, rubber, each	*
A4379	Ostomy pouch, urinary, with faceplate attached, plastic, each	*
A4380	Ostomy pouch, urinary, with faceplate attached, rubber, each	*
A4381	Ostomy Pouch, urinary, for use with faceplate, plastic, each	10 ea/mo
A4382	Ostomy pouch, urinary, for use on faceplate, heavy plastic, each	*
A4383	Ostomy pouch, urinary, for use on faceplate, rubber, each	*
A4384	Ostomy faceplate equivalent, silicone ring, each	*
A4385	Ostomy skin barrier, solid 4x4 or equivalent, extended wear, without built-in convexity, each	*
A4387	Ostomy pouch, closed, with barrier attached, with built-in convexity (1 piece), each	*
A4388	Ostomy pouch, drainable, with extended wear barrier attached, (1 piece), each	*
A4389	Ostomy pouch, drainable, with barrier attached, with built-in convexity (1 piece), each	*
A4390	Ostomy pouch, drainable, with extended wear barrier attached, with built-in convexity (1 piece), each	*
A4391	Ostomy pouch, urinary, with extended wear barrier attached (1 piece), each	*

**No established usual maximum, refer to Medicare Utilization Guidelines.*

Medicare Utilization

HCPCS	Long Description	Quantity
A4392	Ostomy pouch, urinary, with standard wear barrier attached, with built-in convexity (1 piece), each	*
A4393	Ostomy pouch, urinary, with extended wear barrier attached, with built-in convexity (1 piece), each	*
A4394	Ostomy deodorant, with or without lubricant, for use in ostomy pouch, liquid, per fluid ounce	*
A4395	Ostomy deodorant for use in ostomy pouch, solid, per tablet	*
A4396	Ostomy belt with peristomal hernia support	*
A4397	Irrigation supply, sleeve, each	4 ea/mo
A4398	Ostomy irrigation supply, bag, each	2 ea/6 mo
A4399	Ostomy irrigation supply, cone/catheter, including brush	2 ea/6 mo
A4402	Lubricant, per ounce	4 oz/mo
A4404	Ostomy ring, each	10 ea/mo
A4405	Ostomy skin barrier, non-pectin based, paste, per ounce	4 oz/mo
A4406	Ostomy skin barrier, pectin-based, paste, per ounce	4 oz/mo
A4407	Ostomy skin barrier, with flange (solid, flexible or accordion), extended wear, with built-in convexity, 4x4 inches or smaller, each	*
A4408	Ostomy skin barrier, with flange (solid, flexible or accordion), extended wear, with built-in convexity, larger than 4x4 inches, each	*
A4409	Ostomy skin barrier, with flange (solid, flexible or accordion), extended wear, without built-in convexity, 4x4 inches or smaller, each	*
A4410	Ostomy skin barrier, with flange (solid, flexible or accordion), extended wear, without built-in convexity, larger than 4x4 inches, each	*
A4411	Ostomy skin barrier, solid 4x4 or equivalent, extended wear, with built-in convexity, each	*
A4412	Ostomy pouch, drainable, high output, for use on a barrier with flange (2 piece system), without filter, each	*
A4413	Ostomy pouch, drainable, high output, for use on a barrier with flange (2 piece system), with filter, each	*
A4414	Ostomy skin barrier, with flange (solid, flexible or accordion), without built-in convexity, 4x4 inches or smaller, each	20 ea/mo
A4415	Ostomy skin barrier, with flange (solid, flexible or accordion), without built-in convexity, larger than 4x4 inches, each	20 ea/mo
A4416	Ostomy pouch, closed, with barrier attached, with filter (1 piece), each	60 ea/mo
A4417	Ostomy pouch, closed, with barrier attached, with built-in convexity, with filter (1 piece), each	60 ea/mo
A4418	Ostomy pouch, closed, without barrier attached, with filter (1 piece), each	60 ea/mo
A4419	Ostomy pouch, closed, for use on barrier with non-locking flange, with filter (2 piece), each	60 ea/mo
A4420	Ostomy pouch, closed, for use on barrier with locking flange (2 piece), each	60 ea/mo
A4422	Ostomy absorbent material (sheet/pad/crystal packet) for use in ostomy pouch to thicken liquid stomal output, each	*
A4423	Ostomy pouch, closed, for use on barrier with locking flange, with filter (2 piece), each	60 ea/mo
A4424	Ostomy pouch, drainable, with barrier attached, with filter (1 piece), each	20 ea/mo
A4425	Ostomy pouch, drainable, for use on barrier with non-locking flange, with filter (2 piece system), each	20 ea/mo

*No established usual maximum, refer to Medicare Utilization Guidelines.

HCPCS	Long Description	Quantity
A4426	Ostomy pouch, drainable, for use on barrier with locking flange (2 piece system), each	20 ea/mo
A4427	Ostomy pouch, drainable, for use on barrier with locking flange, with filter (2 piece system) each	20 ea/mo
A4428	Ostomy pouch, urinary, with extended wear barrier attached, with faucet-type tap with valve (1 piece), each	*
A4429	Ostomy pouch, urinary, with barrier attached, with built-in convexity, with faucet-type tap with valve (1 piece), each	20 ea/mo
A4430	Ostomy pouch, urinary, with extended wear barrier attached, with built-in convexity, with faucet-type tap with valve (1 piece), each	*
A4431	Ostomy pouch, urinary, with barrier attached, with faucet-type tap, with valve (1 piece), each	20 ea/mo
A4432	Ostomy pouch, urinary, for use on barrier with non-locking flange, with faucet-type tap with valve (2 piece), each	20 ea/mo
A4433	Ostomy pouch, urinary, for use on barrier with locking flange (2 piece), each	20 ea/mo
A4434	Ostomy pouch, urinary, for use on barrier with locking flange, with faucet-type tap with valve (2 piece), each	20 ea/mo
A4450	Tape, non-waterproof, per 18 square inches	40 ea (18sq in)/mo
A4452	Tape, waterproof, per 18 square inches	40 ea (18sq in)/mo
A4455	Adhesive remover or solvent (for tape, cement, or other adhesive), per ounce	16 oz/6mo
A5051	Ostomy pouch, closed, with barrier attached (1 piece), each	60 ea/mo
A5052	Ostomy pouch, closed, without barrier attached, (1 piece), each	60 ea/mo
A5053	Ostomy pouch, closed, for use on faceplate, each	60 ea/mo
A5054	Ostomy pouch, closed, for use on barrier with flange (2 piece), each	60 ea/mo
A5055	Stoma Cap	31 ea/mo
A5061	Ostomy pouch, drainable, with barrier attached (1 piece), each	20 ea/mo
A5062	Ostomy pouch, drainable, without barrier attached (1 piece), each	20 ea/mo
A5063	Ostomy pouch, drainable, for use on barrier with flange (2 piece system), each	20 ea/mo
A5071	Ostomy pouch, urinary, with barrier attached (1 piece), each	20 ea/mo
A5072	Ostomy pouch, urinary, without barrier attached (1 piece), each	20 ea/mo
A5073	Ostomy pouch, urinary, for use on barrier with flange (2 piece), each	20 ea/mo
A5081	Continent device, plug for continent stoma	31 ea/mo
A5082	Continent device, catheter for continent stoma	1 ea/mo
A5093	Ostomy accessory, convex insert	10 ea/mo
A5102	Bedside drainage bottle with or without tubing, rigid or expandable, each	2 ea/6 mo
A5120	Skin barrier, wipes or swabs, each	150/6 mo
A5121	Skin barrier, solid, 6x6 or equivalent, each	20 ea/mo
A5122	Skin barrier, solid, 8x8 or equivalent, each	20 ea/mo
A5126	Adhesive or non-adhesive, disk or foam pad	20 ea/mo
A5131	Appliance cleaner, incontinence and ostomy appliances, per 16 oz	1 ea(16-oz btl)/mo
A6216	Gauze, non-impregnated, non-sterile, pad size 16 sq. in. or less, without adhesive border, each dressing	60 ea/mo

*No established usual maximum, refer to Medicare Utilization Guidelines.

New Image Two-Piece Pouching System

Color Matching System

To help you select compatible New Image components, every box label is color-coded. First, select the proper skin barrier; then select a matching color-coded pouch ("red" skin barrier box label is compatible to "red" pouch box label).



Flat Skin Barriers with Floating Flanges



1460x (CUT-TO-FIT)

Box of 5

- ▶ Flexextend (extended wear) Skin Barrier - Flat
- ▶ With Tape Border

HCPCS: A4409

HCPCS for 14606: A4410

Stock No	Color Match	Flange Size	Barrier Opening
14602	●	1-3/4" (44mm)	up to 1-1/4" (32mm)
14603	●	2-1/4" (57mm)	up to 1-3/4" (44mm)
14604	●	2-3/4" (70mm)	up to 2-1/4" (57mm)
14606	●	4" (102mm)	up to 3-1/2" (89mm)



1560x (CUT-TO-FIT)

Box of 5

- ▶ Flexextend (extended wear) Skin Barrier - Flat
- ▶ Without Tape Border

HCPCS for 15602 & 15603: A4409

HCPCS for 15604: A4410

Stock No	Color Match	Flange Size	Barrier Opening
15602	●	1-3/4" (44mm)	up to 1-1/4" (32mm)
15603	●	2-1/4" (57mm)	up to 1-3/4" (44mm)
15604	●	2-3/4" (70mm)	up to 2-1/4" (57mm)



1420x (CUT-TO-FIT)

Box of 5

- ▶ FlexWear (standard wear) Skin Barrier - Flat
- ▶ With Tape Border

HCPCS: A4414

HCPCS for 14206: A4415

Stock No	Color Match	Flange Size	Barrier Opening
14202	●	1-3/4" (44mm)	up to 1-1/4" (32mm)
14203	●	2-1/4" (57mm)	up to 1-3/4" (44mm)
14204	●	2-3/4" (70mm)	up to 2-1/4" (57mm)
14206	●	4" (102mm)	up to 3-1/2" (89mm)



1520x (CUT-TO-FIT)

Box of 5

- ▶ FlexWear (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border

HCPCS for 15202 & 15203: A4414

HCPCS for 15204: A4415

Stock No	Color Match	Flange Size	Barrier Opening
15202	●	1-3/4" (44mm)	up to 1-1/4" (32mm)
15203	●	2-1/4" (57mm)	up to 1-3/4" (44mm)
15204	●	2-3/4" (70mm)	up to 2-1/4" (57mm)

Flat Skin Barriers with Floating Flanges (continued)



1470x (PRE-SIZED)

Box of 5

- ▶ FlexTend (extended wear) Skin Barrier - Flat
- ▶ With Tape Border

HCPCS: A4409

Stock No	Color Match	Flange Size	Barrier Opening
14701	●	1-3/4" (44mm)	5/8" (16mm)
14702	●	1-3/4" (44mm)	3/4" (19mm)
14703	●	1-3/4" (44mm)	7/8" (22mm)
14704	●	1-3/4" (44mm)	1" (25mm)
14705	●	1-3/4" (44mm)	1-1/8" (29mm)
14706	●	1-3/4" (44mm)	1-1/4" (32mm)
14707	●	2-1/4" (57mm)	1-3/8" (35mm)
14708	●	2-1/4" (57mm)	1-1/2" (38mm)
14709	●	2-1/4" (57mm)	1-3/4" (44mm)



1610x (PRE-SIZED)

Box of 5

- ▶ FlexTend (extended wear) Skin Barrier - Flat
- ▶ Without Tape Border

HCPCS: A4409

Stock No	Color Match	Flange Size	Barrier Opening
16101	●	1-3/4" (44mm)	5/8" (16mm)
16102	●	1-3/4" (44mm)	3/4" (19mm)
16103	●	1-3/4" (44mm)	7/8" (22mm)
16104	●	1-3/4" (44mm)	1" (25mm)
16105	●	1-3/4" (44mm)	1-1/8" (29mm)
16106	●	1-3/4" (44mm)	1-1/4" (32mm)
16107	●	2-1/4" (57mm)	1-3/8" (35mm)
16108	●	2-1/4" (57mm)	1-1/2" (38mm)



1430x (PRE-SIZED)

Box of 5

- ▶ FlexWear (standard wear) Skin Barrier - Flat
- ▶ With Tape Border

HCPCS: A4414

Stock No	Color Match	Flange Size	Barrier Opening
14301	●	1-3/4" (44mm)	5/8" (16mm)
14302	●	1-3/4" (44mm)	3/4" (19mm)
14303	●	1-3/4" (44mm)	7/8" (22mm)
14304	●	1-3/4" (44mm)	1" (25mm)
14305	●	1-3/4" (44mm)	1-1/8" (29mm)
14306	●	1-3/4" (44mm)	1-1/4" (32mm)
14307	●	2-1/4" (57mm)	1-3/8" (35mm)
14308	●	2-1/4" (57mm)	1-1/2" (38mm)
14309	●	2-1/4" (57mm)	1-3/4" (44mm)



1640x (PRE-SIZED)

Box of 5

- ▶ FlexWear (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border

HCPCS: A4414

Stock No	Color Match	Flange Size	Barrier Opening
16401	●	1-3/4" (44mm)	5/8" (16mm)
16402	●	1-3/4" (44mm)	3/4" (19mm)
16403	●	1-3/4" (44mm)	7/8" (22mm)
16404	●	1-3/4" (44mm)	1" (25mm)
16405	●	1-3/4" (44mm)	1-1/8" (29mm)
16406	●	1-3/4" (44mm)	1-1/4" (32mm)
16407	●	2-1/4" (57mm)	1-3/8" (35mm)
16408	●	2-1/4" (57mm)	1-1/2" (38mm)

New Image Two-Piece Pouching System

Convex Skin Barriers with Floating Flanges



1480x (CUT-TO-FIT)

Box of 5

- ▶ Flexlend (extended wear) Skin Barrier - Convex
- ▶ With Tape Border

HCPCS: A4407

Stock No	Color Match	Flange Size	Barrier Opening
14802	●	1-3/4" (44mm)	up to 1" (25mm)
14803	●	2-1/4" (57mm)	up to 1-1/2" (38mm)
14804	●	2-3/4" (70mm)	up to 2" (51mm)



1580x (CUT-TO-FIT)

Box of 5

- ▶ Flexlend (extended wear) Skin Barrier - Convex
- ▶ Without Tape Border

HCPCS for 15802 & 15803: A4407

HCPCS for 15804: A4408

Stock No	Color Match	Flange Size	Barrier Opening
15802	●	1-3/4" (44mm)	up to 1" (25mm)
15803	●	2-1/4" (57mm)	up to 1-1/2" (38mm)
15804	●	2-3/4" (70mm)	up to 2" (51mm)



1440x (CUT-TO-FIT)

Box of 5

- ▶ FlexWear (standard wear) Skin Barrier - Convex
- ▶ With Tape Border

HCPCS: A4373

Stock No	Color Match	Flange Size	Barrier Opening
14402	●	1-3/4" (44mm)	up to 1" (25mm)
14403	●	2-1/4" (57mm)	up to 1-1/2" (38mm)
14404	●	2-3/4" (70mm)	up to 2" (51mm)



1540x (CUT-TO-FIT)

Box of 5

- ▶ FlexWear (standard wear) Skin Barrier - Convex
- ▶ Without Tape Border

HCPCS: A4373

Stock No	Color Match	Flange Size	Barrier Opening
15402	●	1-3/4" (44mm)	up to 1" (25mm)
15403	●	2-1/4" (57mm)	up to 1-1/2" (38mm)
15404	●	2-3/4" (70mm)	up to 2" (51mm)

New!



149xx (PRE-SIZED)

Box of 5

- ▶ Flexlend (extended wear) Skin Barrier - Convex
- ▶ With Tape Border

HCPCS: A4407

Stock No	Color Match	Flange Size	Barrier Opening
14901	●	1-3/4" (44mm)	5/8" (16mm)
14902	●	1-3/4" (44mm)	3/4" (19mm)
14903	●	1-3/4" (44mm)	7/8" (22mm)
14904	●	1-3/4" (44mm)	1" (25mm)
14905	●	2-1/4" (57mm)	1-1/8" (29mm)
14906	●	2-1/4" (57mm)	1-1/4" (32mm)
14907	●	2-1/4" (57mm)	1-3/8" (35mm)
14908	●	2-1/4" (57mm)	1-1/2" (38mm)
14909	●	2-3/4" (70mm)	1-5/8" (41mm)
14910	●	2-3/4" (70mm)	1-3/4" (44mm)
14911	●	2-3/4" (70mm)	2" (51mm)

CONVEX PRODUCT

Convex Skin Barriers with Floating Flanges (continued)



1590x (PRE-SIZED)

Box of 5

- ▶ Flextend (extended wear) Skin Barrier - Convex
- ▶ Without Tape Border

HCPCS: A4407

Stock No	Color Match	Flange Size	Barrier Opening
15903	●	1-3/4" (44mm)	7/8" (22mm)
15904	●	1-3/4" (44mm)	1" (25mm)
15905	●	2-1/4" (57mm)	1-1/8" (29mm)
15906	●	2-1/4" (57mm)	1-1/4" (32mm)



1450x (PRE-SIZED)

Box of 5

- ▶ FlexWear (standard wear) Skin Barrier - Convex
- ▶ With Tape Border

HCPCS: A4373

Stock No	Color Match	Flange Size	Barrier Opening
14503	●	1-3/4" (44mm)	7/8" (22mm)
14504	●	1-3/4" (44mm)	1" (25mm)
14505	●	2-1/4" (57mm)	1-1/8" (29mm)
14506	●	2-1/4" (57mm)	1-1/4" (32mm)



1550x (PRE-SIZED)

Box of 5

- ▶ FlexWear (standard wear) Skin Barrier - Convex
- ▶ Without Tape Border

HCPCS: A4373

Stock No	Color Match	Flange Size	Barrier Opening
15503	●	1-3/4" (44mm)	7/8" (22mm)
15504	●	1-3/4" (44mm)	1" (25mm)
15505	●	2-1/4" (57mm)	1-1/8" (29mm)
15506	●	2-1/4" (57mm)	1-1/4" (32mm)

Drainable Pouches



1818x/1819x

Box of 10

- ▶ Drainable Pouch
- ▶ With Integrated Filter
- ▶ Lock 'n Roll Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels
 - Transparent - Body Side Only
 - Beige - Both Sides

HCPCS: A4425

Transparent Stock No	Beige Stock No	Color Match	Flange Size	Pouch Length
18192	18182	●	1-3/4" (44mm)	12" (30cm)
18193	18183	●	2-1/4" (57mm)	12" (30cm)
18194	18184	●	2-3/4" (70mm)	12" (30cm)

For Your Information

- ▶ Each box of New Image Pouches includes one disposal bag per every pouch for discreet disposal.

New Image Two-Piece Pouching System

Drainable Pouches (continued)



1828x/1829x

Box of 20

- ▶ Drainable Mini-Pouch
- ▶ With Integrated Filter
- ▶ Lock 'n Roll Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels
 - Transparent - Body Side Only
 - Beige - Both Sides

HCPCS: A4425

Transparent Stock No	Beige Stock No	Color Match	Flange Size	Pouch Length
18292	18282	●	1-3/4" (44mm)	7" (18cm)
18293	18283	●	2-1/4" (57mm)	7" (18cm)
18294	18284	●	2-3/4" (70mm)	7" (18cm)



1811x/1813x/1800x

Box of 10

- ▶ Drainable Pouch
- ▶ Without Integrated Filter
- ▶ Lock 'n Roll Closure
- ▶ Odor-Barrier Rustle-Free Film -Transparent and Beige
- ▶ Odor-Barrier Film -Ultra-Clear
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels
 - Transparent and Ultra-Clear - Body Side Only
 - Beige - Both Sides

HCPCS: A5063

Transparent Stock No	Beige Stock No	Ultra-Clear Stock No	Color Match	Flange Size	Pouch Length
18132	18112	18002*	●	1-3/4" (44mm)	12" (30cm)
18133	18113	18003	●	2-1/4" (57mm)	12" (30cm)
18134	18114	18004	●	2-3/4" (70mm)	12" (30cm)
—	—	18006	●	4" (102mm)	12" (30cm)



1814x/1816x

Box of 10

- ▶ Drainable Pouch
- ▶ With Integrated Filter
- ▶ Pouch Clamp Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels
 - Transparent - Body Side Only
 - Beige - Both Sides

HCPCS: A4425

Transparent Stock No	Beige Stock No	Color Match	Flange Size	Pouch Length
18162	18142	●	1-3/4" (44mm)	12" (30cm)
18163	18143	●	2-1/4" (57mm)	12" (30cm)
18164	18144	●	2-3/4" (70mm)	12" (30cm)

*coming soon

Drainable Pouches (continued)



1810x/1812x/1817x

Box of 10

- ▶ Drainable Pouch
- ▶ Without Integrated Filter
- ▶ Pouch Clamp Closure
- ▶ Odor-Barrier Rustle-Free Film -Transparent and Beige
- ▶ Odor-Barrier Film -Ultra-Clear
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels
 - Transparent and Ultra-Clear - Body Side Only
 - Beige - Both Sides

HCPCS: A5063

Transparent Stock No	Beige Stock No	Ultra-Clear Stock No	Color Match	Flange Size	Pouch Length
18102	18122	18172*	●	1-3/4" (44mm)	12" (30cm)
18103	18123	18173	●	2-1/4" (57mm)	12" (30cm)
18104	18124	18174	●	2-3/4" (70mm)	12" (30cm)
—	—	18176	●	4" (102mm)	12" (30cm)



1820x

Box of 10

- ▶ Drainable Mini-Pouch
- ▶ Without Integrated Filter
- ▶ Pouch Clamp Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels - Both Sides

HCPCS: A5063

Beige Stock No	Color Match	Flange Size	Pouch Length
18202	●	1-3/4" (44mm)	9" (23cm)
18203	●	2-1/4" (57mm)	9" (23cm)
18204	●	2-3/4" (70mm)	9" (23cm)

High Output Drainable Pouches



1801x

Box of 10

- ▶ Drainable Pouch
- ▶ Without Integrated Filter
- ▶ Soft Tap
- ▶ Odor-Barrier Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panel - Body-Side Only

HCPCS: A4412

Ultra-Clear Stock No	Color Match	Flange Size	Pouch Length
18012	●	1-3/4" (44mm)	12" (30cm)
18013	●	2-1/4" (57mm)	12" (30cm)
18014	●	2-3/4" (70mm)	12" (30cm)
18016	●	4" (102mm)	12" (30cm)

*coming soon

New Image Two-Piece Pouching System

Closed Pouches



1836x/1837x

Box of 60

- ▶ Closed Pouch
- ▶ With Integrated Filter
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels
 - Transparent - Body Side Only
 - Beige - Both Sides

HCPSCS: A4419

Transparent Stock No	Beige Stock No	Color Match	Flange Size	Pouch Length
18362	18372	●	1-3/4" (44mm)	9" (23cm)
18363	18373	●	2-1/4" (57mm)	9" (23cm)
18364	18374	●	2-3/4" (70mm)	9" (23cm)



1838x/1839x

Box of 60

- ▶ Closed Mini-Pouch
- ▶ With Integrated Filter
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels
 - Transparent - Body Side Only
 - Beige - Both Sides

HCPSCS: A4419

Transparent Stock No	Beige Stock No	Color Match	Flange Size	Pouch Length
18382	18392	●	1-3/4" (44mm)	7" (18cm)
18383	18393	●	2-1/4" (57mm)	7" (18cm)
18384	18394	●	2-3/4" (70mm)	7" (18cm)



1833x 1873x

New!

Box of 30
Box of 60

- ▶ Closed Pouch
- ▶ Without Integrated Filter
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels - Both Sides

HCPSCS: A5054

Stock No Box of 30	Stock No Box of 60	Color Match	Flange Size	Pouch Length
18332	18732	●	1-3/4" (44mm)	9" (23cm)
18333	18733	●	2-1/4" (57mm)	9" (23cm)
18334	18734	●	2-3/4" (70mm)	9" (23cm)

New!



1835x 1875x

Box of 30
Box of 60

- ▶ Closed Mini-Pouch
- ▶ Without Integrated Filter
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels - Both Sides

HCPSCS: A5054

Stock No Box of 30	Stock No Box of 60	Color Match	Flange Size	Pouch Length
18352	18752	●	1-3/4" (44mm)	7" (18cm)
18353	18753	●	2-1/4" (57mm)	7" (18cm)
18354	18754	●	2-3/4" (70mm)	7" (18cm)

Urostomy Pouches



1840x/1841x/1842x

Box of 10

- ▶ Urostomy Pouch
- ▶ Odor-Barrier Rustle-Free Film -Transparent
- ▶ Odor-Barrier Film -Ultra-Clear
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels
 - Transparent and Ultra-Clear - Body Side Only
 - Beige - Both Sides

HPCPS: A4432

Transparent Stock No	Beige Stock No	Ultra-Clear Stock No	Color Match	Flange Size	Pouch Length
18402	18412	18422	●	1-3/4" (44mm)	9" (23cm)
18403	18413	18423	●	2-1/4" (57mm)	9" (23cm)
18404	18414	18424*	●	2-3/4" (70mm)	9" (23cm)

Irrigator Sleeve



1815x

Box of 5

- ▶ Drainable Irrigator Sleeve
- ▶ Odor-Barrier Film
- ▶ With Belt Tabs

HPCPS: A4397

Transparent Stock No	Color Match	Flange Size	Pouch Length
18152	●	1-3/4" (44mm)	30" (76cm)
18153	●	2-1/4" (57mm)	30" (76cm)
18154	●	2-3/4" (70mm)	30" (76cm)

Stoma Cap



1850x

Box of 25

- ▶ With Integrated Filter
- ▶ Odor-Barrier, Rustle-Free Film
- ▶ With ComfortWear Panels

HPCPS: A5055

Beige Stock No	Color Match	Flange Size	Pouch Length
18502	●	1-3/4" (44mm)	4-1/2" (11mm)
18503	●	2-1/4" (57mm)	4-1/2" (11mm)
18504	●	2-3/4" (70mm)	4-1/2" (11mm)

Non-Sterile Kits - Drainable



1900x

Box of 5

- Colostomy/Ileostomy Kit includes:
- ▶ (1) 1420x (FlexWear) with Tape Border
 - ▶ (1) 1800x (Lock 'n Roll Closure), Ultra-Clear
 - ▶ (1) 7767 Loop Ostomy Bridge (in 19006 only)
 - ▶ Single Use Package

HPCPS for 19003 & 19004: A5063, A4414

HPCPS for 19006: A5063, A4415

Ultra-Clear Stock No	Color Match	Flange Size	Barrier Opening
19002	●	1-3/4" (44mm)	up to 1-14" (32mm)
19003	●	2-1/4" (57mm)	up to 1-3/4" (44mm)
19004	●	2-3/4" (70mm)	up to 2-1/4" (57mm)
19006	●	4" (102mm)	up to 3-1/2" (89mm)

*coming soon

New Image Two-Piece Pouching System

Non-Sterile Kits - Drainable (continued)



1910x

Box of 5

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1420x (FlexWear) with Tape Border
- ▶ (1) 1817x (Pouch Clamp Closure), Ultra-Clear
- ▶ (1) 7767 Loop Ostomy Bridge (in 19106 only)
- ▶ Single Use Package

HCPCS for 19103 & 19104: A5063, A4414

HCPCS for 19106: A5063, A4415

Ultra-Clear Stock No	Color Match	Flange Size	Barrier Opening
19102*	●	1-3/4" (44mm)	up to 1-1/4" (32mm)
19103	●	2-1/4" (57mm)	up to 1-3/4" (44mm)
19104	●	2-3/4" (70mm)	up to 2-1/4" (57mm)
19106	●	4" (102mm)	up to 3-1/2" (89mm)

Non-Sterile Kits - Urostomy



1920x

Box of 5

Urostomy Kit includes:

- ▶ (1) 1460x (Flexextend) with Tape Border
- ▶ (1) 1842x (Ultra-Clear Pouch Film)
- ▶ (1) 7331
- ▶ Single Use Package

HCPCS: A4432, A4409

Ultra-Clear Stock No	Color Match	Flange Size	Barrier Opening
19202	●	1-3/4" (44mm)	up to 1-1/4" (32mm)
19203	●	2-1/4" (57mm)	up to 1-3/4" (44mm)
19204*	●	2-3/4" (70mm)	up to 2-1/4" (57mm)

Sterile Kits - Drainable



1905x

Box of 5

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1420x (FlexWear) with Tape Border
- ▶ (1) 1800x (Lock 'n Roll Closure), Ultra-Clear
- ▶ (1) 7767 Loop Ostomy Bridge (in 19056 only)
- ▶ Sterile, Single Use Package

HCPCS for 19053 & 19054: A5063, A4414

HCPCS for 19056: A5063, A4415

Ultra-Clear Stock No	Color Match	Flange Size	Barrier Opening
19052*	●	1-3/4" (44mm)	up to 1-1/4" (32mm)
19053	●	2-1/4" (57mm)	up to 1-3/4" (44mm)
19054	●	2-3/4" (70mm)	up to 2-1/4" (57mm)
19056	●	4" (102mm)	up to 3-1/2" (89mm)

*coming soon

Sterile Kits - Drainable (continued)



1915x

Box of 5

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1420x (FlexWear) with Tape Border
- ▶ (1) 1817x (Pouch Clamp Closure), Ultra-Clear
- ▶ (1) 7767 Loop Ostomy Bridge (in 19156 only)
- ▶ Sterile, Single Use Package

HCPCS for 19153 & 19154: A5063, A4414

HCPCS for 19156: A5063, A4415

Ultra-Clear Stock No	Color Match	Flange Size	Barrier Opening
19152*	●	1-3/4" (44mm)	up to 1-1/4" (32mm)
19153	●	2-1/4" (57mm)	up to 1-3/4" (44mm)
19154	●	2-3/4" (70mm)	up to 2-1/4" (57mm)
19156	●	4" (102mm)	up to 3-1/2" (89mm)

New Image Sterile Kits - Urostomy



1925x

Box of 5

Urostomy Kit includes:

- ▶ (1) 1460x (Flexextend) with Tape Border
- ▶ (1) 1842x (Ultra-Clear Pouch Film)
- ▶ (1) 7331
- ▶ Sterile, Single Use Package

HCPCS: A4432, A4410

Ultra-Clear Stock No	Color Match	Flange Size	Barrier Opening
19252	●	1-3/4" (44mm)	up to 1-1/4" (32mm)
19253	●	2-1/4" (57mm)	up to 1-3/4" (44mm)
19254*	●	2-3/4" (70mm)	up to 2-1/4" (57mm)

*coming soon

Premier One-Piece Pouching System

Premier Drainable Pouches with Flat Barriers

New!



88401 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border
- ▶ Without Integrated Filter
- ▶ Lock 'n Roll Closure
- ▶ Odor-Barrier Film
- ▶ With ComfortWear Panels - Body Side Only

HCPCS: A5061

Ultra-Clear Stock No	Barrier Opening	Pouch Length
88401	up to 2-1/2" (64mm)	12" (30cm)



8531/8551 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ Flexextend (extended wear) Skin Barrier - Flat
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Lock 'n Roll Closure
- ▶ Odor-Barrier Rustle-Free Film - Beige
- ▶ Odor-Barrier Film -Ultra-Clear
- ▶ With ComfortWear Panels
 - Ultra-Clear - Body Side Only
 - Beige - Both Sides

HCPCS: A4388

Ultra-Clear Stock No	Beige Stock No	Barrier Opening	Pouch Length
8531	8551	up to 2-1/2" (64mm)	12" (30cm)



8631/8651 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ Flexextend (extended wear) Skin Barrier - Flat
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Pouch Clamp Closure
- ▶ Odor-Barrier Rustle-Free Film -Transparent
- ▶ Odor-Barrier Film - Ultra-Clear
- ▶ With ComfortWear Panels - Body Side Only

HCPCS: A4388

Transparent Stock No	Ultra-Clear Stock No	Barrier Opening	Pouch Length
8631	8651	up to 2-1/2" (64mm)	12" (30cm)



8641 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Mini-Pouch
- ▶ Flexextend (extended wear) Skin Barrier - Flat
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Pouch Clamp Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels - Body Side Only

HCPCS: A4388

Transparent Stock No	Barrier Opening	Pouch Length
8641	up to 2-1/2" (64mm)	9" (23cm)

Premier Drainable Pouches with Flat Barriers (continued)



88300/88400 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border
- ▶ With Integrated Filter
- ▶ Lock 'n Roll Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels
 - Transparent - Body Side Only
 - Beige - Both Sides

HCPCS: A4424

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
88400	88300	up to 2-1/8" (55mm)	12" (30cm)



88800/88700 (CUT-TO-FIT)

Box of 20

- ▶ Drainable Mini-Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border
- ▶ With Integrated Filter
- ▶ Lock 'n Roll Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels
 - Transparent - Body Side Only
 - Beige - Both Sides

HCPCS: A4424

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
88800	88700	up to 2-1/8" (55mm)	7" (18cm)

New!



88501 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border
- ▶ Without Integrated Filter
- ▶ Pouch Clamp Closure
- ▶ Odor-Barrier Film
- ▶ With ComfortWear Panels - Body Side Only

HCPCS: A5061

Ultra-Clear Stock No	Barrier Opening	Pouch Length
88501	up to 2-1/2" (64mm)	12" (30cm)

For Your Information

- ▶ Cut-to-fit products have a 1/2" (13mm) starter hole and can be used as a pre-sized product for some small stomas.
- ▶ Cutting guides are provided on all cut-to-fit barriers.
- ▶ Drainable Pouches without Lock 'n Roll include one pouch clamp per box.
- ▶ Do not use Skin Gel Wipes (Stock No 7917) or other brands of skin preps with Flexextend skin barriers.

Premier Drainable Pouches with Flat Barriers (continued)



883xx/884xx (PRE-SIZED)

Box of 10

- ▶ Drainable Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border
- ▶ With Integrated Filter
- ▶ Lock 'n Roll Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels
 - Transparent - Body Side Only
 - Beige - Both Sides

HCPCS: A4424

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
88425	88325	1" (25mm)	12" (30cm)
88430	88330	1-3/16" (30mm)	12" (30cm)
88435	88335	1-3/8" (35mm)	12" (30cm)
88440	88340	1-9/16" (40mm)	12" (30cm)



888xx/887xx (PRE-SIZED)

Box of 20

- ▶ Drainable Mini-Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border
- ▶ With Integrated Filter
- ▶ Lock 'n Roll Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels
 - Transparent - Body Side Only
 - Beige - Both Sides

HCPCS: A4424

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
88825	88725	1" (25mm)	7" (18cm)
88830	88730	1-3/16" (30mm)	7" (18cm)
88835	88735	1-3/8" (35mm)	7" (18cm)
88840	88740	1-9/16" (40mm)	7" (18cm)



853x/855x (PRE-SIZED)

Box of 10

- ▶ Drainable Pouch
- ▶ FlexTend (extended wear) Skin Barrier - Flat
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Lock 'n Roll Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels
 - Transparent - Body Side Only
 - Beige - Both Sides

HCPCS: A4388

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
8532	8552	1" (25mm)	12" (30cm)
8538	8558	1-1/4" (32mm)	12" (30cm)
8533	8553	1-1/2" (38mm)	12" (30cm)
8539	8559	1-3/4" (44mm)	12" (30cm)

Drainable Pouches with Flat Barriers (continued)



863x (PRE-SIZED)

Box of 10

- ▶ Drainable Pouch
- ▶ Flexextend (extended wear) Skin Barrier - Flat
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Pouch Clamp Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels - Body Side Only

HPCPS: A4388

Transparent
Stock No

Barrier Opening

Pouch Length

8637	3/4" (19mm)	12" (30cm)
8632	1" (25mm)	12" (30cm)
8638	1-1/4" (32mm)	12" (30cm)
8633	1-1/2" (38mm)	12" (30cm)
8639	1-3/4" (44mm)	12" (30cm)
8634	2" (51mm)	12" (30cm)
8635	2-1/2" (64mm)	12" (30cm)



864x (PRE-SIZED)

Box of 10

- ▶ Drainable Mini-Pouch
- ▶ Flexextend (extended wear) Skin Barrier - Flat
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Pouch Clamp Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels - Both Sides

HPCPS: A4388

Beige
Stock No

Barrier Opening

Pouch Length

8647	3/4" (19mm)	9" (23cm)
8642	1" (25mm)	9" (23cm)
8648	1-1/4" (32mm)	9" (23cm)
8643	1-1/2" (38mm)	9" (23cm)
8649	1-3/4" (44mm)	9" (23cm)
8644	2" (51mm)	9" (23cm)

Drainable Pouches with Convex Barriers



852x/850x (CUT-TO-FIT)

Box of 5

- ▶ Drainable Pouch
- ▶ Flexextend (extended wear) Skin Barrier - Convex
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Lock 'n Roll Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels
 - Transparent - Body Side Only
 - Beige - Both Sides

HPCPS: A4390

Transparent
Stock No

Beige
Stock No

Barrier Opening

Pouch Length

8524	8504	up to 1" (25mm)	12" (30cm)
8528	8508	up to 1-1/2" (38mm)	12" (30cm)
85211	85011	up to 2" (51mm)	12" (30cm)



Premier One-Piece Pouching System

Drainable Pouches with Convex Barriers (continued)



851x/859x (PRE-SIZED)

Box of 5

- ▶ Drainable Pouch
- ▶ Flexextend (extended wear) Skin Barrier - Convex
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Lock 'n Roll Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels
 - Transparent - Body Side Only
 - Beige - Both Sides

HCPCS: A4390

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
8510	8590	3/4" (19mm)	12" (30cm)
8511	8591	7/8" (22mm)	12" (30cm)
8512	8592	1" (25mm)	12" (30cm)
8513	8593	1-1/8" (29mm)	12" (30cm)
8514	8594	1-1/4" (32mm)	12" (30cm)
8515	8595	1-3/8" (35mm)	12" (30cm)



862x (CUT-TO-FIT)

Box of 5

- ▶ Drainable Pouch
- ▶ Flexextend (extended wear) Skin Barrier - Convex
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Pouch Clamp Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels - Body Side Only

HCPCS: A4390

Transparent Stock No	Barrier Opening	Pouch Length
8624	up to 1" (25mm)	12" (30cm)
8628	up to 1-1/2" (38mm)	12" (30cm)
86211	up to 2" (51mm)	12" (30cm)



861x (PRE-SIZED)

Box of 5

- ▶ Drainable Pouch
- ▶ Flexextend (extended wear) Skin Barrier - Convex
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Pouch Clamp Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels - Body Side Only

HCPCS: A4390

Transparent Stock No	Barrier Opening	Pouch Length
8610	3/4" (19mm)	12" (30cm)
8611	7/8" (22mm)	12" (30cm)
8612	1" (25mm)	12" (30cm)
8613	1-1/8" (29mm)	12" (30cm)
8614	1-1/4" (32mm)	12" (30cm)
8615	1-3/8" (35mm)	12" (30cm)
8616	1-1/2" (38mm)	12" (30cm)
8617	1-5/8" (41mm)	12" (30cm)
8618	1-3/4" (44mm)	12" (30cm)
8619	2" (51mm)	12" (30cm)



Post-Operative Drainable Pouches



9763/9764 (CUT-TO-FIT)

Box of 5

- ▶ Drainable Pouch
- ▶ FlexWear (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border
- ▶ Without Integrated Filter
- ▶ Integrated Wire Closure
- ▶ Odor-Barrier Film
- ▶ Without ComfortWear Panels

HCPCS: A5061

**Ultra-Clear
Stock No**

Barrier Opening

Pouch Length

9763 (Non-Sterile)

up to 4" (102mm)

15-1/2" (39cm)

9764 (Sterile)

up to 4" (102mm)

15-1/2" (39cm)

High Output Drainable Pouches



80070/80110 (CUT-TO-FIT)

Box of 10

- ▶ Non-Sterile, Drainable Pouch
- ▶ Flexextend (extended wear) Skin Barrier - Flat
- ▶ With Soft Tap
- ▶ With Access Window
- ▶ Odor-Barrier Film
- ▶ With ComfortWear Panel - Body Side Only

HCPCS: A6154

**Ultra-Clear
Stock No**

Barrier Opening

Pouch Length

80070

up to 2-3/4" (70mm)

12" (30cm)

80110

up to 4-1/3" (110mm)

12" (30cm)



81070/81010 (CUT-TO-FIT)

Box of 5

- ▶ Sterile, Drainable Pouch
- ▶ Flexextend (extended wear) Skin Barrier - Flat
- ▶ With Soft Tap
- ▶ With Access Window
- ▶ Odor-Barrier Film
- ▶ With ComfortWear Panel - Body Side Only

HCPCS: A6154

**Ultra-Clear
Stock No**

Barrier Opening

Pouch Length

81070

up to 2-3/4" (70mm)

12" (30cm)

81010

up to 4-1/3" (110mm)

12" (30cm)

Premier One-Piece Pouching System

Closed Pouches



82300/82400 (CUT-TO-FIT)

Box of 30

- ▶ Closed Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border
- ▶ With Integrated Filter
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels
 - Transparent - Body Side Only & Beige - Both Sides

HCPCS: A4416

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
82400	82300	up to 2-1/8" (55mm)	9" (23cm)



823xx (PRE-SIZED)

Box of 30

- ▶ Closed Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border
- ▶ With Integrated Filter
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels - Both Sides

HCPCS: A4416

Beige Stock No	Barrier Opening	Pouch Length
82325	1" (25mm)	9" (23cm)
82330	1-3/16" (30mm)	9" (23cm)
82335	1-3/8" (35mm)	9" (23cm)



82100/82500 (CUT-TO-FIT)

Box of 30

- ▶ Closed Mini-Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border
- ▶ With Integrated Filter
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels
 - Transparent - Body Side Only
 - Beige - Both Sides

HCPCS: A4416

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
82500	82100	up to 2-1/8" (55mm)	7" (18cm)



821xx (PRE-SIZED)

Box of 30

- ▶ Closed Mini-Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border
- ▶ With Integrated Filter
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels - Both Sides

HCPCS: A4416

Beige Stock No	Barrier Opening	Pouch Length
82125	1" (25mm)	7" (18cm)
82130	1-3/16" (30mm)	7" (18cm)
82135	1-3/8" (35mm)	7" (18cm)

Stoma Cap



1796 (PRE-SIZED)

- ▶ Stoma Cap
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border
- ▶ With Integrated Filter
- ▶ Odor-Barrier Film

HCPCS: A5055

Beige Stock No	Barrier Opening	Size
1796	1-15/16" (50mm)	4" (102mm)

Urostomy Pouches with Flat Barriers



8450/8460/8440 (CUT-TO-FIT)

Box of 10

- ▶ Urostomy Pouch
- ▶ Flexextend (extended wear) Skin Barrier - Flat
- ▶ With Tape Border
- ▶ Odor-Barrier Rustle-Free Film - Transparent
- ▶ Odor-Barrier Film - Ultra-Clear
- ▶ With ComfortWear Panels - Body Side Only

**Beige
Coming
Soon!**

HCPCS: A4428

Transparent Stock No	Ultra-Clear Stock No	Beige Stock No	Barrier Opening	Pouch Length
8460	8450	8440*	up to 2-1/2" (64mm)	9" (23cm)



846x/844x (PRE-SIZED)

Box of 10

- ▶ Urostomy Pouch
- ▶ Flexextend (extended wear) Skin Barrier - Flat
- ▶ With Tape Border
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels - Body Side Only

**Beige
Coming
Soon!**

HCPCS: A4428

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
8467	8447*	3/4" (19mm)	9" (23cm)
8462	8442*	1" (25mm)	9" (23cm)
8468	8448*	1-1/4" (32mm)	9" (23cm)
8463	8443*	1-1/2" (38mm)	9" (23cm)
8469	8449*	1-3/4" (44mm)	9" (23cm)
8464	8444*	2" (51mm)	9" (23cm)

Urostomy Pouches with Convex Barriers



847x/843x (CUT-TO-FIT)

Box of 5

- ▶ Urostomy Pouch
- ▶ Flexextend (extended wear) Skin Barrier - Convex
- ▶ With Tape Border
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels - Body Side Only

**Beige
Coming
Soon!**

HCPCS: A4430

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
8474	8434*	up to 1" (25mm)	9" (23cm)
8478	8438*	up to 1-1/2" (38mm)	9" (23cm)
84711	8434*	up to 2" (51mm)	9" (23cm)

*coming soon

Premier One-Piece Pouching System

Urostomy Pouches with Convex Barriers (continued)



**Beige
Coming
Soon!**

848x/849x (PRE-SIZED)

Box of 5

- ▶ Urostomy Pouch
- ▶ FlexTend (extended wear) Skin Barrier - Convex
- ▶ With Tape Border
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels - Body Side Only

HCPCS: A4430

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
8480	8490*	1/2" (13mm)	9" (23cm)
8481	8491*	5/8" (16mm)	9" (23cm)
8482	8492*	3/4" (19mm)	9" (23cm)
8483	8493*	7/8" (22mm)	9" (23cm)
8484	8494*	1" (25mm)	9" (23cm)
8485	8495*	1-1/8" (29mm)	9" (23cm)
8486	8496*	1-1/4" (32mm)	9" (23cm)
8487	8497*	1-3/8" (35mm)	9" (23cm)
8488	8498*	1-1/2" (38mm)	9" (23cm)
8489	8499*	1-3/4" (44mm)	9" (23cm)

*coming soon

Non-Sterile Kits - Drainable

New!



89001 (CUT-TO-FIT)

Box of 5 Kits

Colostomy/Ileostomy Kit Includes:

- ▶ (1) 88401 (SoftFlex)
- ▶ Without Tape Border
- ▶ Without Integrated Filter
- ▶ Lock 'n Roll Closure
 - ▶ Single Use Package

HCPCS: A5061

Ultra-Clear Stock No	Barrier Opening	Pouch Length
89001	up to 2-1/2" (64mm)	12" (30cm)

New!



89002 (CUT-TO-FIT)

Box of 5 Kits

Colostomy/Ileostomy Kit Includes:

- ▶ (1) 88501 (SoftFlex)
- ▶ Without Tape Border
- ▶ Without Integrated Filter
- ▶ Pouch Clamp Closure
- ▶ Single Use Package

HCPCS: A5061

Ultra-Clear Stock No	Barrier Opening	Pouch Length
89002	up to 2-1/2" (64mm)	12" (30cm)

New!



89003 (CUT-TO-FIT)

Box of 5 Kits

Colostomy/Ileostomy Kit Includes:

- ▶ (1) 8651 (FlexTend)
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Pouch Clamp Closure
- ▶ Single Use Package

HCPCS: A4388

Ultra-Clear Stock No	Barrier Opening	Pouch Length
89003	up to 2-1/2" (64mm)	12" (30cm)

 **CONVEX PRODUCT**

Non-Sterile Kits - Drainable (continued)

New!



- 89004 (CUT-TO-FIT)**
Colostomy/Ileostomy Kit Includes:
- ▶ (1) 8531 (Flexextend)
 - ▶ With Tape Border
 - ▶ Without Integrated Filter
 - ▶ Lock 'n Roll Closure
 - ▶ Single Use Package

Box of 5 Kits

HCPCS: A4388

Ultra-Clear Stock No	Barrier Opening	Pouch Length
89004	up to 2-1/2" (64mm)	12" (30cm)

New!



- 89005 (CUT-TO-FIT)**
Urostomy Kit Includes:
- ▶ (1) 8450 (Flexextend)
 - ▶ With Tape Border
 - ▶ (1) 7331
 - ▶ Single Use Package

Box of 5 Kits

HCPCS: A4428

Ultra-Clear Stock No	Barrier Opening	Pouch Length
89005	up to 2-1/2" (64mm)	9" (23cm)

Wound Drainage Collectors

Wound Drainage Collectors



- 970x/977x**
- ▶ For Management of Draining Wounds, Fistulae, or Draining Tube Sites
 - ▶ Available With and Without FlexWear Skin Barriers
 - ▶ Sterile and Non-Sterile Drainable Pouches
 - ▶ Odor-Barrier Film
 - ▶ With Access Window
 - ▶ Reclosable Cap
 - ▶ Without Starter Hole

With Barrier: Non-Sterile Stock No	HCPCS	Description	Quantity/Box
9773	A6154	Small (for wounds up to 3")	10
9776	A6154	Medium (for wounds up to 3-3/4")	10
9778	A6154	Large (for wounds up to 4" x 8")	3

With Barrier: Sterile Stock No	HCPCS	Description	Quantity/Box
9701	A6154	Small (for wounds up to 3")	5
9702	A6154	Medium (for wounds up to 3-3/4")	5
9703	A6154	Large (for wounds up to 4" x 8")	3

Without Barrier: Sterile Stock No	HCPCS	Description	Quantity/Box
9772	A6154	Small (for wounds up to 3")	10
9775	A6154	Medium (for wounds up to 3-3/4")	10
9777	A6154	Large (for wounds up to 4" x 8")	3

Hollister Universal Catheter Access Port



- 9779**
- ▶ Graduated nipple and punch allows introduction of catheter or drain through Wound Drainage Collector
 - ▶ Sterile

Stock No	Quantity/Box
9779	8

Pouchkins Pediatric Pouching System

Premie Pouch



3777 (CUT-TO-FIT)

Box of 15

- ▶ Use as a Premie Closed or Drainable Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border
- ▶ Without Starter Hole
- ▶ Without Integrated Filter
- ▶ Odor-Barrier Film
- ▶ Soft Wire Closure
- ▶ Without ComfortWear Panels

HCPCS: A5052, A4362

**Ultra-Clear
Stock No**

Barrier Opening

Pouch Length

3777

up to 5/8" (15mm)

4.5" (11.5cm)

Newborn Pouch



3778 (CUT-TO-FIT)

Box of 15

- ▶ Newborn Pouch; Use With or Without the Drain Valve
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border
- ▶ Without Integrated Filter
- ▶ Soft Wire Closure
- ▶ Odor-Barrier Film
- ▶ Without ComfortWear Panels

HCPCS: A5062, A4362

**Ultra-Clear
Stock No**

Barrier Opening

Pouch Length

3778

up to 1-3/8" x 7/8" (35mm x 22mm)

6.5" (16.5cm)

One-Piece Pouches

New!



3796 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Starter Hole
- ▶ Without Tape Border
- ▶ Without Integrated Filter
- ▶ Lock 'n Roll Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels - Body Side Only

HCPCS: A5061

**Transparent
Stock No**

Barrier Opening

Pouch Length

3796

up to 1-1/2" (38mm)

6-1/2" (16cm)

New!



3795 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Starter Hole
- ▶ Without Tape Border
- ▶ With Integrated Filter
- ▶ Lock 'n Roll Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels - Body Side Only

HCPCS: A5061

**Ultra-Clear
Stock No**

Barrier Opening

Pouch Length

3795

up to 2" (51mm)

7" (18cm)

One-Piece Pouches (continued)

New!



3797

Box of 10

- ▶ Urostomy Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Starter Hole
- ▶ Without Tape Border
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels - Body Side Only

HCPCS: A5071

**Transparent
Stock No**

Barrier Opening

Pouch Length

3797

up to 1-1/2" (38mm)

8-3/4" (22cm)

Two-Piece Skin Barrier

New!



3761

Box of 5

- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Starter Hole
- ▶ Without Tape Border

HCPCS: A4414

Stock No

**Transparent
Color Match**

Flange Size

Barrier Opening

3761

●

1-3/4" (44mm)

up to 1-1/4" (32mm)

Two-Piece Pouches

New!



3799

Box of 10

- ▶ Drainable Pouch
- ▶ Without Integrated Filter
- ▶ Lock 'n Roll Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels - Body Side Only

HCPCS: A5063

**Ultra-Clear
Stock No**

Color Match

Flange Size

Pouch Length

3799

●

1-3/4" (44mm)

6-1/2" (16cm)

New!



18292/18282

Box of 20

- ▶ Drainable Pouch
- ▶ With Integrated Filter
- ▶ Lock 'n Roll Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panel on Body Side

HCPCS: A4425

**Transparent
Stock No**

**Beige
Stock No**

**Color
Match**

**Flange
Size**

**Pouch
Length**

18292

18282

●

1-3/4" (44mm)

7" (18cm)

Pouchkins Pediatric Pouching System



18382/18392

Box of 60

- ▶ Closed Mini-Pouch
- ▶ With Integrated Filter
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panels
 - Transparent - Body Side Only
 - Beige - Both Sides

HCPCS: A4419

Transparent Stock No	Beige Stock No	Color Match	Flange Size	Pouch Length
18382	18392	●	1-3/4" (44mm)	7" (18cm)



18402/18412

Box of 10

- ▶ Urostomy Pouch
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ With ComfortWear Panel on Body Side

HCPCS: A4432

Transparent Stock No	Beige Stock No	Color Match	Flange Size	Pouch Length
18402	18412	●	1-3/4" (44mm)	9" (23cm)

Accessories



3774

Single Unit

- ▶ Pediatric Ostomy Belt

HCPCS: A4367

Stock No	Length
3774	10" - 17" (25-43cm)



3779

Single Unit

- ▶ Drain Tube
- ▶ For use with Newborn Ostomy Pouch #3778

HCPCS: A4331

Ultra-Clear Stock No	Length
3779	61" (155cm)

Other Pouchkins Products

One-Piece Pouches



9766*

Box of 15

- ▶ Drainable Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Starter Hole
- ▶ Without Tape Border
- ▶ Without Integrated Filter
- ▶ Soft Wire Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ Without ComfortWear Panels

HCPCS: A5061

Transparent Stock No	Barrier Opening	Pouch Length
9766	up to 2" (51mm)	8-1/2" (21cm)

*Plan to Discontinue 12-31-07

**9768*****Box of 15**

- ▶ Drainable Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Starter Hole
- ▶ Without Tape Border
- ▶ Without
- ▶ Soft Wire Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With ComfortWear Panels - Body Side Only

HCPCS: A5061**Ultra-Clear
Stock No****Barrier Opening****Pouch Length**

9768

up to 1-1/2" (38mm)

6-1/2" (16cm)

**9767*****Box of 15**

- ▶ Urostomy Pouch
- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border
- ▶ Odor-Barrier Rustle-Free Film
- ▶ Without ComfortWear Panels

HCPCS: A5071**Transparent
Stock No****Barrier Opening****Pouch Length**

9767

up to 1-1/2" (38mm)

8-1/2" (21cm)

Two-Piece Skin Barrier

**3770*** (CUT-TO-FIT)**Box of 5**

- ▶ SoftFlex (standard wear) Skin Barrier - Flat
- ▶ Without Tape Border

HCPCS: A4414**Transparent
Stock No****Letter Match****Flange Size****Pouch Length**

3770

①

1-3/4" (44mm)

up to 1-1/4" (32mm)

Two-Piece Pouches

**3771*****Box of 10**

- ▶ Drainable Pouch
- ▶ Without Integrated Filter
- ▶ Soft Wire Closure
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ Without ComfortWear Panels

HCPCS: A5063**Transparent
Stock No****Letter Match****Flange Size****Pouch Length**

3771

①

1-3/4" (44mm)

up to 8-1/2" (21cm)

**3772*****Box of 10**

- ▶ Urostomy Pouch
- ▶ Odor-Barrier Rustle-Free Film
- ▶ With Belt Tabs
- ▶ Without ComfortWear Panels

HCPCS: A5073**Transparent
Stock No****Letter Match****Flange Size****Pouch Length**

3772

①

1-3/4" (44mm)

up to 8-1/2" (21cm)

**Plan to Discontinue 12-31-07*

Adapt Accessories

Skin Barrier Accessories



780x
▶ Adapt Barrier Rings (Flat)
HCPCS: A4385

Box of 10

Stock No	Outer Diameter
7805	2" (48mm)
7806	4" (98mm)



795xx
▶ Adapt Convex Barrier Rings
HCPCS: A4411

Box of 10

Stock No	Inner Diameter
79520	13/16" (20mm) - can be stretched to 1" (25mm)
79530	1-3/16" (30mm) - can be stretched to 1-3/8" (35mm)
79540	1-9/16" (40mm) - can be stretched to 1-3/4" (45mm)



79400
▶ Adapt Barrier Strips
HCPCS: A4385

Box of 10

Stock No	Size	Length
79400	60 grams	5-1/8" (13cm)

Skin Barrier Paste



79300/79301
▶ Adapt Paste
HCPCS: A4406

Stock No	Description	Quantity
79300	2 oz. (57gm) tube	Box of 1
79301	0.5 oz (14gm) tube	Box of 20

Skin Barrier Powder



7906
▶ Adapt Stoma Powder
HCPCS: A4371

Box of 1

Stock No	Description
7906	1 oz. (28.3g) puff bottle

Skin Barrier



8800
8801
▶ FlexTend (extended wear) Skin Barrier
HCPCS for 8800: A4385
HCPCS for 8801: A5122

Box of 5
Box of 3

Stock No	Size
8800	4" x 4" (10cm x 10cm)
8801	8" x 8" (20cm x 20cm)

Skin Gel Products



7917
▶ Skin Gel Protective Dressing Wipes
HCPCS: A5120

12 Cartons of
50 Wipes Each

Stock No
7917

Adhesives and Removers



7760
▶ Universal Remover Wipes
▶ For Adhesives and Barriers
HCPCS: A4365

12 Cartons of
50 Wipes Each

Stock No
7760



7730

Box of 4

- ▶ Medical Adhesive
- ▶ Improves adhesive contact between skin and barrier

HCPCS: A4364

Stock No

Description

7730 3.2 oz. (90g) spray can



7731

Box of 4

- ▶ Medical Adhesive Remover

HCPCS: A4455

Stock No

Description

7731 2.7 oz. (76g) spray can

Odor Eliminators and Lubricants



78500/78501

- ▶ Adapt Lubricating Deodorant

HCPCS Code: A4394

Stock No

Description

Quantity

78500	8 oz. (236mL) bottle	Box of 1
78501	0.27 oz. (8mL) packet	50 Packets/Box



773x

- ▶ m9 Odor Eliminator Spray

HCPCS: No Code

Stock No (unscented)

Description

Quantity

7732	2 oz. (60mL) pump spray	Box of 12
7733	8 oz. (240mL) pump spray	Box of 6

Stock No (scented)

Description

Quantity

7734	2 oz. (60mL) pump spray	Box of 12
7735	8 oz. (240mL) pump spray	Box of 6

Cleaner/Decrystallizer



7736

Box of 12

- ▶ m9 Cleaner/Decrystallizer
- ▶ Cleans Urinary Drainage Systems
- ▶ pH balanced and non-acidic
- ▶ Free wash bottle included

HCPCS: A5131

Stock No

Description

7736 16 oz. (480mL) bottle

Ostomy Belts



7300/7299

Box of 10

- ▶ Ostomy Belts

HCPCS: A4367

Stock No

Size

7300	26"-43" (66 - 109cm) Medium
7299	29"-49" (74 - 125cm) Large

Drainable Pouch Clamp



8770

Box of 20

- ▶ Drainable Pouch Clamp

HCPCS: A4363

Stock No

Color

8770 Beige

Urostomy Drain Tube Adapter



7331

Box of 10

- ▶ Urostomy Drain Tube Adapter

HCPCS: A4421

Stock No

7331

CenterPointLock Two-Piece Pouching System

Letter Matching System

To help you select compatible CenterPointLock Two-Piece System components, every CenterPointLock System product box label features a letter designation: Ⓜ, Ⓛ, Ⓧ, Ⓚ or Ⓦ. Once you've selected a CenterPointLock Two-Piece System skin barrier, select a pouch, irrigator drain or stoma cap with a matching letter on the box label.

Flat Skin Barriers with Floating Flanges



872x (CUT-TO-FIT) **Box of 5**
 ▶ Flexend Skin Barrier, Flat, Tape Border

HCPCS: A4409

Stock No	Letter Match	Flange Size	Barrier Opening
8727	Ⓜ	1-1/2" (38mm)	up to 1" (25mm)
8722	Ⓛ	1-3/4" (44mm)	up to 1-1/4" (32mm)
8723	Ⓧ	2-1/4" (57mm)	up to 1-3/4" (44mm)
8724	Ⓚ	2-3/4" (70mm)	up to 2-1/4" (57mm)



875x (CUT-TO-FIT) **Box of 5**
 ▶ Flexend Skin Barrier, Flat, Without Tape Border

HCPCS FOR 8752, 8753 & 8757: A4409

HCPCS FOR 8754: A4410

Stock No	Letter Match	Flange Size	Barrier Opening
8752	Ⓛ	1-3/4" (44mm)	up to 1-1/4" (32mm)
8753	Ⓧ	2-1/4" (57mm)	up to 1-3/4" (44mm)
8754	Ⓚ	2-3/4" (70mm)	up to 2-1/4" (57mm)



372x (CUT-TO-FIT) **Box of 5**
 ▶ FlexWear Skin Barrier, Flat, Tape Border

HCPCS FOR 3722-3724 & 3727: A4414

HCPCS FOR 3726: A4415

Stock No	Letter Match	Flange Size	Barrier Opening
3727	Ⓜ	1-1/2" (38mm)	up to 1" (25mm)
3722	Ⓛ	1-3/4" (44mm)	up to 1-1/4" (32mm)
3723	Ⓧ	2-1/4" (57mm)	up to 1-3/4" (44mm)
3724	Ⓚ	2-3/4" (70mm)	up to 2-1/4" (57mm)
3726	Ⓦ	4" (102mm)	up to 3-1/2" (89mm)



376x (CUT-TO-FIT) **Box of 5**
 ▶ SoftFlex Skin Barrier, Flat, Without Tape Border

HCPCS FOR 3762-3764 & 3767: A4414

HCPCS FOR 3766: A4415

Stock No	Letter Match	Flange Size	Barrier Opening
3767	Ⓜ	1-1/2" (38mm)	up to 1" (25mm)
3762	Ⓛ	1-3/4" (44mm)	up to 1-1/4" (32mm)
3763	Ⓧ	2-1/4" (57mm)	up to 1-3/4" (44mm)
3764	Ⓚ	2-3/4" (70mm)	up to 2-1/4" (57mm)
3766	Ⓦ	4" (102mm)	up to 3-1/2" (89mm)



370x (CUT-TO-FIT) **Box of 5**
 ▶ HoliHesive Skin Barrier, Flat, Tape Border

HCPCS FOR 3702-3704 & 3707: A4414

HCPCS FOR 3706: A4415

Stock No	Letter Match	Flange Size	Barrier Opening
3707	Ⓜ	1-1/2" (38mm)	up to 3/4" (19mm)
3702	Ⓛ	1-3/4" (44mm)	up to 1" (25mm)
3703	Ⓧ	2-1/4" (57mm)	up to 1-1/2" (38mm)
3704	Ⓚ	2-3/4" (70mm)	up to 2" (51mm)
3706	Ⓦ	4" (102mm)	up to 3-1/2" (89mm)

Flat Skin Barriers with Floating Flanges (continued)



874x (PRE-SIZED)

Box of 5

► Flexextend Skin Barrier, Flat, Tape Border

HCPCS: A4409

Stock No	Letter Match	Flange Size	Barrier Opening
8743	Ⓜ	1-1/2" (38mm)	7/8" (22mm)
8744	Ⓜ	1-1/2" (38mm)	1" (25mm)
8745	Ⓛ	1-3/4" (44mm)	1-1/8" (29mm)
8746	Ⓛ	1-3/4" (44mm)	1-1/4" (32mm)
8747	Ⓝ	2-1/4" (57mm)	1-3/8" (35mm)
8748	Ⓝ	2-1/4" (57mm)	1-1/2" (38mm)



374x (PRE-SIZED)

Box of 5

► FlexWear Skin Barrier, Flat, Tape Border

HCPCS: A4414

Stock No	Letter Match	Flange Size	Barrier Opening
3747	Ⓜ	1-1/2" (38mm)	3/4" (19mm)
3742	Ⓜ	1-1/2" (38mm)	1" (25mm)
3748	Ⓛ	1-3/4" (44mm)	1-1/4" (32mm)
3743	Ⓝ	2-1/4" (57mm)	1-1/2" (38mm)
3749	Ⓝ	2-1/4" (57mm)	1-3/4" (44mm)

Convex Skin Barriers with Floating Flanges



877x (CUT-TO-FIT)

Box of 5

► Flexextend Skin Barrier, Convex, Tape Border

HCPCS: A4407

Stock No	Letter Match	Flange Size	Barrier Opening
8774	Ⓛ	1-3/4" (44mm)	up to 1" (25mm)
8778	Ⓝ	2-1/4" (57mm)	up to 1-1/2" (38mm)
87711	Ⓚ	2-3/4" (70mm)	up to 2" (51mm)



873x (PRE-SIZED)

Box of 5

► Flexextend Skin Barrier, Convex, Tape Border

HCPCS: A4407

Stock No	Letter Match	Flange Size	Barrier Opening
8731	Ⓜ	1-1/2" (38mm)	5/8" (15mm)
8732	Ⓜ	1-1/2" (38mm)	3/4" (19mm)
8733	Ⓛ	1-3/4" (44mm)	7/8" (22mm)
8734	Ⓛ	1-3/4" (44mm)	1" (25mm)
8735	Ⓝ	2-1/4" (57mm)	1-1/8" (29mm)
8736	Ⓝ	2-1/4" (57mm)	1-1/4" (32mm)
8737	Ⓝ	2-1/4" (57mm)	1-3/8" (35mm)
8738	Ⓝ	2-1/4" (57mm)	1-1/2" (38mm)



373x (PRE-SIZED)

Box of 5

► FlexWear Skin Barrier, Convex, Tape Border

HCPCS: A4373

Stock No	Letter Match	Flange Size	Barrier Opening
3730	Ⓜ	1-1/2" (38mm)	1/2" (13mm)
3731	Ⓜ	1-1/2" (38mm)	5/8" (15mm)
3732	Ⓜ	1-1/2" (38mm)	3/4" (19mm)
3733	Ⓛ	1-3/4" (44mm)	7/8" (22mm)
3734	Ⓛ	1-3/4" (44mm)	1" (25mm)
3735	Ⓝ	2-1/4" (57mm)	1-1/8" (29mm)
3736	Ⓝ	2-1/4" (57mm)	1-1/4" (32mm)
3737	Ⓝ	2-1/4" (57mm)	1-3/8" (35mm)
3738	Ⓝ	2-1/4" (57mm)	1-1/2" (38mm)
3739	Ⓚ	2-3/4" (70mm)	1-5/8" (41mm)
37310	Ⓚ	2-3/4" (70mm)	1-3/4" (44mm)
37311	Ⓚ	2-3/4" (70mm)	2" (51mm)



CenterPointLock Two-Piece Pouching System

Drainable Pouches



389x/383x

Box of 10

- ▶ Drainable Pouch, Replaceable Filter, Clamp Closure, Odor-Barrier Rustle-Free Film, Belt Tabs, Without ComfortWear Panels. Includes 1 Bag of 20 Hollister Replacement Filters and 2 Plugs, #7766

HPCPS: A4425, A4368

Transparent Stock No	Beige Stock No	Letter Match	Flange Size	Pouch Length
3897	3837	Ⓜ	1-1/2" (38mm)	12" (30cm)
3892	3832	Ⓛ	1-3/4" (44mm)	12" (30cm)
3893	3833	Ⓢ	2-1/4" (57mm)	12" (30cm)
3894	3834	Ⓚ	2-3/4" (70mm)	12" (30cm)



881x

Box of 10

- ▶ Drainable Pouch, Without Integrated Filter, Clamp Closure, Odor-Barrier Rustle-Free Film, Belt Tabs, ComfortWear Panels - Both Sides

HPCPS: A5063

Beige Stock No	Letter Match	Flange Size	Pouch Length
8817	Ⓜ	1-1/2" (38mm)	12" (30cm)
8812	Ⓛ	1-3/4" (44mm)	12" (30cm)
8813	Ⓢ	2-1/4" (57mm)	12" (30cm)
8814	Ⓚ	2-3/4" (70mm)	12" (30cm)



380x/381x

Box of 10

- ▶ Drainable Pouch, Without Integrated Filter, Clamp Closure, Odor-Barrier Rustle-Free Film, Belt Tabs, Without ComfortWear Panels

HPCPS: A5063

Transparent Stock No	Beige Stock No	Letter Match	Flange Size	Pouch Length
3807	3817	Ⓜ	1-1/2" (38mm)	11" (28cm)
3802	3812	Ⓛ	1-3/4" (44mm)	11" (28cm)
3803	3813	Ⓢ	2-1/4" (57mm)	11" (28cm)
3804	3814	Ⓚ	2-3/4" (70mm)	11" (28cm)
3806	-	Ⓛ	4" (102mm)	11" (28cm)



387x

Box of 10

- ▶ Drainable Pouch, Without Integrated Filter, Clamp Closure, Odor-Barrier Film, Belt Tabs Without ComfortWear Panels

HPCPS: A4412

Transparent Stock No	Letter Match	Flange Size	Pouch Length
3872	Ⓛ	1-3/4" (44mm)	16" (40cm)
3873	Ⓢ	2-1/4" (57mm)	16" (40cm)
3874	Ⓚ	2-3/4" (70mm)	16" (40cm)



384x

Box of 10

- ▶ Drainable Mini-Pouch, Without Integrated Filter, Clamp Closure, Odor-Barrier Rustle-Free Film, Belt Tabs, Without ComfortWear Panels

HPCPS: A5063

Beige Stock No	Letter Match	Flange Size	Pouch Length
3847	Ⓜ	1-1/2" (38mm)	9" (23cm)
3842	Ⓛ	1-3/4" (44mm)	9" (23cm)
3843	Ⓢ	2-1/4" (57mm)	9" (23cm)
3844	Ⓚ	2-3/4" (70mm)	9" (23cm)

Closed Pouches



834x

Box of 15

- ▶ Closed Pouch, Integrated Filter, Odor-Barrier Rustle-Free Film, Belt Tabs, ComfortWear Panels - Both Sides

HCPCS: A4419

Beige Stock No	Letter Match	Flange Size	Pouch Length
8347	Ⓜ	1-1/2" (38mm)	7-1/2" (19cm)
8342	Ⓛ	1-3/4" (44mm)	7-1/2" (19cm)
8343	Ⓝ	2-1/4" (57mm)	7-1/2" (19cm)
8344	Ⓚ	2-3/4" (70mm)	7-1/2" (19cm)



334x

Box of 15

- ▶ Closed Pouch, Integrated Filter, Odor-Barrier Rustle-Free Film, Belt Tabs, Without ComfortWear Panels

HCPCS: A4419

Beige Stock No	Letter Match	Flange Size	Pouch Length
3347	Ⓜ	1-1/2" (38mm)	9" (23cm)
3342	Ⓛ	1-3/4" (44mm)	9" (23cm)
3343	Ⓝ	2-1/4" (57mm)	9" (23cm)
3344	Ⓚ	2-3/4" (70mm)	9" (23cm)

Urostomy Pouches



890x

Box of 10

- ▶ Urostomy Pouch, Odor-Barrier Rustle-Free Film, Belt Tabs, ComfortWear Panels - Both Sides

HCPCS: A4432

Beige Stock No	Letter Match	Flange Size	Pouch Length
8902	Ⓛ	1-3/4" (44mm)	9" (23cm)
8903	Ⓝ	2-1/4" (57mm)	9" (23cm)



390x

Box of 10

- ▶ Urostomy Pouch, Odor-Barrier Rustle-Free Film, Belt Tabs, Without ComfortWear Panels

HCPCS: A4432

Transparent Stock No	Letter Match	Flange Size	Pouch Length
3907	Ⓜ	1-1/2" (38mm)	9" (23cm)
3902	Ⓛ	1-3/4" (44mm)	9" (23cm)
3903	Ⓝ	2-1/4" (57mm)	9" (23cm)
3904	Ⓚ	2-3/4" (70mm)	9" (23cm)

Stoma Cap



340x

Box of 25

- ▶ Stoma Cap, Integrated Filter, Odor-Barrier Rustle-Free Film Without ComfortWear Panels

HCPCS: A5055

Beige Stock No	Letter Match	Flange Size	Pouch Length
3407	Ⓜ	1-1/2" (38mm)	4-1/2" (11cm)
3402	Ⓛ	1-3/4" (44mm)	4-1/2" (11cm)
3403	Ⓝ	2-1/4" (57mm)	4-1/2" (11cm)
3404	Ⓚ	2-3/4" (70mm)	4-1/2" (11cm)

CenterPointLock Two-Piece Pouching System

Irrigator Sleeve



382x

▶ Irrigator Sleeve, Belt Tabs

Box of 5

HCPCS: A4397

Transparent Stock No	Letter Match	Flange Size	Pouch Length
3827	Ⓧ	1-1/2" (38mm)	30" (76cm)
3822	Ⓛ	1-3/4" (44mm)	30" (76cm)
3823	Ⓢ	2-1/4" (57mm)	30" (76cm)
3824	Ⓚ	2-3/4" (70mm)	30" (76cm)
3826	Ⓛ	4" (102mm)	30" (76cm)

*For Two-Piece Pouching System Kits,
refer to pages 27 -29.*

Karaya 5 One-Piece Pouching System

Drainable Pouches



366x (PRE-SIZED)

Box of 10

- ▶ Drainable Pouch, Karaya 5 Skin Barrier, Tape Border, Replaceable Filter, Pouch Clamp Closure, Odor-Barrier Rustle-Free Film, Belt Tabs, Without ComfortWear Panels, Includes 1 Bag of 20 Hollister Replacement Filters and 2 Plugs, #7766

HPCPS: A4424, A4368

Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3662	1" (25mm)	5/8" (15mm)	12" (30cm)
3668	1-1/4" (32mm)	7/8" (22mm)	12" (30cm)
3663	1-1/2" (38mm)	1-1/8" (29mm)	12" (30cm)
3669	1-3/4" (44mm)	1-3/8" (35mm)	12" (30cm)
3664	2" (51mm)	1-1/2" (38mm)	12" (30cm)
3665	2-1/2" (64mm)	2" (51mm)	12" (30cm)



360x (PRE-SIZED)

Box of 10

- ▶ Drainable Pouch, Karaya 5 Skin Barrier, Tape Border, Without Integrated Filter, Pouch Clamp Closure, Odor-Barrier Rustle-Free Film, Belt Tabs, Without ComfortWear Panels

HPCPS: A5061

Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3608	1-1/4" (32mm)	7/8" (22mm)	12" (30cm)
3603	1-1/2" (38mm)	1-1/8" (29mm)	12" (30cm)
3609	1-3/4" (44mm)	1-3/8" (35mm)	12" (30cm)
3604	2" (51mm)	1-1/2" (38mm)	12" (30cm)
3605	2-1/2" (64mm)	2" (51mm)	12" (30cm)
3606	3" (76mm)	2-1/2" (64mm)	12" (30cm)



323x (PRE-SIZED)

Box of 30

- ▶ Drainable Mini-Pouch, Karaya 5 Skin Barrier, Tape Border, Without Integrated Filter, Pouch Clamp Closure, Odor-Barrier Film, Belt Tabs, Without ComfortWear Panels

HPCPS: A5061

Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3232	1" (25mm)	5/8" (15mm)	9" (23cm)
3238	1-1/4" (32mm)	7/8" (22mm)	9" (23cm)
3233	1-1/2" (38mm)	1-1/8" (29mm)	9" (23cm)
3239	1-3/4" (44mm)	1-3/8" (35mm)	9" (23cm)
3234	2" (51mm)	1-1/2" (38mm)	9" (23cm)



322x (PRE-SIZED)

Box of 30

- ▶ Drainable Pouch, Karaya 5 Skin Barrier, Tape Border, Without Integrated Filter, Pouch Clamp Closure, Odor-Barrier Film, Belt Tabs, Without ComfortWear Panels

HPCPS: A5061

Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3222	1" (25mm)	5/8" (15mm)	12" (30cm)
3228	1-1/4" (32mm)	7/8" (22mm)	12" (30cm)
3223	1-1/2" (38mm)	1-1/8" (29mm)	12" (30cm)
3229	1-3/4" (44mm)	1-3/8" (35mm)	12" (30cm)
3224	2" (51mm)	1-1/2" (38mm)	12" (30cm)
3225	2-1/2" (64mm)	2" (51mm)	12" (30cm)
3226	3" (76mm)	2-1/2" (64mm)	12" (30cm)

Karaya 5 One-Piece Pouching System

Drainable Pouches (continued)



327x (PRE-SIZED)

Box of 30

- ▶ Drainable Pouch, Karaya 5 Skin Barrier, Tape Border, Without Integrated Filter, Pouch Clamp Closure, Odor-Barrier Film, Belt Tabs, Without ComfortWear Panels

HPCPS: A5061

Transparent Stock No	Pouch Barrier Size Opening	Pouch Length	
3278	1-1/4" (32mm)	7/8" (22mm)	16" (40cm)
3273	1-1/2" (38mm)	1-1/8" (29mm)	16" (40cm)
3279	1-3/4" (44mm)	1-3/8" (35mm)	16" (40cm)
3274	2" (51mm)	1-1/2" (38mm)	16" (40cm)
3275	2-1/2" (64mm)	2" (51mm)	16" (40cm)
3276	3" (76mm)	2-1/2" (64mm)	16" (40cm)

Closed Pouches



332x (PRE-SIZED)

Box of 30

- ▶ Closed Pouch, Karaya 5 Skin Barrier, Tape Border, Integrated Filter, Odor-Barrier Film, Belt Tabs, Without ComfortWear Panels

HPCPS: A4416

Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3328	1-1/4" (32mm)	7/8" (22mm)	9" (23cm)
3323	1-1/2" (38mm)	1-1/8" (29mm)	9" (23cm)
3329	1-3/4" (44mm)	1-3/8" (35mm)	9" (23cm)
3324	2" (51mm)	1-1/2" (38mm)	9" (23cm)
3325	2-1/2" (64mm)	2" (51mm)	9" (23cm)
3326	3" (76mm)	2-1/2" (64mm)	9" (23cm)



716x (PRE-SIZED)

Box of 30

- ▶ Closed Pouch, Karaya 5 Skin Barrier, Without Tape Border, Without Integrated Filter, Odor-Barrier Film, Belt Tabs, Without ComfortWear Panels

HPCPS: A5051

Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
7168	1-1/4" (32mm)	7/8" (22mm)	9" (23cm)
7163	1-1/2" (38mm)	1-1/8" (29mm)	9" (23cm)
7169	1-3/4" (44mm)	1-3/8" (35mm)	9" (23cm)
7164	2" (51mm)	1-1/2" (38mm)	9" (23cm)
7165	2-1/2" (64mm)	2" (51mm)	9" (23cm)
7166	3" (76mm)	2-1/2" (64mm)	9" (23cm)

Urostomy Pouches



143x (PRE-SIZED)

Box of 10

- ▶ Urostomy Pouch, Karaya 5 Skin Barrier, Tape Border, Odor-Barrier Film, Belt Tabs, Without ComfortWear Panels

HPCPS: A4431

Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
1432	1" (25mm)	5/8" (16mm)	9" (23cm)
1438	1-1/4" (32mm)	7/8" (22mm)	9" (23cm)
1433	1-1/2" (38mm)	1-1/8" (29mm)	9" (23cm)
1439	1-3/4" (44mm)	1-3/8" (35mm)	9" (23cm)
1434	2" (51mm)	1-1/2" (38mm)	9" (23cm)

Other One-Piece Products

Closed Pouches



314x (PRE-SIZED)

Box of 30

- ▶ Closed Mini-Pouch, Tape Adhesive, Integrated Filter, Odor-Barrier Film, Without ComfortWear Panels

HCPCS: A4418

Transparent Stock No	Barrier Opening	Pouch Length
3142	1" (25mm)	7" (18cm)
3148	1-1/4" (32mm)	7" (18cm)
3143	1-1/2" (38mm)	7" (18cm)
3149	1-3/4" (44mm)	7" (18cm)
3144	2" (51mm)	7" (18cm)

Stoma Cap



318x (PRE-SIZED)

Box of 30

- ▶ Stoma Cap, Tape Adhesive, Integrated Filter, Odor-Barrier Film, Without ComfortWear Panels

HCPCS: A5055

White Stock No	Opening	Size
3184	2" (51mm)	4-1/4" (11cm)
3186	3" (76mm)	4-1/4" (11cm)

Ostomy Accessories

Skin Care Products



517210

▶ Restore Skin Cleanser

HCPCS: No Code

Box of 12

Stock No

Description

517210

8 oz. (237mL) pump spray



517220

▶ Restore Skin Conditioning Creme

HCPCS: A6250

Box of 12


Stock No

Description

517220

4 oz. (120mL) bottle

Irrigation Products



7721

▶ Cone Irrigator Kit includes:

- ▶ (1) Irrigator
- ▶ (1) Stoma Cone with Connector
- ▶ (1) Cleaning Brush

HCPCS: A4398, A4399

Box 1

Stock No

7721



7723

▶ Stoma Cone with Connector

▶ Replacement Unit

HCPCS: A4399

Box of 10

Stock No

7723



7728/7724

▶ Irrigator Sleeve with Belt Tabs

HCPCS: A4397

Box of 20

Stock No

Opening Size

Length

7728

2" (51mm)

22" (56cm)

7724

3" (76mm)

22" (56cm)



7740

▶ Adapt Stoma Lubricant

HCPCS: A4402

Box of 10

Stock No

Description

7740

4 oz. (120mL) bottle

Loop Ostomy Bridge



7767

▶ Loop Ostomy Bridge - Sterile

HCPCS: A4421

Box of 20

Stock No

7767

Replacement Filters



7766

Box of 5 Packs

- ▶ Hollister Replacement Filters
- ▶ For use with 366x, 383x and 389x
- ▶ 1 pack contains 20 Filters and 2 Plugs

HCPCS: A4368

Stock No

7766

Bedside Drainage Collection System with Anti-Reflux Valve



9839

Box of 20

- ▶ Latex-Free
- ▶ Sample Port
- ▶ Sterile Fluid Pathway

HCPCS: A4357

Stock No

9839

Description

68 oz. (2000ml) Capacity

Skin Barrier Pastes and Powders



7910

Box of 12

- ▶ Karaya Paste

HCPCS: A4405

Stock No

7910

Description

4.5 oz. (135ml) tube



7905

Box of 12

- ▶ Karaya Powder

HCPCS: A4371

Stock No

7905

Description

2.5 oz. (75ml) puff bottle

Skin Barriers



7800

Box of 5

- ▶ Premium (standard wear) Skin Barrier

HCPCS: A4362

Stock No

7800

Size

4" x 4" (10cm x 10cm)



7700

Box of 5

- ▶ HollisHesive (standard wear) Skin Barrier

HCPCS: A4362

Stock No

7700

Size

4" x 4" (10cm x 10cm)

Odor Eliminators



771x

- ▶ m9 Odor Eliminator Drops

HCPCS: A4394

Stock No

7715

Description

1 oz. (30ml) bottle

7717

8 oz. (240ml) bottle

Quantity

Box of 12

Box of 6

Discontinued Products

FirstChoice to Premier — Discontinued Date 12-31-06

FirstChoice One-Piece Cut-to-Fit Drainable Pouch →

FlexWear Standard Wear, without Filter, with Tape, Pouch Clamp Closure

Stock No	Barrier Opening	Length	Color
3631	up to 2-1/2" (64mm)	12"	Transparent

SoftFlex Standard Wear, without Filter, with Tape, Pouch Clamp Closure

Stock No	Barrier Opening	Length	Color
3630	up to 2-1/2" (64mm)	12"	Transparent

FirstChoice One-Piece Pre-Sized Drainable Pouch →

SoftFlex Standard Wear, without Filter, with Tape, Pouch Clamp Closure

Stock No	Barrier Opening	Length	Color
3632	1" (25mm)	12"	Transparent
3638	1-1/4" (32mm)	12"	Transparent
3633	1-1/2" (38mm)	12"	Transparent
3639	1-3/4" (44mm)	12"	Transparent
3634	2" (51mm)	12"	Transparent

FirstChoice One-Piece Convex Pre-Sized Drainable Pouch →

FlexWear Standard Wear, without Filter, with Tape, Pouch Clamp Closure

Stock No	Barrier Opening	Length	Color
3610	3/4" (19mm)	12"	Transparent
3611	7/8" (22mm)	12"	Transparent
3612	1" (25mm)	12"	Transparent
3613	1-1/8" (29mm)	12"	Transparent
3614	1-1/4" (32mm)	12"	Transparent
3615	1-3/8" (35mm)	12"	Transparent
3616	1-1/2" (38mm)	12"	Transparent

3617	1-5/8" (41mm)	12"	Transparent
3680	3/4" (19mm)	12"	Beige
3681	7/8" (22mm)	12"	Beige
3682	1" (25mm)	12"	Beige
3683	1-1/8" (29mm)	12"	Beige
3684	1-1/4" (32mm)	12"	Beige

FirstChoice One-Piece Cut-to-Fit Closed Pouch →

SoftFlex Standard Wear, with Filter, without Tape

Stock No	Barrier Opening	Length	Color
3520	up to 2-1/2" (64mm)	7.5"	Transparent
3450	up to 2-1/2" (64mm)	7.5"	Beige

FirstChoice One-Piece Pre-Sized Closed Pouch →

SoftFlex Standard Wear, with Filter, without Tape

Stock No	Barrier Opening	Length	Color
3452	1" (25mm)	7.5"	Beige
3458	1-1/4" (32mm)	7.5"	Beige
3453	1-1/2" (38mm)	7.5"	Beige
3459	1-3/4" (44mm)	7.5"	Beige
3454	2" (51mm)	7.5"	Beige
3455	2-1/2" (64mm)	7.5"	Beige

FirstChoice One-Piece Cut-to-Fit Urostomy Pouch →

FlexWear Standard Wear, with Tape

Stock No	Barrier Opening	Length	Color
1460	up to 2-1/2" (64mm)	9"	Transparent

Premier One-Piece Cut-to-Fit Drainable Pouch

SoftFlex Standard Wear Barrier, without Filter, without Tape

Lock 'n Roll	Pouch Clamp	Barrier Opening	Length	Color	Comments
88401	88501	up to 2-1/2" (64mm)	12"	Ultra-Clear	

SoftFlex Standard Wear Barrier, without Filter, without Tape

Lock 'n Roll	Pouch Clamp	Barrier Opening	Length	Color	Comments
88401	88501	up to 2-1/2" (64mm)	12"	Ultra-Clear	

Premier One-Piece Pre-Sized Drainable Pouch

Flexextend Extended Wear Barrier, without Filter, with Tape

Lock 'n Roll	Pouch Clamp	Barrier Opening	Length	Color	Comments
8532	8632	1" (25mm)	12"	Transparent	Change in skin barrier
8538	8638	1-1/4" (32mm)	12"	Transparent	Change in skin barrier
8533	8633	1-1/2" (38mm)	12"	Transparent	Change in skin barrier
8539	8639	1-3/4" (44mm)	12"	Transparent	Change in skin barrier
8531*	8634	2" (51mm)	12"	Transparent	Change in skin barrier, *Change to Cut-to-Fit

Premier One-Piece Convex Pre-Sized Drainable Pouch

Flexextend Extended Wear Barrier, without Filter, with Tape

Lock 'n Roll	Pouch Clamp	Barrier Opening	Length	Color	Comments
8510	8610	3/4" (19mm)	12"	Transparent	Change in skin barrier
8511	8611	7/8" (22mm)	12"	Transparent	Change in skin barrier
8512	8612	1" (25mm)	12"	Transparent	Change in skin barrier
8513	8613	1-1/8" (29mm)	12"	Transparent	Change in skin barrier
8514	8614	1-1/4" (32mm)	12"	Transparent	Change in skin barrier
8515	8615	1-3/8" (35mm)	12"	Transparent	Change in skin barrier
8528*	8616	1-1/2" (38mm)	12"	Transparent	Change in skin barrier, *Change to Cut-to-Fit
85211	8617	1-5/8" (41mm)	12"	Transparent	Change in skin barrier
8590	n/a	3/4" (19mm)	12"	Beige	Change in skin barrier
8591	n/a	7/8" (22mm)	12"	Beige	Change in skin barrier
8592	n/a	1" (25mm)	12"	Beige	Change in skin barrier
8593	n/a	1-1/8" (29mm)	12"	Beige	Change in skin barrier
8594	n/a	1-1/4" (32mm)	12"	Beige	Change in skin barrier

Premier One-Piece Cut-to-Fit Closed Mini Pouch

SoftFlex Standard Wear, with Filter, without Tape

Stock No	Barrier Opening	Length	Color	Comments
82500	up to 2-1/8" (55mm)	7"	Transparent	Change in Max Cut-to-Fit size
82100	up to 2-1/8" (55mm)	7"	Beige	Change in Max Cut-to-Fit size

Premier One-Piece Cut-to-Fit Closed Mini Pouch

SoftFlex Standard Wear, with Filter, without Tape

Stock No	Barrier Opening	Length	Color	Comments
82325	1" (25mm)	7"	Beige	
82300	up to 2-1/8" (55mm)	7"	Beige	Change to Cut-to-Fit
82300	up to 2-1/8" (55mm)	7"	Beige	Change to Cut-to-Fit
82300	up to 2-1/8" (55mm)	7"	Beige	Change to Cut-to-Fit
82300	up to 2-1/8" (55mm)	7"	Beige	Change to Cut-to-Fit
82300	up to 2-1/8" (55mm)	7"	Beige	Change to Cut-to-Fit

Premier One-Piece Cut-to-Fit Urostomy Pouch

Flexextend Extended Wear, with Tape

Stock No	Barrier Opening	Length	Color	Comments
8460	up to 2-1/2" (64mm)	9"	Transparent	Change in skin barrier

FirstChoice to Premier — Discontinued Date 12-31-06 (continued)

FirstChoice One-Piece Convex Pre-Sized Urostomy Pouch →

FlexWear Standard Wear, with Tape

Stock No	Barrier Opening	Length	Color
1481	5/8" (16mm)	9"	Transparent
1482	3/4" (19mm)	9"	Transparent
1483	7/8" (22mm)	9"	Transparent
1484	1" (25mm)	9"	Transparent
1485	1-1/8" (29mm)	9"	Transparent
1486	1-1/4" (32mm)	9"	Transparent
1487	1-3/8" (35mm)	9"	Transparent
1488	1-1/2" (38mm)	9"	Transparent

FirstChoice Drainable Kit - Non-Sterile →

FlexWear Standard Wear, with Tape, Pouch Clamp Closure

Stock No	Barrier Opening	Length	Color
3636	up to 2-1/2" (64mm)	12"	Transparent

FirstChoice Urostomy Kit - Non-Sterile →

FlexWear Standard Wear, with Tape

Stock No	Barrier Opening	Length	Color
1466	up to 2-1/2" (64mm)	9"	Transparent
3455	2-1/2" (64mm)	7.5"	Beige

Premier One-Piece Convex Pre-Sized Urostomy Pouch

FlexWear Extended Wear, with Tape

Stock No	Barrier Opening	Length	Color	Comments
8481	5/8" (16mm)	9"	Transparent	Change to extended wear barrier
8482	3/4" (19mm)	9"	Transparent	Change to extended wear barrier
8483	7/8" (22mm)	9"	Transparent	Change to extended wear barrier
8484	1" (25mm)	9"	Transparent	Change to extended wear barrier
8485	1-1/8" (29mm)	9"	Transparent	Change to extended wear barrier
8486	1-1/4" (32mm)	9"	Transparent	Change to extended wear barrier
8487	1-3/8" (35mm)	9"	Transparent	Change to extended wear barrier
8488	1-1/2" (38mm)	9"	Transparent	Change to extended wear barrier

Premier Single Use Kit Drainable - Non-Sterile

FlexWear Extended Wear, with Tape

Lock 'n Roll	Pouch Clamp	Barrier Opening	Length	Color	Comments
89004	89003	up to 2-1/2" (64mm)	12"	Ultra-Clear	Change to extended wear

Premier Single Use Kit Urostomy - Non-Sterile

FlexWear Extended Wear, with Tape

Stock No	Barrier Opening	Length	Color	Comments
89005	up to 2-1/2" (64mm)	9"	Ultra-Clear	Change to extended wear
82300	up to 2-1/8" (55mm)	7"	Beige	Change to Cut-to-Fit

Pouchkins Pediatric Pouching Systems — Planned Discontinuation Date 12-31-07

One-Piece Drainable Products →

Stock No	Barrier Opening	Length	Color
9766	up to 2" (51mm)	8-1/2" (21cm)	Transparent
9768	up to 1-1/2" (38mm)	6-1/2" (16cm)	Ultra-Clear

One-Piece Urostomy Product →

Stock No	Barrier Opening	Length	Color
9767	up to 1-1/2" (38mm)	8-1/2" (21cm)	Transparent

Two-Piece Barrier →

Stock No	Flange Size	Barrier Opening
3770	1-3/4" (44mm)	up to 1-1/4" (32mm)

Two-Piece Pouch →

Stock No	Flange Size	Length	Type
3771	1-3/4" (44mm)	8-1/2" (21cm)	Drainable
3772	1-1/3" (44mm)	8-1/2" (21cm)	Urostomy

One-Piece Drainable Products

Stock No	Barrier Opening	Length	Color
3795	up to 2" (51mm)	7" (18cm)	Ultra-Clear
3796	up to 1-1/2" (38mm)	6-1/2" (16cm)	Transparent
3796	up to 1-1/2" (38mm)	6-1/2" (16cm)	Transparent

One-Piece Urostomy Product

Stock No	Barrier Opening	Length	Color
3797	up to 1-1/2" (38mm)	8-3/4" (22cm)	Transparent

Two-Piece Barrier

Stock No	Flange Size	Barrier Opening
3761	1-3/4" (44mm)	up to 1-1/4" (32mm)

Two-Piece Pouch

Stock No	Flange Size	Length	Type
3799	1-3/4" (44mm)	6-1/2" (16cm)	Drainable
18292	1-3/4" (44mm)	7" (18cm)	Drainable
18402	1-3/4" (44mm)	9" (23cm)	Urostomy

New Image Closed Pouch (Filter) — Planned Discontinuation Date 12-31-07

New Image Closed Pouch (Old Filter) →

Stock No	Flange Size	Length	Color
18302	1-3/4" (44mm)	7-1/2" (19cm)	Transparent
18304	2-3/4" (70mm)	7-1/2" (19cm)	Transparent
18312	1-3/4" (44mm)	7-1/2" (19cm)	Beige
18313	2-1/4" (57mm)	7-1/2" (19cm)	Beige
18314	2-3/4" (70mm)	7-1/2" (19cm)	Beige

New Image Closed Pouch (New Filter)

Stock No	Flange Size	Length	Color
18362	1-3/4" (44mm)	9" (23cm)	Transparent
18364	2-3/4" (70mm)	9" (23cm)	Transparent
18372	1-3/4" (44mm)	9" (23cm)	Beige
18373	2-1/4" (44mm)	9" (23cm)	Beige
18374	2-3/4" (70mm)	9" (23cm)	Beige

Odor-Barrier Pouch Odor-Barrier Pouch is highly resistant in allowing odors to pass through the pouch film. The Odor-Barrier quality increases confidence when wearing a pouch.

Pre-Sized Skin Barrier A skin barrier that is already cut to the proper size if stoma is round.

Rustle-Free Pouch Rustle-Free Pouch provides for ultimate discretion since the pouch does not rustle (make noise) during normal wear.

Skin Barrier Portion of the pouching system that protects the skin from contact with the discharge. It also provides an adhesive attachment for the pouching system to the body. Examples include FlexTend, FlexWear and SoftFlex.

SoftFlex Skin Barrier SoftFlex is a synthetic skin barrier, which provides gentle adhesive attachment to the skin. It's designed to limit skin damage with frequent removal and can be changed daily (or more often).

Standard Wear Skin Barrier Provides a gentle adhesive attachment to the skin. These skin barriers absorb less and disintegrate faster than extended wear skin barriers. Examples include FlexWear and SoftFlex.

Stoma The place where the drainage comes out of the body when you have an ostomy. The stoma is red and moist and doesn't have feeling.

Urostomy An ostomy (surgical opening) created to drain urine.

Urostomy Pouch Pouch with a drain valve at the bottom so urine can be emptied quickly and easily.

Wear Time Wear time is the length of time a product can be worn before failure. Wear time varies widely based upon a variety of factors.

Closed Pouch A pouch without a spout or a clamp. Quick and easy to discard. It must be removed to be emptied.

Colostomy An ostomy (surgical opening) created in the colon, part of the large intestine.

ComfortWear Pouch Panels Lining over the plastic of the pouch to enhance comfort and discretion. The soft pouch panels may be on one or both sides of the pouch.

Convex Skin Barrier May be used when the stoma does not protrude from the abdomen. The skin barrier is shaped to press more deeply around the stoma.

Cut-to-Fit Skin Barrier If the stoma is not round, measure where it is the largest. Select a skin barrier that can be cut to the right size and shape.

Drainable Pouch A pouch with an opening at the bottom. An integrated closure or clamp is used to keep the pouch closed until it's time to empty it.

Extended Wear Skin Barriers Skin barriers with special additives that achieve a stronger adhesive seal and are more resistant to breakdown. Skin barrier of choice for urostomies and ileostomies. FlexTend is an example. It is not recommended to use skin-prepping agents under extended wear skin barriers.

Filter Lets gas but not odor out of pouch.

Flange The plastic ring in a two-piece ostomy system that is used to connect the two parts together.

Flat Skin Barrier May be used when the stoma does not protrude from the abdomen.

FlexTend Skin Barrier FlexTend is an extended wear, synthetic skin barrier, which is designed to be the most resistant to discharge.

FlexWear Skin Barrier FlexWear is a standard wear, synthetic skin barrier, which is designed to be worn more than one day. FlexWear barrier is more resistant to breakdown than a SoftFlex skin barrier.

Floating Flange Allows for fingers to be put under the flange for support while attaching the pouch to the barrier.

Ileostomy An ostomy (surgical opening) created in the small intestine.

Irrigation A procedure used to control the time of elimination from colostomy. It is not for everyone and requires special equipment and teaching.

Karaya Skin Barrier Karaya is a natural barrier used for sensitive skin or in some cases, where there is sensitivity to synthetic products.

Lock 'n Roll Closure Patented system for closing the bottom of drainable pouches. Utilizes plastic interlocking fasteners to create a secure closure. Found on New Image and Premier drainable pouches.
