

What goes into our new
CalciCare calcium alginate dressing
helps ensure all the dressing comes out.

We changed the name and reinforced the benefits.

- New strengthening layer provides improved wet strength
 - Excellent absorption capacity
- Same natural calcium alginate dressing

Available with Silver for additional benefits:


- Effective against a broad spectrum of microorganisms through recommended wear time
 - Odor reduction

CalciCare

 **Hollister**

New strengthening layer increases confidence that what goes in comes out.

CalciCare calcium alginate dressings—addressing your wound packing needs.



Strengthening nylon web

- Increased integrity when wet, helping to achieve clean, one-piece removal

Calcium alginate material

- Highly absorbent
- Natural calcium alginate fibers
- Creates conformable protective gel when in contact with exudate

CalciCare calcium alginate dressing

- Natural calcium alginate fibers**
 - Residual fibers are biodegradable
- Protects the wound**
 - Aids autolytic debridement, facilitating wound healing
- High wet strength**
 - Facilitates intact, one-piece removal
- Hemostatic properties**
 - As the product contains alginate, it may assist in supporting the control of minor bleeding in superficial wounds
- Highly absorbent**
 - Uniformly absorbs exudates and protects the peri-wound skin
- Helps form and maintain a moist wound interface**
 - Creates a conformable protective gel when in contact with exudate


CalciCare calcium alginate silver dressing

All the benefits with the addition of sustained silver release


- Silver ions are released in the presence of wound exudate
- Effective against a broad spectrum of microorganisms through recommended wear time

New strengthening layer increases confidence that **what goes in comes out.**

Excellent wet strength¹



Comparable absorbency¹



CalciCare



CalciCare calcium alginate silver dressing

Stock No	Product Size	Quantity/Box	HCPGS
529967R	2" x 2" (5 cm x 5 cm)	10	A6196
529968R	4" x 4.75" (10.2 cm x 12 cm)	10	A6197
529969R	1" x 18" (2.5 cm x 46 cm)	5	A6199



CalciCare calcium alginate dressing

Stock No	Product Size	Quantity/Box	HCPGS
529938R	2" x 2" (5 cm x 5 cm)	10	A6196
529937R	4" x 4.75" (10.2 cm x 12 cm)	10	A6197
529939R	4" x 8" (10.2 cm x 20.3 cm)	5	A6197
529940R	1" x 18" (2.5 cm x 46 cm)	5	A6199



Restore foam with adhesive border

Stock No	Product Size	Quantity/Box
509387 Pad	4" x 4" (10 cm x 10 cm) 2" x 2" (5 cm x 5 cm)	10
509388 Pad	6" x 6" (15 cm x 15 cm) 4" x 4" (10 cm x 10 cm)	10
509389 Pad	6" x 8" (15 cm x 20 cm) 3.3" x 5.8" (8.5 cm x 14.7 cm)	10
509390 Pad	8" x 8" (20 cm x 20 cm) 5" x 5" (14 cm x 14 cm)	10



Restore foam with TRIACT technology

Stock No	Product Size	Quantity/Box
509380	2.5" x 2.5" (6 cm x 6 cm)	10
509381	4" x 4" (10 cm x 10 cm)	10
509382	6" x 6" (15 cm x 15 cm)	10

1. Data on file - LWIDEV030 2014.

Manufactured for:
Hollister Incorporated
 2000 Hollister Drive
 Libertyville, IL 60048 USA

Distributed in Canada by:
 Distribué au Canada par:
Hollister Limited
 95 Mary Street
 Aurora, Ontario L4G 1G3

www.hollisterwoundcare.com

Rx Caution: For dressings which contain silver, Federal (USA) law restricts this device to sale by or on the order of a physician (or properly licensed practitioner). Prior to use, be sure to read the entire Instructions For Use package insert supplied with the product for Device Intended Use, Description, Contraindications, Warnings, Precautions, Adverse events, and Instructions For Use.

For product questions, sampling needs, or detailed clinical questions concerning our products in the US call **1.888.740.8999**. In Canada call **1.800.263.7400**.

The Hollister logo, CalciCare, and Restore logos are trademarks of Hollister Incorporated.

©2014 Hollister Incorporated.

Made and printed in the USA.

922571

