



Side-loaded catheter



Integrated 1500 mL collection bag with tear strip
allows touch free catheterization



Smooth eyelets
help provide comfort upon insertion and removal



Pre-lubricated catheter
top half lubrication reduces friction during insertion



Protective tip
helps the catheter bypass bacteria located within the first 15mm of the distal urethra

Trusted Quality. Reliable Performance. Cost Effective.

Sterile, and self-contained, the Hollister Apogee Plus intermittent catheter system empowers you to maintain your independence and live an active, mobile lifestyle.

Apogee Plus

Intermittent Catheter System



Continence Care

Apogee Plus

Intermittent Catheter System

Touch Free Protective Tip

- Helps protect the catheter tip from exposure to bacteria located within the first 15mm of the distal urethra
- Helps reduce the risk of carrying bacteria up into the urinary system
- Serves as a guide for proper and easy insertion

Side Loaded, Pre-Lubricated Catheter

- Pre-lubricated catheter (top-half)
- Side loaded catheter
- Smooth catheter eyelets help provide comfort upon insertion and removal

1500 mL Urine Collection Bag with Tear Strip

- Allows secure, discreet catheterization; use anywhere, anytime
- Simple to empty the collection bag and dispose with the tear strip

Additional Trusted Features

- Protective cap covers the introducer tip
- Not made with natural rubber latex

Hollister Incorporated

2000 Hollister Drive, Libertyville, IL USA

www.hollister.com



Apogee Plus 16" Intermittent Catheter System

Apogee Plus	Apogee Plus Kit [†]	Size	Description
—	B6FB	6 Fr	Firm
B8F	B8FB	8 Fr	Firm
B10F	B10FB	10 Fr	Firm
B12F	B12FB	12 Fr	Firm
B14F	B14FB	14 Fr	Firm
B16F	B16FB	16 Fr	Firm
B18F	B18FB	18 Fr	Firm
B12C	B12CB	12 Fr	Coudé
B14C	B14CB	14 Fr	Coudé
B16C	B16CB	16 Fr	Coudé
B12S	B12SB	12 Fr	Firm
B14S	B14SB	14 Fr	Firm
B14F-Female	B14FB-Female	14 Fr	Firm

[†]Kit Includes: Apogee Plus intermittent catheter system, gloves, underpad drape, gauze, and PVP swabsticks

Qty: 100/box

For detailed clinical questions concerning our products:

1.888.740.8999

For orders only: **1.800.323.4060**

For more information or to try a free sample, please visit us at www.hollister.com/apogee



Prior to use of Apogee Plus Intermittent Catheters, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions. Rx Only

The Hollister logo and Apogee Plus are trademarks of Hollister Incorporated.
© 2020 Hollister Incorporated. 922326-0420 US-01032

