



## VaPro Plus Pocket

Touch Free Hydrophilic Intermittent Catheter

**You don't have to compromise protection for discretion.**

Touch free catheter with an integrated collection bag in a pocket size package.

**The VaPro Plus Pocket touch free hydrophilic intermittent catheter combines the right balance of:**

### Ease of Use

Prehydrated and ready to use right out of the package

### Protection

The protective tip and sleeve help keep germs away

### Convenience

Both the catheter and collection bag are packaged in a small, discreet pouch that's easy to open



### Why the protective tip?

In many people, bacteria are concentrated in the first part of the urethra. The protective tip is designed to help reduce the risk of carrying these bacteria into your bladder during catheter insertion. Inserting the protective tip into your urethra and then passing the catheter through the end of the protective tip may help the catheter avoid these bacteria.

VaPro Plus Pocket Touch Free Hydrophilic Intermittent Catheter			
Stock No	Length (inches)	Size	Box Qty
71082-30	8	8Fr	30
71102-30	8	10Fr	30
71122-30	8	12Fr	30
71142-30	8	14Fr	30



For more information or to request a product sample, visit [hollisterpeoplefirst.com/vapro](http://hollisterpeoplefirst.com/vapro)



- Ready-to-use, single-use intermittent catheter
- Hydrophilic-coated PVC catheter
- Protective tip and sleeve supports no-touch insertion technique
- 2 smooth catheter eyelets
- Color-coded funnel

This intermittent catheter is a flexible tubular device that is inserted through the urethra by male, female, and pediatric patients who need to drain urine from the bladder.

**WARNING:** To help reduce the potential for infection and/or other complications, do not reuse. If discomfort or any sign of trauma occurs, discontinue use immediately and consult your healthcare practitioner.

**PRECAUTION:** Please consult a medical professional before using this product if any of the following conditions are present: Severed urethra, unexplained urethral bleeding, pronounced stricture, urethritis – inflammation of the urethra. Self-catheterization should follow the plan of care and advice given by your healthcare practitioner and be carried out only in accordance with the instructions for use provided. Because catheterization frequency varies by person, the recommended frequency of your catheterization should be provided by your healthcare practitioner. For any other questions about your catheterization, please contact your healthcare practitioner.

**CAUTION:** Rx Only. Federal (USA) Law restricts this device to sale by or on the order of a physician or other healthcare practitioner licensed under state law to order this product. Prior to use of this device, be sure to read (i) the complete information on how to use this device including Warnings, Cautions, and Instructions for Use, and (ii) all other package inserts and labels supplied with the product.



The Hollister logo and VaPro Plus Pocket are trademarks of Hollister Incorporated.  
©2015 Hollister Incorporated.  
922893-715

