



Continence Care Products



Making a difference in the journey of life.

Table of Contents

Hollister Navigator SM services.....	2-3
Intermittent Catheters	4-16
Infyna Chic TM Intermittent Catheters.....	4-5
VaPro TM No Touch Intermittent Catheters.....	6-9
Onli TM Intermittent Catheters	10-11
Advance Plus TM Intermittent Catheter Systems.....	12-13
Apogee Plus TM Intermittent Catheter Systems	14-16
Male External Catheters	17-18
InView TM Silicone Male External Catheters	17
Extended Wear TM Latex Male External Catheters.....	18
Everyday TM Latex Male External Catheters.....	18
Leg Bag Systems	19-20
Vented Urinary Leg Bag Combination Packs.....	19
Urinary Leg Bag Combination Packs	19
Urinary Leg Bags	20
Leg Bag Accessories	21
Fabric Leg Bag Straps.....	21
Extension Tubings	21
Collection Systems.....	22-24
Bedside Drainage Collection System with Anti-Reflux Valve.....	22
Male Urinary Pouch External Collection Devices	23
Female Urinary Pouch External Collection Device	24
Latex Statement.....	25
HCPCS Codes Disclaimer	25
Index	28



With Navigator services, no one has to journey alone.

Hollister is here to help you navigate the path ahead. From your first prescription to your first order, you don't have to figure it out on your own.



For more information or to enroll in NavigatorSM services visit www.hollister.com/navigator or contact the Navigator services team at: 888.712.6628 or navigator@hollister.com



Hollister offers free support for new intermittent catheter customers through NavigatorSM services.

NavigatorSM services help new Hollister intermittent catheter customers start strong and stay connected on the new journey ahead as an intermittent catheter user. Those enrolled in this exclusive program receive:



Instant DME Placement:

No matter the insurance, no matter the complexity, Navigator services will match Hollister customers with the right durable medical equipment (DME) company as soon as we receive enrollment.



Begin with Confidence:

New customers enrolled in Navigator services have access to free resources and experts to answer their intermittent catheter questions on **day 1 or day 100**.



A Guide Along the Way:

Changes happen, and they can be difficult to navigate alone when you rely on catheters daily. Hollister's Navigator services stays connected with customers and helps them overcome new challenges in their intermittent catheter journey.



Be a VIP:

Members are the first to know of Hollister's latest innovations that can help make intermittent catheter use easier, more comfortable, and improve quality of life.



Better outcomes start in the 'comfort zone' with Infyna Chic™ intermittent catheters

Comfort of clinical confidence: protection from potential UTI and urethral trauma.
Physical comfort: easy and convenient with HydraBalance™ Lubricating Technology.
Emotional comfort: overcoming barriers like fear and embarrassment.

Call to request a sample today!
1.888.808.7456



Infyna Chic™ Intermittent Catheters



Scan to
Request a
Free Sample

Infyna Chic™ Intermittent Catheters

81xx-30

Box of 30

- ▶ Discreet to Carry with Unique Combination of Color, Design and Small Size
- ▶ Flip-top Cap on Catheter Case is Easy to Open with Just One Hand
- ▶ Reseal Case After Use and Dispose of Without Leaks or Odor
- ▶ Ready-to-Use, Prehydrated Hydrophilic Catheter with Enhanced Lubrication
- ▶ Just Right Catheter Stiffness Helps Support Easy, Touch-Free Insertion
- ▶ Catheter Length Helps Give Confidence Bladder is Fully Drained

HCPCS: A4295

Stock No	Length (inches)	Size	Description
8108-30	5.5	8 Fr	Straight
8110-30	5.5	10 Fr	Straight
8112-30	5.5	12 Fr	Straight
8114-30	5.5	14 Fr	Straight



Making a difference in the journey of life.



With 100% No Touch Protection,
whatever you touch won't touch
VaPro™ catheters.


The protective tip and sleeve of VaPro™ catheters have been tested and have been shown to prevent the transmission of bacteria that may cause UTIs. Don't settle for anything less than 100% No Touch Protection.

Call to request a sample today!
1.888.808.7456



VaPro™
No Touch Intermittent Catheter

VaPro™ No Touch Intermittent Catheters



Scan to
Request a
Free Sample

VaPro Pocket™ Coudé No Touch Intermittent Catheters

77xxx-30

Box of 30

- ▶ Full Length No Touch Catheter in a Pocket-Size Package
- ▶ Protective Tip and Sleeve Allow for 100% No Touch Catheterization
- ▶ Coudé Tip
- ▶ Ready-to-Use, Prehydrated Hydrophilic Catheter
- ▶ Smooth Eyelets
- ▶ Not Made with Natural Rubber Latex

HPCPS: A4296

Stock No	Length (inches)	Size	Description
77124-30	16	12Fr	Coudé
77144-30	16	14Fr	Coudé
77164-30	16	16Fr	Coudé



Scan to
Request a
Free Sample

VaPro™ Coudé No Touch Intermittent Catheters

73xxx

Box of 30

- ▶ Protective Tip and Sleeve Allow for 100% No Touch Catheterization
- ▶ Ready-to-Use, Prehydrated Hydrophilic Catheter
- ▶ Smooth Eyelets
- ▶ Not Made with Natural Rubber Latex
- ▶ Coudé Tip

HPCPS: A4296

Stock No	Length (inches)	Size	Description
73124	16	12 Fr	Coudé
73144	16	14 Fr	Coudé
73164	16	16 Fr	Coudé



Making a difference in the journey of life.



Scan to
Request a
Free Sample

VaPro Pocket™ No Touch Intermittent Catheters

70xxx-30

Box of 30

- ▶ Full length No Touch Catheter in a Pocket-Size Package
- ▶ Protective Tip and Sleeve Allow for 100% No Touch Catheterization
- ▶ Ready-to-Use, Prehydrated Hydrophilic Catheter
- ▶ Smooth Eyelets
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4297

Stock No	Length (inches)	Size	Description
70104-30	16	10 Fr	Straight
70124-30	16	12 Fr	Straight
70144-30	16	14 Fr	Straight
70164-30	16	16 Fr	Straight



Scan to
Request a
Free Sample

VaPro Plus Pocket™ No Touch Intermittent Catheters

71xxx-30

Box of 30

- ▶ Full Length No Touch Catheter in a Pocket-Size Package
- ▶ Protective Tip and Sleeve Allow for 100% No Touch Catheterization
- ▶ Integrated 1000 mL Collection Bag
- ▶ Smooth Eyelets
- ▶ Ready-to-Use, Prehydrated Hydrophilic Catheter
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4297

Stock No	Length (inches)	Size	Description
71082-30	8	8Fr	Straight
71102-30	8	10Fr	Straight
71122-30	8	12Fr	Straight
71142-30	8	14Fr	Straight
71084-30	16	8Fr	Straight
71104-30	16	10Fr	Straight
71124-30	16	12Fr	Straight
71144-30	16	14Fr	Straight
71164-30	16	16Fr	Straight



Scan to
Request a
Free Sample

VaPro Plus™ No Touch Intermittent Catheters

74xxx

Box of 30

- ▶ Protective Tip and Sleeve Allows for 100% No Touch Catheterization
- ▶ Smooth Eyelets
- ▶ Integrated 1000 mL Collection Bag
- ▶ Ready-to-Use, Prehydrated Hydrophilic Catheter
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4297

Stock No	Length (inches)	Size	Description
74122	8	12 Fr	Straight
74142	8	14 Fr	Straight
74124	16	12 Fr	Straight
74144	16	14 Fr	Straight



DID YOU KNOW...

VaPro™ No Touch Intermittent Catheters offer the right balance of ease of use and protection:

- Prehydrated, hydrophilic and ready to use
- Protective sleeve allows for the catheter to be gripped anywhere and is proven to prevent the transmission of pathogens that may cause UTIs
- Protective tip helps shield the sterile catheter during insertion from bacteria located within the first 15mm of the distal urethra
- Packaging is designed with finger holes to be easy to open
- VaPro Plus™ and VaPro Plus Pocket™ catheters have an integrated 1000ml collection bag for catheterization at home or on the go
- Offered in pocket-size packaging for easy transportation and discreet out of home usage

Onli™ Intermittent Catheters



Scan to
Request a
Free Sample

Onli™ Intermittent Catheters

82xxx-30

Box of 30

- Ready-to-Use, Prehydrated Hydrophilic Catheter with Enhanced Lubrication
- Easy-to-Open, Peel-Apart Packaging with Two Finger Holes
- Ergonomic Sliding Gripper* and Length of Flexible Funnel Help Provide Control
- Catheter has a Balance of Flexibility and Rigidity that Helps Promote Easy, Touch-Free Insertion
- Clear Tube Provides Increased Visibility
- Doesn't Stain Common Fabrics** if Contents are Spilled

HCPCS: A4295

Stock No	Length (inches)	Size	Description
82081-30	7	8 Fr	Straight
82101-30	7	10 Fr	Straight
82121-30	7	12 Fr	Straight
82141-30	7	14 Fr	Straight
82084-30	16	8 Fr	Straight
82104-30	16	10 Fr	Straight
82124-30	16	12 Fr	Straight
82144-30	16	14 Fr	Straight
82164-30	16	16 Fr	Straight

Onli™ Coudé Intermittent Catheters

85xxx

Box of 30

- Ready-to-Use, Prehydrated Hydrophilic Catheter with Enhanced Lubrication
- Easy-to-Open, Peel-Apart Packaging with Two Finger Holes
- Coudé Tip
- Mess-Free, Doesn't Stain Common Fabrics**
- Ergonomic, Sliding Gripper* and Length of Flexible Funnel Help Provide Control

HCPCS: A4296

Stock No	Length (inches)	Size	Description
85124	16	12 Fr	Coudé
85144	16	14 Fr	Coudé
85164	16	16 Fr	Coudé

*Sliding gripper is only on 16" catheter

**Common fabrics include cotton, denim, linen and polyester



Easier to open:
two finger holes and
peel-apart packaging



Enhanced lubrication:
more lubricious hydrophilic
coating for easier insertion
and withdrawal



Clear catheter tube:
increased visibility and a
more appealing look



Inspired for easier use

Ready when you are

Onli™ intermittent catheters are designed with ease in mind by incorporating features that make a big difference to the self catheterization experience.

Call to request a sample today!
1.888.808.7456



Onli™
Intermittent Catheter



Making a difference in the journey of life.



Scan to
Request a
Free Sample

Advance Plus™ Intermittent Catheter Systems

94xxx/95xxx

Box of 100

- ▶ Gel Reservoir
- ▶ Protective Tip
- ▶ 1500 mL Collection Bag with Large Handle
- ▶ Smooth Eyelets
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4353

Stock No	Length (inches)	Size	Description
94064	16	6 Fr	Straight
94084	16	8 Fr	Straight
94104	16	10 Fr	Straight
94124	16	12 Fr	Straight
94144	16	14 Fr	Straight
94164	16	16 Fr	Straight
94184	16	18 Fr	Straight
95124	16	12 Fr	Coudé
95144	16	14 Fr	Coudé
95164	16	16 Fr	Coudé



Advance Plus™ Intermittent Catheter System Kits

96xxx/97xxx

Box of 100

- ▶ Advance Plus™ Catheter
- ▶ Kit Contains Gloves, Underpad Drape, and Benzalkonium Chloride Swabsticks
- ▶ Not Made with Natural Rubber Latex



Scan to
Request a
Free Sample

HCPCS: A4353

Stock No	Length (inches)	Size	Description
96064	16	6 Fr	Straight
96084	16	8 Fr	Straight
96104	16	10 Fr	Straight
96124	16	12 Fr	Straight
96144	16	14 Fr	Straight
96164	16	16 Fr	Straight
96184	16	18 Fr	Straight
97124	16	12 Fr	Coudé
97144	16	14 Fr	Coudé
97164	16	16 Fr	Coudé



Scan to
Request a
Free Sample

Advance Plus Pocket™ Intermittent Catheter Systems

93xxx

Box of 100

- ▶ Full Length Advance Plus™ Intermittent Catheter System in a Pocket-Size Package

HCPCS: A4353

Stock No	Length (inches)	Size	Description
93104	16	10 Fr	Pocket
93124	16	12 Fr	Pocket
93144	16	14 Fr	Pocket
93164	16	16 Fr	Pocket

**Apogee IC™ Intermittent Catheters****10xx/10xxx/11xxx****Box of 30**

- ▶ Smooth Eyelets
- ▶ Not Made with Natural Rubber Latex



Scan to
Request a
Free Sample

HCPCS: A4351 for all except**HCPCS: A4352 for xxx26**

Stock No	Length (inches)	Size	Description
10806	6	8 Fr	Firm
11006	6	10 Fr	Firm
11206	6	12 Fr	Firm
11406	6	14 Fr	Firm
11446*	6	14 Fr	Firm
10610	10	6 Fr	Firm
10810	10	8 Fr	Firm
11010	10	10 Fr	Firm
10616	16	6 Fr	Firm
10816	16	8 Fr	Firm
11016	16	10 Fr	Firm
11216	16	12 Fr	Firm
11416	16	14 Fr	Firm
11616	16	16 Fr	Firm
11816	16	18 Fr	Firm
10826	16	8 Fr	Coudé
11026	16	10 Fr	Coudé
11226	16	12 Fr	Coudé
11426	16	14 Fr	Coudé
11626	16	16 Fr	Coudé
11826	16	18 Fr	Coudé
1061	16	14 Fr	Soft
1065	16	14 Fr	Firm, Curved Package

*No funnel Qty: 30/box except curved catheters: 50/box



Scan to
Request a
Free Sample

Apogee Plus™ Intermittent Catheter Systems**Bxx/Bxxx****Box of 100**

- ▶ Protective Tip
- ▶ Pre-lubricated
- ▶ Smooth Eyelets
- ▶ 1500 mL Collection Bag
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4353

Stock No	Length (inches)	Size	Description
B8F	16	8 Fr	Firm
B10F	16	10 Fr	Firm
B12F	16	12 Fr	Firm
B14F	16	14 Fr	Firm
B16F	16	16 Fr	Firm
B18F	16	18 Fr	Firm
B12C	16	12 Fr	Coudé
B14C	16	14 Fr	Coudé
B16C	16	16 Fr	Coudé
B12S	16	12 Fr	Soft
B14S	16	14 Fr	Soft
B14F-Female	8	14 Fr	Firm



Apogee Plus™ Intermittent Catheter System Kits

BxxB/BxxxB

Box of 100

- ▶ Apogee IC™ Intermittent Catheter
- ▶ Kit Contains Gloves, Underpad Drape, Gauze, and PVP Swabsticks
- ▶ Not Made with Natural Rubber Latex



Scan to
Request a
Free Sample

HCPCS: A4353

Stock No	Length (inches)	Size	Description
B6FB	16	6 Fr	Firm
B8FB	16	8 Fr	Firm
B10FB	16	10 Fr	Firm
B12FB	16	12 Fr	Firm
B14FB	16	14 Fr	Firm
B16FB	16	16 Fr	Firm
B18FB	16	18 Fr	Firm
B12CB	16	12 Fr	Coudé
B14CB	16	14 Fr	Coudé
B16CB	16	16 Fr	Coudé
B12SB	16	12 Fr	Soft
B14SB	16	14 Fr	Soft
B14FB-Female	8	14 Fr	Firm



Scan to
Request a
Free Sample



InView™ Silicone Male External Catheters

974xx/975xx/976xx

974xx-100/975xx-100/976xx-100

Box of 30

Box of 100

- ▶ Self-Adhesive
- ▶ Clear Material Enables Monitoring of Skin
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4349

SPECIAL: shorter length sheath

Stock No Box of 30	Stock No Box of 100	Color	Size
97425	97425-100	Green	25 mm
97429	97429-100	Purple	29 mm
97432	97432-100	Brown	32 mm
97436	97436-100	Blue	36 mm
97441	97441-100	Orange	41 mm

STANDARD: the daily catheter for general use

Stock No Box of 30	Stock No Box of 100	Color	Size
97525	97525-100	Green	25 mm
97529	97529-100	Purple	29 mm
97532	97532-100	Brown	32 mm
97536	97536-100	Blue	36 mm
97541	97541-100	Orange	41 mm

EXTRA: 70% more adhesive than Standard for added security

Stock No Box of 30	Stock No Box of 100	Color	Size
97625	97625-100	Green	25 mm
97629	97629-100	Purple	29 mm
97632	97632-100	Brown	32 mm
97636	97636-100	Blue	36 mm
97641	97641-100	Orange	41 mm



Scan to
Request a
Free Sample

Extended Wear™ Latex Male External Catheters

920x

Box of 30

- ▶ Self-Adhesive
- ▶ Unique Inner Flap
- ▶ Easy Applicator Collar
- ▶ Contains Natural Rubber Latex

HCPCS: A4349

Stock No Box of 30	Color	Size
9206	■	22 mm-25 mm
9207	■	26 mm-30 mm
9208	■	31 mm-35 mm
9209	■	36 mm-39 mm

Caution: This product contains natural rubber latex which may cause allergic reactions.



Scan to
Request a
Free Sample

Everyday™ Latex Male External Catheters

910x

Box of 30

- ▶ Self-Adhesive
- ▶ Easy Applicator Collar
- ▶ Contains Natural Rubber Latex

HCPCS: A4349

Stock No Box of 30	Color	Size
9106	■	22 mm-25 mm
9107	■	26 mm-30 mm
9108	■	31 mm-35 mm
9109	■	36 mm-39 mm

Caution: This product contains natural rubber latex which may cause allergic reactions.



Scan to
Request a
Free Sample

Vented Urinary Leg Bag Combination Packs

9645/9655

Box of 10

- ▶ Includes:
 - (1) Sterile, Vented Urinary Leg Bag — Medium or Large
 - (1) Sterile, Extension Tubing with Vacuum Relief Valve, (18" Long, Oval, 11/32" Effective Inside Diameter), with Connector
 - (1) Pair, Fabric Leg Bag Straps — Medium (23" Long, 1-1/8" Wide)
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4358

Stock No	Description	Length (inches)	Width
9645	Medium — 18 oz (540 mL)	10"	4-1/4"
9655	Large — 30 oz (900 mL)	11-1/8"	5"



Scan to
Request a
Free Sample

Urinary Leg Bag Combination Packs

9348/9349

Box of 12

- ▶ Includes:
 - (1) Sterile, Urinary Leg Bag — Medium or Large
 - (1) Non-Sterile, Extension Tubing (18" Long, Oval, 11/32" Effective Inside Diameter) with Connector
 - (1) Pair, Fabric Leg Bag Straps — Medium (23" Long, 1-1/8" Wide)
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4358

Stock No	Description	Length (inches)	Width
9348	Medium — 18 oz (540 mL)	10"	4-1/4"
9349	Large — 30 oz (900 mL)	11-1/8"	5"



Scan to
Request a
Free Sample

Urinary Leg Bags

9814/9805

Box of 10

- Includes:
(1) Sterile, Urinary Leg Bag — Medium or Large
- Not Made with Natural Rubber Latex

HPCPS: A9270

HPCPS: A4358 (if straps are added)

Stock No	Description	Length (inches)	Width
9814	Medium — 18 oz (540 mL)	10"	4-1/4"
9805	Large — 30 oz (900 mL)	11-1/8"	5"



Scan to
Request a
Free Sample

Fabric Leg Bag Straps

9342/9343

10 Pair

- Not Made with Natural Rubber Latex

HPCPS: A5114

Stock No	Description	Length (inches)	Width
9342	Small	15"	1-1/8"
9343	Medium	23"	1-1/8"



Scan to
Request a
Free Sample

Extension Tubings

9345/9346

Box of 10

- Kink-resistant with Connector
Oval, 11/32" Effective Inside Diameter
- Thermoplastic Elastomer
- Not Made with Natural Rubber Latex

HPCPS: A4331

Stock No	Description	Length (inches)
9345	Non Sterile	18"
9346	Sterile	18"



DID YOU KNOW...

Hollister Leg Bags Offer Many Benefits:

- Comfort with soft, vinyl material
- Discreet profile with pleats
- Strap security with wide, woven fabric
- Minimal dexterity needed with easy-to-operate drain valve
- Reduced risk of pressure with protected strap connection buttons
- Helps decrease urine back flow with anti-reflux valve



Scan to
Request a
Free Sample

Bedside Drainage Collection System with Anti-Reflux Valve

9839

Box of 20

- ▶ Sample Port
- ▶ Sterile Fluid Pathway
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4357

Stock No	Description	Length (inches)	Width
9839	68 oz (2000 mL)	12"	9-3/4"



Male Urinary Pouch External Collection Devices*

9873

Box of 10

- ▶ Cut-to-Fit, FlexTend™ Skin Barrier
- ▶ Thin and Flexible
- ▶ No Additional Tape Needed
- ▶ Odor-Barrier Film
- ▶ Convenience Drain
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4326

Stock No	Description	Length (inches)	Width
9873	Male Urinary Pouch	7-1/2"	4-1/2"



Scan to
Request a
Free Sample

9811

Box of 10

- ▶ Cut-to-Fit, SoftFlex™ Skin Barrier
- ▶ Odor-Barrier Film
- ▶ Convenience Drain
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4326

Stock No	Description	Length (inches)	Width
9811	Male Urinary Pouch	7-1/2"	4-1/2"

**Meant to be connected to leg bag or gravity drainage*



Scan to
Request a
Free Sample

Female Urinary Pouch External Collection Device*

9840

Box of 10

- ▶ Cut-to-Fit, SoftFlex™ Skin Barrier
- ▶ Odor-Barrier Film
- ▶ Convenience Drain
- ▶ Alternative to Indwelling Catheter for Bedridden Females
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4328

Stock No	Description	Length (inches)	Width
9840	Female Urinary Pouch	7-1/2"	4-1/2"

**Meant to be connected to leg bag or gravity drainage*

HCPCS Codes Disclaimer

The reimbursement information provided herein is intended to provide general information concerning coding of Hollister products only. Hollister does not guarantee coverage or payment for any products. The ultimate responsibility for proper coding, satisfying reimbursement requirements, and obtaining reimbursement remains with the provider. Coding and coverage policies and guidelines are complex, can vary from one carrier or region to another, and are updated frequently. Providers should check with their local carriers or intermediaries often and should consult with counsel, a reimbursement specialist, and/or the DMECS website for any coding, coverage, reimbursement, or billing questions.

HCPCS codes shown in this Product Catalog are current as of January 2026.

Latex Statement

The specifications for the Hollister Continence Care products listed below include natural rubber latex as a component; our suppliers add natural rubber latex as part of the production process. If a Hollister product(s) found in this catalog is not listed below, the specification for that product(s) and packaging does not include natural rubber latex as a component, and our component suppliers do not add natural rubber as part of their production process.

Current Stock Numbers:

9106	9206
9107	9207
9108	9208
9109	9209

Prior to use of ALL continence care products, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

All intermittent catheters are Rx Only.

Stock No	Page
10xx	14
10xxx	14
11xxx	14
70xxx-30	8
71xxx-30	8
73xxx	7
74xxx-30	9
77xxx-30	7
81xx-30	5
82xxx-30	10
85xxx	10
910x	18
920x	18
93xxx	12
9342	21
9343	21
9345	21
9346	21
9348	19
9349	19

Stock No	Page
94xxx	12
95xxx	12
96xxx	13
9645	19
9655	19
97xxx	13
974xx	17
974xx-100	17
975xx	17
975xx-100	17
976xx	17
976xx-100	17
9805	20
9811	23
9814	20
9839	22
9840	24
9873	23
Bxx/Bxxx	15
BxxB/BxxxB	16

Contact Us!

Questions about our products or services?

1.888.808.7456

For orders only:

1.800.323.4060

Fax 1.847.680.1017

Hollister Incorporated

2000 Hollister Drive

Libertyville, Illinois 60048

www.hollister.com

Hollister, the Hollister logo, Advance Plus, Advance Plus Pocket, Apogee, Apogee IC, Apogee Plus, Everyday, Extended Wear, Flexend, HydraBalance, the HydraBalance logo, Infyna Chic, InView, Navigator, Onli, Softflex, VaPro, VaPro Plus, VaPro Plus Pocket and VaPro Pocket are trademarks of Hollister Incorporated.

© 2025 Hollister Incorporated. 922530-1225 US-02346

