

Female Urinary Pouch and Male Urinary Pouches

The Hollister Female Urinary Pouch and Male Urinary Pouches provide effective ways to manage urine output for non-ambulatory patients. The urinary pouches are not made with natural rubber latex and meant to be connected to leg or gravity drainage bags. The male urinary pouches also provide options for men when an external catheter is not an option.

Both products are alternatives to:

- Indwelling catheters, which may be associated with increased risk of urinary tract infections
- Absorbent pads and diapers, which may affect patient dignity, allow for compromised skin integrity, and may require frequent changes

Both products feature:

- Skin barrier that may be cut for a custom fit
- Odor-barrier film that helps manage odor effectively
- Universal connector that fits most leg or night drainage bag systems



Female Urinary Pouch



Male Urinary Pouch

Male Urinary Pouches

Male Urinary Pouch External Collection Device*

9873 **Box of 10**

- Cut-to-Fit, *Flexend* Skin Barrier
- Thin and Flexible
- No Additional Tape Needed
- Odor-Barrier Film
- Convenience Drain
- Not Made with Natural Rubber Latex

HCPCS: A4326

Stock No	Description	Length	Width
9873	Male Urinary Pouch	7-1/2"	4-1/2"

9811 **Box of 10**

- Cut-to-Fit, *SoftFlex* Skin Barrier
- Odor-Barrier Film
- Convenience Drain
- Tape Strips and Skin Gel Wipe Included
- Not Made with Natural Rubber Latex

HCPCS: A4326

Stock No	Description	Length	Width
9811	Male Urinary Pouch	7-1/2"	4-1/2"

*Meant to be connected to leg bag or gravity drainage

Female Urinary Pouch

Female Urinary Pouch External Collection Device*

9840 **Box of 10**

- Cut-to-Fit, *SoftFlex* Skin Barrier
- Odor-Barrier Film
- Convenience Drain
- Alternative to Indwelling Catheter for Bedridden Females
- *Adapt* Paste, Tape Strips, and Skin Gel Wipe Included
- Not Made with Natural Rubber Latex

HCPCS: A4328

Stock No	Description	Length	Width
9840	Female Urinary Pouch	7-1/2"	4-1/2"

*Meant to be connected to leg bag or gravity drainage

Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

Hollister, the Hollister logo, Adapt, Flexend and SoftFlex are trademarks of Hollister Incorporated.
© 2021 Hollister Incorporated. US-01457

