



VaPro Pocket™ Catheters

No Touch Protection now available with a coudé tip

VaPro Pocket™
No Touch Intermittent Catheter



VaPro Pocket™ Coudé No Touch Intermittent Catheters provide:



Protection.

Whatever you touch
won't touch VaPro™ catheters.



Protective tip

- Helps shield the sterile catheter during insertion from bacteria located within the first 15mm of the distal urethra
- Helps reduce the risk of carrying bacteria into the urinary tract

Protective sleeve

- Allows for catheter to be gripped anywhere
- Provides a barrier to the catheter to help keep germs away



Discretion.

Pocket-size packaging for easy transportation and
out-of-home use.

- Discreet packaging design
- Flexible design of the packaging allows catheter to feel comfortable in a pocket
- Does not leave stains on clothes*

*Polyester, cotton, linen, and denim



Ease of use.

- Prehydrated hydrophilic catheter is ready-to-use right out of the package - no extra steps required
- A finger hole makes the packaging easy to open

Packaging designed
with finger holes
to be easy to open.



Coudé tip.

- Designed to help maneuver around tight bends found in a narrow urethra
- Suitable for males who have enlarged prostate or similar condition

Scan to request
product sample



To learn more, call 1.888.808.7456
or visit www.Hollister.com/VaProCoude





“It’s all handy. I just have it in my bag or pocket and can go anywhere.”

— Alan, VaPro™ catheter user

Consider VaPro™ catheters right from the start

NEW! VaPro Pocket™ Coudé
No Touch Intermittent Catheter without Collection Bag

HCPCS A4352	Stock No	Size	Length	Tip	Box Qty
	77124-30	12 Fr	16"	Coudé	30
	77144-30	14 Fr	16"	Coudé	30
	77164-30	16 Fr	16"	Coudé	30



VaPro Pocket™

No Touch Intermittent Catheter



Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions. Rx Only.

Hollister Incorporated

2000 Hollister Drive
Libertyville, Illinois 60048
1.888.808.7456
www.hollister.com

HCPCS A4352: Insurance coverage based on medical necessity.

The testimonials, statements and opinions presented are applicable to the person depicted. These testimonials are representative of that person's experience, but the exact results and experience will be unique and individual to each person.

Hollister, the Hollister logo, VaPro, and VaPro Pocket are trademarks of Hollister Incorporated.
©2023 Hollister Incorporated. 923531-0123 US-01893

