



With 100% No Touch Protection,
whatever you touch won't touch
VaPro coude catheters.

The VaPro touch free hydrophilic intermittent catheters with coude tip are hygienic by design to help protect patients from pathogens throughout the entire catheterization process.




VaPro

Touch Free Hydrophilic
Intermittent Catheter

VaPro coudé tip catheters: 100% No Touch Protection

Protective tip

- Helps shield the sterile catheter during insertion from bacteria located within the first 15 mm of the distal urethra



Protective sleeve

- Allows for catheter to be gripped anywhere
- Provides a barrier to the catheter to help keep pathogens away

Request a sample

Contact your local rep or call **1.888.808.7456**

| Size | Length | Description | Color Code | Stock Number |
|-------|--------|-------------|------------|--------------|
| 12 Fr | 16" | Coudé | ○ White | 73124 |
| 14 Fr | 16" | Coudé | ● Green | 73144 |
| 16 Fr | 16" | Coudé | ● Orange | 73164 |

HCPCS A4352: Insurance coverage based on medical necessity.

Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

Rx Only

Hollister Incorporated

2000 Hollister Drive
Libertyville, Illinois 60048

www.hollister.com

The Hollister logo and VaPro are trademarks of Hollister Incorporated.
©2018 Hollister Incorporated.
923235-1218 US-00611

