



With 100% No Touch Protection, whatever you touch won't touch VaPro™ Coudé Catheters.

VaPro™ Coudé No Touch Intermittent Catheters provide 100% No Touch Protection with a protective tip and sleeve, and they are easy to use, with a coudé tip that is designed to help maneuver around the tight bends found in a narrow urethra.

Learn more by visiting Hollister.com or calling 1.888.808.7456



VaPro™
No Touch Intermittent Catheter




Making a difference in the journey of life.

VaPro™ Coudé Catheters: 100% No Touch Protection

Protective tip

- Helps shield the sterile catheter during insertion from bacteria located within the first 15 mm of the distal urethra



Protective sleeve

- Allows for catheter to be gripped anywhere
- Provides a barrier to the catheter to help keep pathogens away

VaPro™ Coudé No Touch Intermittent Catheter

Stock Number	Size	Length	Color Code	Box Qty
73124	12 Fr	16"	○ White	30
73144	14 Fr	16"	● Green	30
73164	16 Fr	16"	● Orange	30

HCPCS A4296: Insurance coverage based on medical necessity.

Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

Rx Only

Scan to request
product sample



For more information,
visit www.Hollister.com/Coude

Hollister Incorporated

2000 Hollister Drive
Libertyville, Illinois 60048
1.800.323.4060
www.hollister.com



VaPro™
No Touch Intermittent Catheter

