


Calcium Alginate Dressings

For Wound Care

Ordering Information

Stock No	Product Size	Qty/Box	HCPCS
 Restore Calcium Alginate Dressing, Silver Sterile			
529967	2" x 2" (5cm x 5cm)	10	A6196
529968	4" x 4.75" (10cm x 12cm)	10	A6197
529969	12" (30.5cm) rope	5	A6199
Restore Calcium Alginate Dressing Sterile			
529938	2" x 2" (5cm x 5cm)	10	A6196
529937	4" x 4" (10cm x 10cm)	10	A6196
529939	4" x 8" (10cm x 20cm)	5	A6197
529940	12" (30.5cm) rope	5	A6199



Complex Problems. Innovative Solutions.

Manufactured for
Hollister Wound Care
2000 Hollister Drive
Libertyville, IL 60048
1.888.740.8999

Distributed in Canada by
Hollister Limited
95 Mary Street
Aurora, Ontario L4G 1G3
1.800.263.7400

www.hollisterwoundcare.com

Hollisterwoundcare and wave logo are trademarks of Hollister Incorporated. Restore is a trademark of Hollister Wound Care.
Made and printed in USA. ©2013 Hollister Wound Care. 911037-0113



Restore

Calcium Alginate Dressings

High Absorbency without Sacrificing Strength



High wet strength allows intact, one-piece removal

EXUDATE LEVEL*



HEAVY

The highly absorptive composition of the dressing – made from calcium sodium alginate – creates a conformable, protective gel when in contact with exudate.

Benefits

- **Highly absorbent**
Uniformly absorbs exudates and protects the peri-wound skin
- **Forms and helps maintain a moist wound interface**
Creates a conformable, protective gel when in contact with exudate
- **Superior wet strength**
Intact, one-piece removal facilitates ease of removal
- **Less residue than other dressings**
Reduces need for repeated irrigation that might damage new tissue
- **Protects the wound**
Aids autolytic debridement, facilitating wound healing
- **Hemostatic Properties**
As this product contains alginate, it may assist in supporting the control of minor bleeding in superficial wounds

INDICATIONS FOR USE

- Restore Calcium Alginate Dressings can be used on a wide variety of moderate to heavily exuding partial to full-thickness wounds
- Can be used under compression bandages
- Also indicated for: Post-operative surgical wounds
1st and 2nd degree burns

Restore

Calcium Alginate Dressings, Silver

All the Benefits of our Calcium Alginate Dressings *Plus* the Benefits of Silver

Restore Calcium Alginate Dressings, Silver provide all the benefits of Restore Calcium Alginate Dressings *plus* the proven antimicrobial benefits of silver.

Benefits

- **Sustained Silver Release**
Silver ions are released in the presence of wound exudate
- **Antimicrobial Effect Protects Dressing**
Effective against a broad spectrum of micro-organisms over a period of up to 14 days, based on in-vitro testing (which is beyond indicated wear time) including: MRSA, VRE, Fungi
- **Odor reduction**
Results from the antibacterial effect in the dressing
- **Effective barrier** to bacterial penetration

INDICATIONS FOR USE

- Restore Calcium Alginate Dressings, Silver can be used on a wide variety of moderate to heavily exuding partial to full-thickness wounds
- Can be used under compression bandages
- Also indicated for: Post-operative surgical wounds
1st and 2nd degree burns
Wounds that are critically colonized




¹Data on file with Hollister Incorporated.


See IFU for complete information.

Caution: Federal (US) Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

RESTORE® Calcium Alginate Dressing, Silver absorbs 43% more than the leading silver hydrofiber dressing.¹



High Wet Strength



EXUDATE LEVEL*



HEAVY