


Advance Plus

Touch Free Intermittent Catheter System

Confidence. Pure and Simple.



Touch free 1500 mL collection bag with large handle

Helps you freely handle the catheter while helping to shield it from potential contamination

Protective tip

Helps the catheter bypass environmental bacteria

Ultra-smooth eyelets

For gentle insertion and removal

Gel reservoir

Lets you control catheter lubrication without additional supplies



The Advance Plus intermittent catheter system features a Touch Free technique enabling users to catheterize with confidence – virtually anytime, anywhere.

People First:

Raquel Escudero Villaverde

"Okay, I'm in a wheelchair. But when I travel and see the conditions in which many people live, I think, lucky me."




Advance Plus

Touch Free Intermittent Catheter System

Touch Free protective tip

- Helps protect catheter tip from exposure to bacteria in the initial 15 mm of the urethral opening
- Helps reduce the risk of carrying bacteria up into the urinary system
- Serves as a guide for proper and easy insertion



User-Regulated Gel Reservoir

- No additional lubrication required
- Stabilizes catheter for easy insertion and withdrawal

Latex- and phthalates-free PVC material

- Safe for users with latex-related allergies
- Environmentally responsible

Advance Plus Pocket Touch Free Intermittent Catheter System

- Designed for confidence and convenience
- The compact packaging helps make the catheter easier to carry and store
- The same detailed features as the Advance Plus intermittent catheter



Advance Plus Pocket Touch Free Intermittent Catheter System			
Stock No	Tip	Size (16 in, 40cm)	Box Qty
93104	Straight	10Fr	100
93124	Straight	12Fr	100
93144	Straight	14Fr	100
93164	Straight	16Fr	100

Advance Plus Touch Free Intermittent Catheter System			
Stock No	Tip	Size (16 in, 40cm)	Box Qty
94064	Straight	6Fr	100
94084	Straight	8Fr	100
94104	Straight	10Fr	100
94124	Straight	12Fr	100
94144	Straight	14Fr	100
94164	Straight	16Fr	100
94184	Straight	18Fr	100
95124	Coudé	12Fr	100
95144	Coudé	14Fr	100
95164	Coudé	16Fr	100

Advance Plus Kit Touch Free Intermittent Catheter System			
Stock No	Tip	Size (16 in, 40cm)	Box Qty
96064	Straight	6Fr	100
96084	Straight	8Fr	100
96104	Straight	10Fr	100
96124	Straight	12Fr	100
96144	Straight	14Fr	100
96164	Straight	16Fr	100
96184	Straight	18Fr	100
97124	Coudé	12Fr	100
97144	Coudé	14Fr	100
97164	Coudé	16Fr	100

Kit includes: Advance Plus Intermittent Catheter system, gloves, underpad drape, and Benzalkonium swabsticks

WARNING: To help reduce the potential for infection and/or other complications, do not reuse. If discomfort or any sign of trauma occurs, discontinue use immediately and consult your healthcare practitioner.

- Ready-to-use, single-use intermittent catheter
- Hydrophilic-coated PVC catheter
- Protective tip and sleeve supports no-touch insertion technique
- 2 smooth catheter eyelets
- Color-coded funnel
- Adhesive tab for adhering packaging to a surface

This intermittent catheter is a flexible tubular device that is inserted through the urethra by male, female, and pediatric patients who need to drain urine from the bladder.

PRECAUTION: Please consult a medical professional before using this product if any of the following conditions are present: Severed urethra, unexplained urethral bleeding, pronounced stricture, false passage, urethritis – inflammation of the urethra, prostatitis – inflammation of the prostate gland, epididymitis – inflammation of the epididymis (testicle tube). Self-catheterization should follow the plan of care and advice given by your healthcare practitioner and be carried out only in accordance with the instructions for use provided. Because catheterization frequency varies by person, the recommended frequency of your catheterization should be provided by your healthcare practitioner. For any other questions about your catheterization, please contact your healthcare practitioner.

CAUTION: Rx Only. Federal (USA) Law restricts this device to sale by or on the order of a physician or other healthcare practitioner licensed under state law to order this product. Prior to use of this device, be sure to read (i) the complete information on how to use this device including Warnings, Cautions, and Instructions for Use, and (ii) all other package inserts and labels supplied with the product.



Hollister Incorporated

Libertyville, IL USA

For detailed clinical questions concerning our products:

1.888.740.8999

For orders only: 1800.323.4060

**Distributed by
Hollister Limited**

Aurora, Ontario

1.800.263.7400

www.hollisterpeoplefirst.com

The Hollister logo, Advance Plus and Advance Plus Pocket are trademarks of Hollister Incorporated.

©2014 Hollister Incorporated.

922517-114