

## Apogee IC™ Intermittent Catheter



Apogee Plus™ Intermittent Catheter System



# Trusted Quality. Reliable Performance. Cost Effective.

## Apogee IC

### Intermittent Catheter

Apogee IC catheters are designed to help you catheterize easily.

- Smooth eyelets to aid comfortable insertion and removal **1**
- Flexible funnel helps promote easier handling **2**
- Guide notch on funnel indicates coude tip position

## Apogee Plus

### Intermittent Catheter System

The Apogee Plus catheter is a gel pre-lubricated, touch free catheter with integrated collection bag.

- Integrated 1500 mL collection bag supports touch free catheterization **3**
- Protective tip is designed to help bypass bacteria in the first 15mm of the distal urethra and reduce the risk of carrying bacteria up into the urinary system **4**
- Smooth catheter eyelets aid comfortable insertion and removal
- Available with and without insertion supplies: gloves, underpad drape, gauze and PVP swabsticks **5**

**Apogee™**  
Intermittent Catheters

### Apogee IC Intermittent Catheters

Stock No.	Length (inches)	Size	Description
10806	6	8 Fr	Firm
11006	6	10 Fr	Firm
11206	6	12 Fr	Firm
11406	6	14 Fr	Firm
11446	6	14 Fr	Firm
10610	10	6 Fr	Firm
10810	10	8 Fr	Firm
11010	10	10 Fr	Firm
10616	16	6 Fr	Firm
10816	16	8 Fr	Firm
11016	16	10 Fr	Firm
11216	16	12 Fr	Firm
11416	16	14 Fr	Firm
11616	16	16 Fr	Firm
11816	16	18 Fr	Firm
10826	16	8 Fr	Coudé
11026	16	10 Fr	Coudé
11226	16	12 Fr	Coudé
11426	16	14 Fr	Coudé
11626	16	16 Fr	Coudé
11826	16	18 Fr	Coudé
1061	16	14 Fr	Soft
1065	16	14 Fr	Firm, Curved Package

No funnel.

#### Hollister Incorporated

2000 Hollister Drive, Libertyville, IL USA

[www.hollister.com](http://www.hollister.com)

For detailed clinical questions concerning our products: **1.888.740.8999**

For orders only: **1.800.323.4060**

For more information or to try a free sample, please visit us at [www.hollister.com/apogee](http://www.hollister.com/apogee)

### Apogee Plus 16" Intermittent Catheter Systems

Apogee Plus	Apogee Plus Kit <sup>†</sup>	Size	Description
—	B6FB	6 Fr	Firm
B8F	B8FB	8 Fr	Firm
B10F	B10FB	10 Fr	Firm
B12F	B12FB	12 Fr	Firm
B14F	B14FB	14 Fr	Firm
B16F	B16FB	16 Fr	Firm
B18F	B18FB	18 Fr	Firm
B12C	B12CB	12 Fr	Coudé
B14C	B14CB	14 Fr	Coudé
B16C	B16CB	16 Fr	Coudé
B12S	B12SB	12 Fr	Soft
B14S	B14SB	14 Fr	Soft
B14F-Female	B14FB-Female	14 Fr, 6"	Firm

<sup>†</sup>Kit Includes: Apogee Plus intermittent catheter system, gloves, underpad drape, gauze, and PVP swabsticks

Qty: 100/box

Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions. Rx Only

Hollister, the Hollister logo, Apogee, Apogee IC and Apogee Plus are trademarks of Hollister Incorporated. © 2021 Hollister Incorporated. US-01467

